making work easy

de

en

fr

it

es

pt

zh

ja

ko

SIMPLEX

FAQ - frequently asked questions

Made in Germany

21-7150 10042024

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


FAQ – English

In the case of chapter references, see the full instructions for use available online.

en

Question / Error / Problem	Cause	Troubleshooting / Solution
<p>Unit displays a filament flow error. Either too little filament or no filament at all is extruded from the nozzle.</p>	<ul style="list-style-type: none"> Nozzle temperature is too low. The filament does not melt and cannot be extruded from the nozzle. When using TEMP1 and TEMP2 group filaments, always load / unload the filament with the TEMP2 setting. 	<ul style="list-style-type: none"> Only use the temperature specified for the filament. Relevant information is provided on the spool and on the packaging of the filament. Objects printed using SIMPLEX aligner model Filament must be sliced using the "SIMPLEX aligner model" profile.
<p>Nozzle is blocked</p>	<ul style="list-style-type: none"> If the temperature of the nozzle is too high, the filament can vitrify on the inside wall of the nozzle and harden. The flow of material may gradually come to a stop and the nozzle may become blocked. Ensure that the temperature for the filament is used for printing: <ul style="list-style-type: none"> TEMP1: SIMPLEX study model / working model / multi-use TEMP2: SIMPLEX aligner model. 	<ul style="list-style-type: none"> Step 1: <ul style="list-style-type: none"> Heat the nozzle to 255°. Select the EXTRUDER menu. Wait until the extruder has reached the required temperature. Press the "Load filament" button. This automatically results in extrusion of 15 cm of filament. Perform this action three times. Now select the nozzle temperature required for the filament used. Wait until it is reached. Press the "Load filament" button and check whether the blockage has been resolved. If step 1 is not successful, replace the nozzle.
<p>Nozzle is blocked</p>	<ul style="list-style-type: none"> If the filament is transported via the extruder into the hotend and then through the nozzle and dust and dirt also accumulate, there is also a risk of a blockage. Store the filament in a dry location away from dust. 	<ul style="list-style-type: none"> Step 1: <ul style="list-style-type: none"> Heat the nozzle to 255°. Select the EXTRUDER menu. Wait until the extruder has reached the required temperature. Press the "Load filament" button. This automatically results in extrusion of 15 cm of filament. Perform this action three times. Now select the nozzle temperature required for the filament used. Wait until it is reached. Press the "Load filament" button and check whether the blockage has been resolved. If step 1 is not successful, replace the nozzle.
<p>How long is the average service life of the nozzle?</p>		<ul style="list-style-type: none"> The service life depends on the filament used. <ul style="list-style-type: none"> As it contains plaster, SIMPLEX multi-use filament causes greater wear of the nozzle. SIMPLEX multi-use filament approx. 400 operating hours. SIMPLEX study model, SIMPLEX working model, SIMPLEX aligner model approx. 800 operating hours.
<p>How should the outside of the nozzle be cleaned?</p>		<ul style="list-style-type: none"> The outside of the nozzle can be cleaned using a thick, temperature-resistant cotton cloth. Do not use alcohol, pins, drills, or sharp-edged tools! Caution: Danger of burns! The dark color of the nozzle is a normal reaction due to the temperature of up to 250 °C.
<p>Filament is extruded between the nozzle and the hotend.</p>	<ul style="list-style-type: none"> The nozzle is not screwed on securely. 	<ul style="list-style-type: none"> The nozzle must be screwed on securely to the hotend.
<p>The nozzle has become blocked after replacing SIMPLEX study model filament with SIMPLEX aligner model filament or vice-versa.</p>	<ul style="list-style-type: none"> The temperature at which the filament was replaced was too low. 	<ul style="list-style-type: none"> TEMP2 must always be selected when replacing the filament. This is the only way to ensure that the previous filament is completely removed from the nozzle.
<p>Print object does not adhere to the print bed.</p>	<ul style="list-style-type: none"> The distance between the nozzle and the print bed is too large. 	<ul style="list-style-type: none"> Level the print bed.
	<ul style="list-style-type: none"> Kapton film of the print bed is dirty. 	<ul style="list-style-type: none"> Thoroughly clean the Kapton film with alcohol. The Kapton film may need to be replaced.

Question / Error / Problem	Cause	Troubleshooting / Solution
<p>Model loses adhesion to the print bed.</p> 	<ul style="list-style-type: none"> • Print bed not aligned / leveled. • The distance between the nozzle and the print bed has changed. Distance too large. If the nozzle is too far away from the print bed, the filament does not have any contact to the print bed during extrusion and, as a result, no adhesion to the print bed either. The model slips during printing. 	<ul style="list-style-type: none"> • Level the print bed.
<p>Unwanted lines are visible on the underside.</p>	<ul style="list-style-type: none"> • The print bed is not clean, for example as a result of fingerprints. Grease residue prevents adhesion of the filament to the print bed. 	<ul style="list-style-type: none"> • Clean the print bed with pure alcohol / spirits. • Please wear suitable gloves to clean the print bed. • Afterwards carefully examine the first print layer.
<p>The print job is interrupted.</p>	<ul style="list-style-type: none"> • Unfavorable ambient conditions, a draft, or cold air from air-conditioning. • Model not flatly aligned in the CAD software. • Contact surface to the print bed is insufficient for optimal adhesion. 	<ul style="list-style-type: none"> • Check the ambient conditions. <ul style="list-style-type: none"> - The SIMPLEX 3D filament printer may not be exposed to direct sunlight. - Protect against drafts. - Observe an ambient temperature of 15 - 32 °C. • At a room temperature of ≤ 21 °C, it is recommended that the print bed and nozzle are preheated before printing begins in order to compensate for differences in temperature between the room temperature and the printer. • Model must be aligned again to the zero position in the CAD software.
<p>The first layers do not adhere properly to the print bed and loosen. Unwanted lines are visible on the underside.</p>  	<ul style="list-style-type: none"> • The model was not placed flat on the virtual print bed in the SIMPLEX sliceware. Model only has limited contact with the print bed and therefore insufficient adhesion. 	<ul style="list-style-type: none"> • Ensure that the models are dropped correctly. • When enlarging, a gap should not be visible between the model and the print bed. The model must be placed flat and fitted on the virtual print bed. • Tip: Use drop mode in Renfert mode.
<p>Poor adhesion between the individual layers, gaps in the layers. Display may show filament feed errors.</p>	<ul style="list-style-type: none"> • Not enough filament is fed. 	<ul style="list-style-type: none"> • Clean the cogwheel of the extruder.

Question / Error / Problem	Cause	Troubleshooting / Solution
<p>Print error Side shift / layer shift The upper part of the print object has shifted laterally.</p> 	<ul style="list-style-type: none"> • Print bed was not repositioned correctly in the printer after removal and repositioning. • Latch on the build plate was not correctly closed. • The print bed moved during printing, thus impairing printing. 	<ul style="list-style-type: none"> • Before printing, ensure that the print bed is installed correctly. • Latch clasps must be closed to secure the print bed.
<p>Printing is interrupted with a filament flow error.</p>	<ul style="list-style-type: none"> • Filament does not roll off during printing. 	<ul style="list-style-type: none"> • The filament has become knotted. Roll off the filament up to the location where it is knotted and cut it off. Load the filament again and start printing.
<p>The extruder rattles when printing the first layer (first few minutes of printing).</p>	<ul style="list-style-type: none"> • Incorrect calibration of the print bed (nozzle too close to the print bed, filament cannot be extruded from the nozzle). 	<ul style="list-style-type: none"> • Level the print bed.
<p>Incomplete print / layer missing.</p>	<ul style="list-style-type: none"> • USB flash drive is faulty. • GCode files are not saved completely. 	<ul style="list-style-type: none"> • Test with another USB flash drive or replace USB flash drive.
	<ul style="list-style-type: none"> • USB stick removed during the save process. 	<ul style="list-style-type: none"> • To avoid data loss, it is important to safely remove the USB flash drive.
<p>Printer displays an error message "Short circuit of 12 V power supply" or printer crashes when the jack plug of the extractor fan is inserted.</p>	<ul style="list-style-type: none"> • If the plug of the extractor fan is not fully inserted or is inserted or disconnected during operation, the socket causes a short circuit. 	<ul style="list-style-type: none"> • The jack plug of the extractor fan may only be inserted or disconnected when it is switched off.
<p>What does the message "1000 operating hours" mean?</p>		<ul style="list-style-type: none"> • When 1000 operating hours have been reached, servicing is required in accordance with Chapter 11 of the instructions for use.
<p>Unit displays "Short circuit detected at the motor driver" as well as an ErrorCode.</p>	<ul style="list-style-type: none"> • Overload detected in a motor or motor is faulty. 	<ul style="list-style-type: none"> • Unit must be repaired by the service team.
<p>Printer displays "Error in GCode" if the height of the object to be printed is more than 100 mm.</p>		<ul style="list-style-type: none"> • Update sliceware to min. version 1.0.3.
<p>Printer displays "Error in GCode".</p>	<ul style="list-style-type: none"> • USB flash drive is faulty. G-Code files were not correctly saved. 	<ul style="list-style-type: none"> • Test with another USB flash drive or replace USB flash drive.
<p>Display remains dark, does not function after pressing the reset button.</p>	<ul style="list-style-type: none"> • If the reset button is pressed on the display or the circuit board, the firmware is deleted. 	<ul style="list-style-type: none"> • Firmware must be installed again.
<p>What happens if the customer loses or damages the spacer?</p>		<ul style="list-style-type: none"> • The spacer card is available as a spare part. Art. no. 900219048
<p>When does the Kapton film need to be replaced?</p>		<ul style="list-style-type: none"> • The Kapton film must be replaced if it is scratched / damaged / faulty or if adhesion of the models no longer occurs during printing.
<p>Printing via USB cable is aborted.</p>	<ul style="list-style-type: none"> • If the computer is operating at capacity with other tasks (e.g., scan software is open at the same time), USB transfer may be interrupted, and printing aborted. 	<ul style="list-style-type: none"> • Only print via USB cable if the computer is permanently available and not operating to 100 % capacity. Otherwise, use the USB flash drive.
<p>Is it possible to connect the printer using a LAN cable?</p>		<ul style="list-style-type: none"> • No, a LAN connection is not possible. The printer does not have a corresponding port.

Question / Error / Problem	Cause	Troubleshooting / Solution
How does the print data get from the PC to the printer?		<ul style="list-style-type: none"> • We generally recommend printing via USB flash drive: the G-Codes are saved from the PC to the USB flash drive. The USB flash drive must then be inserted into the USB port of the printer. • If the printer is located adjacent to the PC, it is possible to print using a USB cable: the USB flash drive is not required. During printing, the connected computer must not be shut down or switched to stand-by / energy-saving mode. • Printing via WiFi: Only recommended for WiFi networks with good coverage and few outages. The USB flash drive must be inserted in the USB port. The USB flash drive is used as an external storage medium.
Why is a WiFi connection not recommended?	<ul style="list-style-type: none"> • In the case of a poor WiFi connection, the connection between the PC and the SIMPLEX may be interrupted with the result that it cannot be ensured that SIMPLEX print can be used securely and reliably. 	<ul style="list-style-type: none"> • If data transfer via WiFi is not stable, we recommend printing using a USB flash drive.
What frequency band is used by the Simplex WiFi module?		<ul style="list-style-type: none"> • The WiFi module of the printer only operates at 2.4 Ghz.
Can acetone be used to clean the print bed?		<ul style="list-style-type: none"> • No, acetone would damage the Kapton. Corresponding cleaning agents are listed in the IFU.
What grease should be used to lubricate the printer?		<ul style="list-style-type: none"> • Grease is not included in the scope of delivery of the printer. • Generally, PTFE lubrication grease should always be used.
Connectivity of the Renfert CONNECT app and SIMPLEX 3D filament printer.		<ul style="list-style-type: none"> • The SIMPLEX printer and the CONNECT app cannot communicate directly with each other due to different software languages. • No notifications (e.g., print stop) are sent via the app. • For news or software updates, please visit our website at myRenfert.
How many models can be printed on average using one reel?		<ul style="list-style-type: none"> • On average, 40 - 50 models depending on the size and type of model.
Can the models be cleaned using acetone?		<ul style="list-style-type: none"> • Over time, acetone degrades PLA. If briefly wiped with acetone, models become matt.
Process for replacing filaments.		<ul style="list-style-type: none"> • If TEMP1 or TEMP2 filament is used, the filament must always be loaded or unloaded using TEMP2 settings. • For more information, see Chapter 7.3 in the instructions for use.
Does every filament have an expiry date and storage conditions?		<ul style="list-style-type: none"> • Our filament does not have an expiry date, but it should be stored in a closed container and dark location with silica gel in order to protect the filament against damp.
Is it possible to bond SIMPLEX filament models?		<ul style="list-style-type: none"> • Yes <ul style="list-style-type: none"> - using fast-acting adhesive: Using fast-acting adhesive can cause white spots to develop. - Hot-melt adhesive: Ensure that the resin begins to soften from 60 - 70 °C. It must be ensured that the areas to be bonded do not become too hot. - or using 2-part adhesive

Question / Error / Problem	Cause	Troubleshooting / Solution
Print quality not satisfactory, small bubbles visible on the print surface.	<ul style="list-style-type: none"> Filament has absorbed humidity from the ambient air. Slight stringing may be visible. Humidity should not exceed 60 % rH. 	<ul style="list-style-type: none"> Keep your filament dry. In case it has become humid, use a commercially available filament dryer box with the following properties: <ul style="list-style-type: none"> - active heating, - adjustable temperature, - filament outlet ideally on the top.
Why is a different filament required for a study model than for a working model?		<ul style="list-style-type: none"> The only difference between the SIMPLEX study model and the SIMPLEX working model is the color. Otherwise, the filament is the same. Why do we offer customers two colors? In orthodontics, study or diagnostic models are primarily white. These are used for planning and for communication with the patient. The appliances are fabricated by technicians on the working model. Working on a white model is tiring for the eyes. Green has a calming effect and feels natural. The green shade has proven its worth for fatigue-free vision and work. That's why we chose a green shade for working models.
Can resin be applied to the model? For example, for orthodontic reworking or in order to fabricate an orthodontic appliance on the printed model?		<ul style="list-style-type: none"> Yes, this is possible. Ensure that before use, a resin/resin separator is used. If a separator is not used, the applied resin bonds with the filament model and can no longer be lifted after curing.
Can the multi-use filament be milled? Would it have a different surface structure afterwards?		<ul style="list-style-type: none"> As it contains plaster, SIMPLEX multi-use can be easily worked with rotary instruments and a scalpel.
If the splints are very close to each other, is working with an intermediate film permitted?		<ul style="list-style-type: none"> Yes, this is possible.
If an intermediate film is used, is there a risk that the splints will not be seated tightly enough?		<ul style="list-style-type: none"> Observe the manufacturer recommendations of the splint material / system used.
Worry that in the case of slender or narrow tooth tips, the splint will not reach this point, in other words that it does not reach low enough.		<ul style="list-style-type: none"> Observe the manufacturer recommendations of the splint material / system used.
Does the model require insulation before thermoforming with a film?		<ul style="list-style-type: none"> That depends on what film is used for thermoforming. Recommendations for model insulation are provided in the instructions for use of the respective film. Observe the manufacturer recommendations.
Small holes in the film, heating.		<ul style="list-style-type: none"> Possible cause: Small particles of dust were present on the film that burned during heating. Ensure that the film is free from dust.

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