# Stabiloplast

No. 448 0004

PORTUGUÊS

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ENGLISH

## 1. Description/Applications

Stabiloplast is a highly viscous resin which quickly hardens after being sprayed with the Stabiloplast accelerator. Stabiloplast bonds to metal, wax, plastic and plaster proving useful for various techniques in dental technology. Stabiloplast is particularly suitable for stabilising the fine wax mesh and for reinforcing wax casting sprues. It is just as effective in adhering metals, e.g. preparation of soldering objects. Stabiloplast can be used as undercut wax, it is heat resistant, dimensionally stable, hard, can be ground and burns without residues. It is applied via the dose tip and can be modelled with the plastic spatula right up to the moment when it is sprayed with the accelerator.

## 2. Safety Instructions

Cyanoacrylate bonds skin and eyes in seconds. If the adhesive should come into contact with the skin, wash the area thoroughly with soap and water. Rubbing the affected area of skin repeatedly with oil or butter will accelerate the detachment of the glue. The adhesive should not be scratched or scraped off. This will cause damage to the skin. The glue will fall off by itself within 1–2 days due to natural skin exfoliation. It is not advisable to use a thinner to remove glue from skin. If the adhesive should come into contact with the eyes, immediately rinse with lukewarm water. It is then necessary to seek medical help for an ophthalmological examination.

## 3. Instructions for Use

To use, remove the protective lid and place the ready-for-use atomizer head upon it. Using the plastic spatula only apply the Stabiloplast in the desired form and thickness. Spray the accelerator with a distance of 15–20 cm on the Stabiloplast as otherwise the resin could form bubbles. Always spray in a well-ventilated area and away from your body. Stabiloplast then hardens within 10 to 60 seconds according to the layer thickness. After use keep the bottle upright to allow the resin to flow back into the bottle. Please do not cut the dose tip, as otherwise the seal by the closing mandrel can not be assured.

**Important:** After first use, always store the resin and accelerator separately to avoid any reaction between the accelerator and resin.

## 4. Storage

**Prior to first use**, keep in a cool and dark place at a constant temperature between 6–15°C. **Before opening for the first time**, allow the bottle to reach full room temperature slowly. **After opening the bottle** store at room temperature and keep away from direct sunlight.

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