making work easy

SILENT compact

FAQ - frequently asked questions

Made in Germany

- [FAQ - Deutsch](#)
- [FAQ - English](#)
- [FAQ - Français](#)
- [FAQ - Italiano](#)
- [FAQ - Español](#)
- [FAQ - Português](#)
- [FAQ - Türkçe](#)
- [FAQ - Русский](#)
- [FAQ - Polski](#)
- [FAQ - Dansk](#)
- [FAQ - Čeština](#)
- [FAQ - Українська](#)
- [FAQ – Svenska](#)
- [FAQ - Suomi](#)
- [FAQ - Românesc](#)
- [FAQ - 中文](#)
- [FAQ - 日本語](#)
- [FAQ - 한국어](#)
- [اللغة العربية - FAQ](#)


Fig. 1


Fig. 2

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FAQ - English

For information on operating elements, see [Fig. 1](#) and [Fig. 2](#) on page 2 of this document.

en

For chapter indications, refer to the complete instruction manual on the internet.

Trouble	Possible Cause	Corrective Action
When the device is switched on, there is a loud vibration noise for approx. 8 seconds.	<ul style="list-style-type: none"> The suction unit performs an automatic filter cleaning sequence. 	<ul style="list-style-type: none"> This is functional and requires no intervention.
During an interruption in operation, an acoustic signal occurs 3 times and then the filter cleaning function is carried out.	<ul style="list-style-type: none"> After 16 hrs. constant operation (turbine running time) without having been switched off in between, the filter cleaning function is carried out. 	<ul style="list-style-type: none"> Always switch appliances off at the on / off switch at the end of the working day (9, Fig. 1).
When the unit is switched on, the display dust drawer (24, Fig. 2) lights up and an acoustic signal occurs 3 times.	<ul style="list-style-type: none"> The time period for emptying the dust drawer has elapsed and the dust drawer has not yet been emptied. The emptying of the dust drawer has not been confirmed. 	<ul style="list-style-type: none"> Empty the dust drawer and confirm by pressing the key (25, Fig. 2). Empty the dust drawer and confirm by pressing the key (25, Fig. 2).
Display error (26, Fig. 2) lights up and an acoustic signal occurs 3 times.	<ul style="list-style-type: none"> The automatic switch-on function cannot be set. The difference between the chosen rpm at the handpiece and the stand-by current is too low. The automatic switch-on function cannot be set. 	<ul style="list-style-type: none"> Adjust the switch-on threshold to a higher speed on the handpiece. Use the unit in continuous operating mode.
The error message (26, Fig. 2) blinks.	<ul style="list-style-type: none"> The electronic device has become too hot. 	<ul style="list-style-type: none"> Switch the device off and allow it to cool down. Ensure sufficient cooling, e.g. with: <ul style="list-style-type: none"> - Chapter 4.2 installation instructions. - External ventilation (see chapter 4.6). - Change fine filter.
The suction performance is insufficient.	<ul style="list-style-type: none"> The set suction performance is too low. There is a blockage or leak in the suction tube. Dust drawer is not airtight. The fine filter is full. 	<ul style="list-style-type: none"> Select higher suction level. Check suction tube. Please observe the points in chapter 4.4 Check the dust drawer is in the correct position (chapter 6.2). Switch the unit off and back on in order to initiate the filter cleaning function. Change the fine filter (see chapter 6.3) (if there is no improvement in suction performance after the filter cleaning process).
Dust drawer is over full.	<ul style="list-style-type: none"> The time period for "empty dust drawer" has been set too high. 	<ul style="list-style-type: none"> Adjust to a lower time period (see chapter 6.2.1).
The signal to empty the dust drawer shows even though the dust drawer is not yet full.	<ul style="list-style-type: none"> The time period for "empty dust drawer" has been set too low. 	<ul style="list-style-type: none"> Adjust to a higher time period (see chapter 6.2.1).
The electronic device connected to the unit power socket cannot be operated.	<ul style="list-style-type: none"> Suction is not switched on. 	<ul style="list-style-type: none"> Switch the suction on.

Trouble	Possible Cause	Corrective Action
The suction starts even though the electronic device connected to the unit is not being used.	<ul style="list-style-type: none"> The automatic switch-on function has been set too low. Malfunction due to another device using the same power circuit. 	<ul style="list-style-type: none"> Setting the automatic switch-on (calibration). Connect the extraction unit and the device causing interference to different power circuits (phases).
The suction does not stop when the electronic device is switched off.	<ul style="list-style-type: none"> The automatic switch-on function has been set too low. Malfunction due to another device using the same power circuit. 	<ul style="list-style-type: none"> Setting the automatic switch-on (calibration) (see chapter 5.5). Connect the extraction unit and the device causing interference to different power circuits (phases).
The suction does not start even though the electronic device is connected to the extraction unit.	<ul style="list-style-type: none"> The automatic switch-on function has been set too high. Malfunction due to another device using the same power circuit. 	<ul style="list-style-type: none"> Setting the automatic switch-on (calibration) (see chapter 5.5). Connect the extraction unit and the device causing interference to different power circuits (phases).
The suction suddenly stops during automatic operation, even though the device connected at the appliance socket is still in use. The signal for automatic operation (20, Fig. 2) is still on.	<ul style="list-style-type: none"> Malfunction due to another device using the same power circuit. 	<ul style="list-style-type: none"> Connect the extraction unit and the device causing interference to different power circuits (phases).
In the automatic operation mode, the suction starts after an unusually long delay time.	<ul style="list-style-type: none"> Malfunction due to another device using the same power circuit. 	<ul style="list-style-type: none"> Connect the extraction unit and the device causing interference to different power circuits (phases).
The suction does not start after long transportation or after change of turbine.	<ul style="list-style-type: none"> The connection plug behind the service flap (13, Fig. 1) is not completely inserted. 	<ul style="list-style-type: none"> Open the service flap (13, Fig. 1) and push the connection plug fully up to the top.
When the unit is switched on, not all of the displays light up.	<ul style="list-style-type: none"> The display which does not light up is defect. 	<ul style="list-style-type: none"> Contact Renfert / Service.
The fuse is released when a connected electronic device is switched on.	<ul style="list-style-type: none"> The total capacity for the connected devices is too high. 	<ul style="list-style-type: none"> Observe the max. connected capacity (see chapter 8).
The suction process in the continuous operation or automatic operation stops abruptly and the relevant displays are still on (20 / 22, Fig. 2).	<ul style="list-style-type: none"> The suction turbine is overheated. The suction turbine is defect. 	<ul style="list-style-type: none"> Switch the unit off and allow to cool down for at least 60 min. Check whether the suction tube is blocked. Eliminate blockage. Switch the unit off and then on again, so that the filter cleaning process can be carried out. Change the fine filter (see chapter 6.3) (if the filter cleaning function does not improve the suction performance). Change the suction turbine.

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