making work easy

Twister / Twister venturi

Made in Germany

de

making work easy

Twister / Twister venturi

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Made in Germany

Inhalt

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making work easy

Twister / Twister venturi

TRANSLATION OF THE ORIGINAL INSTRUCTIONS FOR USE

Made in Germany

21-6603 18122023

Contents

1	Introduction	3
1.1	Employed Symbols	3
2	Safety	3
2.1	Intended use	3
2.2	Improper use	3
2.3	Ambient conditions (according to DIN EN 61010-1)	4
2.4	Hazards and warnings	4
2.5	Authorised persons	4
2.6	Disclaimer	4
3	Product description	5
3.1	Components and functional elements	5
3.2	Standard Delivery	5
3.3	Accessories	6
4	Commissioning	6
4.1	Setup Recommendations	6
4.2	Wall Mounting	6
4.3	Benchtop Models	7
4.4	Electrical connection	7
4.5	Compressed air supply	7
5	Operation	7
5.1	Key Symbols	7
5.2	Display	8
5.3	On / Off	8
5.3.1	Standby mode	8
5.4	Setting Mixing Parameters	8
5.5	Mixing Process	9
5.5.1	Abort mixing process	10
5.5.2	Power cut	10
5.6	Adjustment During the Mixing Process	10
6	Cleaning and Maintenance	10
6.1	Housing	10
6.2	Seal Surfaces	10
6.3	Changing suction filter system	11
6.3.1	Sponge filter	11
6.3.2	Sinter filter	11
6.4	Changing fuses	11
6.5	Changing inlet filter	12
6.6	Changing the silencer	12
6.7	Spare Parts	12
7	Troubleshooting	13
7.1	Twister (No. 1826)	13
7.2	Twister venturi (No. 1827)	13
7.3	Error codes	14
8	Technical data	14
9	Guarantee	15
10	Disposal instructions	15
10.1	Disposal instructions for countries in the EU	15

1 Introduction

1.1 Employed Symbols

This manual contains symbols also found on the appliance itself with the following meaning:

Danger!
Immediate risk of injury. Please observe accompanying documentation!

High Voltage!
Danger – high voltage.

Caution
Risk of damage to the appliance if the instruction is not followed.

Notice
A general notice that improves and eases use.

For indoor use only.

Disconnect from mains supply before opening the appliance.

The appliance complies with applicable EU directives.

**This product complies with the relevant UK legislation.
See UKCA Declaration of Conformity in the Internet under www.renfert.com.**

**Within the EU, this device is subject to the provisions of the directive 2002/96/EG
(WEEE directive).**

**There is a FAQ video available on this topic in our video gallery under
www.renfert.com/P49.**

► **List, particular attention should be paid**

- List
- List

⇒ **Instructions / appropriate action / input / operational sequence:**

You will be asked to carry out the action in a specified order.

♦ **Result of an action / reaction of the device / reaction of the program:**

The unit or program reacts as a result of your actions or when a specific incident occurs.

Other symbols are explained when displayed.

2 Safety

Please use this information to brief operators on how to use this appliance and about possible applications and potential hazards during use.

This user information should be made available to the operator at all times.

2.1 Intended use

The *Twister / Twister venturi* vacuum mixer is designed solely for homogeneous, bubble-free mixing of dental casting and modelling materials such as plasters, investments, and silicones.

2.2 Improper use

Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the device, increase the risk of serious injury and lead to damage to the environment or the device itself.

2.3 Ambient conditions (according to DIN EN 61010-1)

The appliance should only be operated:

- indoors
- up to 2,000 metres above sea level
- at an ambient temperature of 5 °C - 40 °C [41 °F - 104 °F] *)
- at a maximum relative humidity of 80 % at 31 °C [87.8 °F], linear reduction up to 50 % relative humidity at 40 °C [104 °F] *)
- with mains electricity supply provided that the voltage fluctuation is within 10 % of the rated value
- with Pollution Degree 2
- with Overvoltage Category II.

*) The appliance can be used at a temperature of 5 °C - 30 °C [41 °F - 86 °F] and at a humidity of up to 80 %. At temperatures of 31 °C - 40 °C [87.8 °F - 104 °F] the humidity must reduce proportionately to ensure that the appliance can be operated (e.g. at 35 °C [95 °F] = 65 % humidity, at 40 °C [104 °F] = 50 % humidity). The appliance should not be operated at temperatures above 40 °C [104 °F].

2.4 Hazards and warnings

- Granted protection is not given any longer in case the appliance is not operated according to the instruction manual on hand.
- Only to be used indoors. The appliance is only intended for dry operation and should not be used or stored outdoors or in wet conditions.
- The appliance should only be operated using a mains cable with a country-specific plug system. Any modification required should only be carried out by an electrician.
- The appliance should only be operated if the data on the rating plate corresponds with the data of regional mains voltage.
- The appliance should only be plugged into earthed sockets.
- The mains socket should be easily accessible.
- Regularly check connecting cables and hoses (e.g. mains cable) for damage (e.g. kinks, tears, porosity) or deterioration. Do not operate appliances with damaged mains wires, hoses or other defects.
- The appliance should always be operated under supervision.
- Risk of injury!
Risk of injury if the components used are not approved. Always use original Renfert components.
- Always read the safety data sheets supplied by the manufacturers when mixing investment materials! Wear the appropriate personal protective equipment as and when required!
- Never mix flammable or explosive materials.
- Never connect the paddle without the mixing bowl.
- Manipulating the automatic bowl coupler and the intake opening can result in damage to the device and personal injury.
- Switch off the appliance when work is complete.
- Switch off and unplug the appliance before carrying out repairs and servicing work on electrical components.
- Switch off and unplug the appliance before cleaning or servicing it.
- It is the responsibility of the operator that national regulations during operation and regarding a repeated safety inspection of electrical equipment are complied with. For Germany these are the regulation 3 by DGUV (German Statutory Accident Insurance) in relation with VDE 0701-0702 (Association for Electrical & Electronic Technology).
- Information on REACH and SVHC is available on our website www.renfert.com, in the Support area.

2.5 Authorised persons

The appliance may only be operated and serviced by trained personnel.

2.6 Disclaimer

Renfert GmbH is not liable for claims for compensation or claims under guarantee if:

- the product is used for purposes other than those stated in the operating instructions.
- the product is modified in any way – apart from modifications described in the operating instructions.
- the product has not been repaired by a specialist firm or original Renfert replacement parts have not been used.
- there is continued use of the product despite obvious safety defects or damage.
- the product has been subjected to mechanical knocks or has been dropped.

3 Product description

3.1 Components and functional elements

- | | |
|---|--|
| 1 Control knob (set mixing parameters, start, stop, aerate) | 9 Knurled nut |
| 2 Parameter keys | 10 Silencer (<i>Twister venturi only</i>) |
| 3 Display | 11 Compressed air connection (<i>Twister venturi only</i>) |
| 4 Sponge filter | 12 Inlet filter (<i>Twister venturi only</i>) |
| 5 Appliance switch | 13 Sinter filter |
| 6 Fuse | 14 Bowl incl. blunger |
| 7 Appliance socket | 15 Pneumatic connection set |
| 8 Cable safety catch | 16 Mains cable |
| | 17 Compressed air hose |


Fig. 1

3.2 Standard Delivery

- 1 Vacuum mixer
- 1 500 ml bowl, incl. blunger
- 1 Mains cable
- 1 Compressed air hose, 2 m (*Twister venturi only*)
- 1 Pneumatic connection set (*Twister venturi only*)
- 1 Operating instructions
- 1 Drill template
- 1 Fastening set
- 1 Replacement sponge filter

3.3 Accessories

- 1821-0102 Stand for benchtop device
W x H x D: 227 x 622 x 294 mm [8.937 x 24.488 x 11.575 inch]
- 1821-0200 Mixing spatula
- 1820-6500 Bowl incl. blunger, 65 ml
- 1820-6520 Bowl, 65 ml
- 1820-0200 Bowl incl. blunger, 200 ml
- 1820-0220 Bowl, 200 ml
- 1820-0500 Bowl incl. blunger, 500 ml
- 1820-0520 Bowl, 500 ml
- 1820-0700 Bowl incl. blunger, 700 ml
- 1820-0720 Bowl, 700 ml
- 1820-1001 Bowl incl. blunger, 1000 ml
- 1820-1020 Bowl, 1000 ml
- 1823-0500 Alginate Mixing bowl 500 ml, incl. paddle
- 2929-0000 Filter pressure regulator

4 Commissioning

4.1 Setup Recommendations

Operate the appliance at room temperature 15 - 30 °C [59 – 86 °F].

When setting up the appliance, please note the following:

- Do not place the appliance under a heat source.
- Do not place the appliance in front of open windows.
- Do not set the appliance up where it will be subjected to direct sunlight.
- Do not subject the appliance to high humidity.

4.2 Wall Mounting

Please have the following items on hand:

- pencil,
- drill template,
- Phillips head screwdriver,
- drill
- drill bit 8 mm [5/16th inch] - appropriate for the wall material.

Make sure the wall on which you intend to mount the device is capable of bearing its weight!

Please ensure that there are no electric cables or water pipes running through the wall at the place of installation!

- ⇒ Establish optimum working height.
- ⇒ Align the drill template and mark the drill holes.
- ⇒ Drill the required holes >>> depth at least 55 mm [2.2 inch].
- ⇒ Insert the dowels.
- ⇒ Screw the top screws into the dowels. Allow the screws to protrude 9 - 10 mm [0.35 - 0.4 inch].
- ⇒ Screw the lower mounting bracket into place.
- ⇒ Hang the appliance on the wall.
- ⇒ Secure the device by tightening the knurled nut (9, Fig. 3).


Fig. 2

4.3 Benchtop Models

The stand (18) can be used to convert the vacuum mixer to a floor-mounted appliance.

- ⇒ Position the stand on a level surface.
- ⇒ Hang the appliance in the stand.
- ⇒ Secure the appliance by tightening the knurled nut (9).

Tighten the knurled nut securely.

The stand is not part of the standard delivery and must be ordered as an optional accessory (refer to section „3.3 Accessories“).


Fig. 3

4.4 Electrical connection

Ensure the mains supply and voltage marked on the nameplate is identical.

- ⇒ Secure the mains cable (16) through the cable safety catch (8).
- ⇒ Plug mains cable into socket (7) on appliance.
- ⇒ Plug in at the mains.


Fig. 4

4.5 Compressed air supply

Twister venturi only

- ⇒ Plug the compressed air hose (17) into the tube attachment on the compressed air connection (11) as far as it will go. There will be a slight resistance to overcome when attaching.
- ⇒ Select the proper coupling from the set (15, Fig. 1) provided and install it on the other hose end to connect it to the compressed air system.
- ⇒ Connect to compressed air line system.


Fig. 5-a

In previous devices, the inlet filter (12) was situated outside the device (Fig. 5-b). In these cases, push the compressed air hose (17) onto the inlet filter (12). There will be a slight resistance to overcome when attaching.

NEVER operate the device without the inlet filter!

The best vacuum performance is attained at the operating pressure specified in the technical data (see section 8). Use filter pressure regulator if necessary (see section „3.3 Accessories“).


Fig. 5-b

The vacuum mixer is now ready for operation.

5 Operation

5.1 Key Symbols

Mixing parameters	Symbol	Adjustment range	Factory setting
Mixing time		0:00 - 9:55 min:sec	1:00
Speed		100 - 450 1/min	350
Vacuum		70 - 100 *) 80 / 100 **) %	100

*) Twister, no. 1826: in 5 % increments

**) Twister venturi, no. 1827: Choice of 80 % and 100 % vacuum levels only.

5.2 Display

The following will appear on the display:

- 19 Large-scale display of parameter (here: mixing time)
- 20 Parameter symbols; The symbol „vacuum“ only shows when less than 100 % vacuum has been set
- 21 Selected speed
- 22 Selected vacuum level
- 23 Current vacuum level (bar graph)


Fig. 6

5.3 On / Off

The appliance is switched on and off at the switch (5, Fig. 1).

The mixing parameters last used will reappear in the display when the appliance is switched on.

The factory settings will appear the first time the appliance is switched on (see section 5.1).

5.3.1 Standby mode

The appliance will go into standby mode if no button is pressed or it stands idle for more than three minutes. The display will be dimmed when this happens.

To leave the standby mode:

- ⇒ Press any parameter key.
 - ⇒ Turn the control knob.
 - ⇒ Press the control knob.
 - ⇒ Couple the mixing bowl to the appliance.
- ◆ The mixing parameters last used will reappear in the display.

5.4 Setting Mixing Parameters

Procedure for changing a mixing parameter:

- ⇒ Press parameter key (2).
- ◆ The parameter key symbol is shown in the display.
- ◆ The parameter value is displayed in large-scale characters.
- ⇒ Turn the control knob (1) to change the value.
- ◆ The modified values are saved when the mixing bowl is coupled to the device.

i In initial state the mixing time is always displayed in large-scale characters and can also be changed directly without pressing the parameter key.

TIP - Setting the Vacuum:

Tests with various investment materials have shown that the best mixing results (smooth, homogeneous casting results) are obtained at the maximum vacuum setting. As a rule, this also applies to plasters. In individual instances, independent tests have shown that - under extremely high vacuum – the partial pressure in the mixing bowl can sink to a point where bubbles due to boiling can form in certain plasters.

Then reduce the vacuum setting.


Fig. 7

5.5 Mixing Process

Always read the safety data sheets supplied by the manufacturers when mixing investment materials! Wear the appropriate personal protective equipment as and when required!

Observe maximum mark on mixing bowl. Do not fill the mixing bowl over the maximum mark. This limit is valid for powders and liquids measured before mixing. If the mixing bowl is over-filled the suction filter system can become soiled.

Attempting to mix small quantities in large bowls will result in inadequate mixing.

- ⇒ Set mixing parameters as instructed by the manufacturer.
- ⇒ Select an appropriate mixing bowl size.
- ⇒ Mix powder and liquid as instructed by the manufacturer.
Allow plaster to sit for approx. 15 - 20 sec.
- ⇒ Using a spatula, lightly mix the materials together.
- ⇒ Install the appropriate paddle.
- ⇒ Make sure the bowl rim and cover rim (25, Fig. 10) are clean.
- ⇒ Couple the mixing bowl (14) to the appliance.
 - ◆ When the bowl is added the motor runs briefly in order to engage the paddle.
 - ◆ Vacuum generation starts automatically.

- ⇒ Do not let go until the bowl stays in place alone and vacuum level bar (23) has passed the halfway mark!

If the set vacuum is < 100 %, the pump will switch off when it is achieved. The pump can restart briefly, in order to reach the set vacuum level exactly.

Twister venturi only

The reduced vacuum level (80 %) is achieved through an opened bypass. This is why it sounds similar to when 100 % vacuum is achieved.

- ⇒ Start the mixing process: Press the control knob (1).
- ◆ The remaining mixing time (countdown) is displayed during the mixing process.

By delaying the start of the mixing process, a prevacuum can be achieved.

- ◆ The appliance beeps when the mixing cycle has finished.
- ◆ The time which has elapsed since the end of the mixing cycle is shown in the display.
- ⇒ Hold the bowl!
- ⇒ Aerate bowl: Press the control knob (1).
 - ◆ The bowl will be released in a few seconds.
- ⇒ Remove the mixing bowl.

Tip:

Only mix the same type of material in a given bowl. Residue from previous mixing processes can have a negative influence (e.g., silicon fails to harden, etc.).

We recommend a separate bowl for each type of material (plaster, investment, silicon).

Use the adhesive labels enclosed to the mixing bowl.


Fig. 8


Fig. 9

5.5.1 Abort mixing process

⇒ Stop the mixing process: Press the control knob once.

- ◆ The appliance will beep.
- ◆ Mixing process being aborted

⇒ Hold the bowl!

⇒ Aerate bowl: Press the control knob.

- ◆ The bowl will be released in a few seconds.

⇒ Remove the mixing bowl.

5.5.2 Power cut

Twister only

If a power cut should occur or the appliance is switched off during a mixing operation, the vacuum will be maintained and the bowl will stay on the device.

! When the power returns or the device is switched back on again, the bowl will become ventilated and then fall off.

Twister venturi only

! In case of a power failure or if the appliance is switched off, the bowl will be vented and be uncoupled.

5.6 Adjustment During the Mixing Process

All the mixing parameters can be displayed in large-scale characters during the mixing process if you press the relevant parameter key.

i The parameters can also be changed during the mixing process by pressing the relevant parameter key and turning the control knob.

i Changes made to the mixing parameter settings during the mixing process will be applied to that mixing process only and will not be saved.

6 Cleaning and Maintenance

Switch off and unplug the appliance before cleaning or servicing.

6.1 Housing

Do not steam clean the appliance.

Do not use solvent-based, aggressive or scouring cleaners.

⇒ Switch the appliance off.

⇒ Disconnect from mains.

⇒ Wipe housing with a damp cloth.

6.2 Seal Surfaces

The following seal surfaces must always be kept clean in order to ensure correct vacuum built up and the secure retention of the mixing bowl during the mixing process:

- Seal between the paddle and the bowl (25).
- Seal between the appliance and the paddle (26).

Tip

You should regularly coat the rubber seals with Vaseline.

This will significantly increase their service life and optimize to appliance's function


Fig. 10

6.3 Changing suction filter system

The suction filter system consists of a sponge filter (4) and a sinter filter (13).

⚠ When using investment materials, sponge filters and sintered filters can contain substances (dusts) that are hazardous to health. Please refer to the information in the safety data sheet. Use appropriate personal protective equipment.

! Vacuum can only build-up in a clean suction filter system.

! Never operate the device without the complete suction filter!


Fig. 11

6.3.1 Sponge filter

Pull the sponge filter (4) out in the direction of the base and replace with a new sponge filter.


Fig. 12

6.3.2 Sinter filter

- ⇒ Pull the sponge filter out in the direction of the base (Fig. 12).
- ⇒ Loosen the sinter filter (13) by turning counter-clockwise and remove.
- ⇒ Insert new sinter filter and tighten clockwise. The sinter filter should fit tightly
- ⇒ Replace the sponge filter.

! In previous devices (serial number index A) the sinter filter was only plugged in. In this case, pull out the sinter filter downwards and push the new sintered filter in until it stops.

**i The sinter filter can also be cleaned in the ultrasonic bath.
(Recommendation: Plaster solvent GO-2011, Art.no. 2011-0000).**


Fig. 13

6.4 Changing fuses

⚠ Before replacing a fuse, disconnect the appliance from the mains.

⚠ Never use a fuse of larger impedance.

- ⇒ Switch the appliance off.
- ⇒ Disconnect from mains.
- ⇒ Unlatch the fuse holder (6) on both sides and remove.
- ⇒ Change blown fuse.
- ⇒ Replace the fuse holder by sliding back into the socket until both sides lock into place.


Fig. 14

6.5 Changing inlet filter

Twister venturi only

In the Twister venturi units there are two different types of inlet filters:

► **Current version:**

- Inlet filter internally, behind the compressed air connection (Fig. 15-a).**
- ⇒ Disconnect the appliance from the compressed air supply.
 - ⇒ Press the compressed air connection ring (11) upwards and remove the compressed air tube (17).
 - ⇒ Loosen the compressed air connection screw fitting (11) and remove the filter (12).
 - ⇒ Screw the new filter into the compressed air connection, replace the screw fitting and hand-tighten.
 - ⇒ Plug the compressed air tube (17) into the tube attachment on the compressed air connection (11) as far as it will go, passing the point of slight resistance.
 - ⇒ Reconnect the appliance to the compressed air supply.


Fig. 15-a

► **Previous version:**

Inlet filter hidden, behind the compressed air connection (Fig. 15-b).

- ⇒ Disconnect the appliance from the compressed air supply.
- ⇒ Press the ring on the inlet filter (12) up and pull the compressed air hose (17) off.
- ⇒ Press the ring on the compressed air coupling (11) up and pull the inlet filter (12) off.
- ⇒ Insert a new inlet filter to the compressed air coupling (11), passing the point of slight resistance.
- ⇒ Insert the compressed air hose to the new intake filter, passing the point of slight resistance.
- ⇒ Reconnect the appliance to the compressed air supply.


Fig. 15-b

Frequent inlet filter blockage is an indication of contaminated compressed air. In such cases we recommend the upstream installation of a filter pressure regulator (see section „3.3 Accessories“).

6.6 Changing the silencer

Twister venturi only

- ⇒ Unscrew and remove silencer (10).
- ⇒ Fit and screw in new silencer.


Fig. 16

6.7 Spare Parts

You can find components subject to wear and the spare parts on the spare part list in the internet at www.renfert.com/p918.

Enter the following item number: 18260000 or 18270000.

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

7 Troubleshooting

7.1 Twister (No. 1826)

en

Error	Cause	Solution
ON / OFF switch fails to work.	<ul style="list-style-type: none">• Power not connected.• Device fuse blown.	<ul style="list-style-type: none">• Check the power supply.• Replace the fuse (Sec. 6.4).
Motor fails to start.	<ul style="list-style-type: none">• Motor fault.	<ul style="list-style-type: none">• Have the device repaired.
No, or insufficient vacuum (not in the green range on the scale).	<ul style="list-style-type: none">• Suction filter system plugged.• Seal surfaces dirty.• Solenoid valve fault.• Vacuum pump fault.	<ul style="list-style-type: none">• Have filter cleaned or changed (Sec. 6.3).• Clean the seal surfaces (Sec. 6.2).• Have the device repaired.• Have the device repaired.
Vacuum release / venting, too slow.	<ul style="list-style-type: none">• Suction filter system plugged.• Solenoid valve fault.	<ul style="list-style-type: none">• Have filter cleaned or changed (Sec. 6.3).• Have the device repaired.
Vacuum pump doesn't start.	<ul style="list-style-type: none">• Bowl detection pin sticking.	<ul style="list-style-type: none">• Clean area around pin.• Have the device repaired.
The vacuum pump switches on and off intermittently during the mixing process.	<ul style="list-style-type: none">• At < 100 % vacuum, the vacuum pump will switch off once the set vacuum has been reached.	<ul style="list-style-type: none">• None, if the vacuum should drop due to seepage, the vacuum pump will automatically switch back on until the set vacuum level has been reached.

7.2 Twister venturi (No. 1827)

Error	Cause	Solution
ON / OFF switch fails to work.	<ul style="list-style-type: none">• Power not connected.• Device fuse blown.	<ul style="list-style-type: none">• Check the power supply.• Replace the fuse (Sec. 6.4).
Motor fails to start.	<ul style="list-style-type: none">• Motor fault.	<ul style="list-style-type: none">• Have the device repaired.
No, or insufficient vacuum.	<ul style="list-style-type: none">• Compressor switched off.• Operating pressure too low or too high.• Compressed air hose not connected, leaky or bent.• Air hose too long.• Square are of air hose too small.• Suction filter system plugged.• Seal surfaces dirty.• Inlet filter plugged.• Silencer plugged.• Solenoid valve fault.• Venturi nozzle plugged.	<ul style="list-style-type: none">• Switch the compressor on.• Check dynamic operating pressure (see technical data, Sec. 8).• Switch off other consumers.• Check the hose.• Maximum length 2 m.• Minimum inner diameter 4 mm.• Have filter cleaned or changed (Sec. 6.3).• Clean the seal surfaces (Sec. 6.2).• Change inlet filter (Sec. 6.5).• Replace the silencer (Sec. 6.6).• Connect filter pressure regulator in series.• Have the device repaired.• Have the device repaired.
Vacuum release / venting, too slow.	<ul style="list-style-type: none">• Suction filter system plugged.• Silencer plugged.• Solenoid valve fault.	<ul style="list-style-type: none">• Have filter cleaned or changed (Sec. 6.3).• Replace the silencer. (Sec. 6.6).• Connect filter pressure regulator in series.• Have the device repaired.
No vacuum being generated.	<ul style="list-style-type: none">• Bowl detection pin sticking.	<ul style="list-style-type: none">• Clean area around pin.• Have the device repaired.
Permanent venting noise.	<ul style="list-style-type: none">• Solenoid valve fault.• Bowl detection pin sticking.	<ul style="list-style-type: none">• Have the device repaired.• Clean area around pin.• Have the device repaired.
Fluttering noise from the silencer.	<ul style="list-style-type: none">• Working pressure too low.	<ul style="list-style-type: none">• Check dynamic operating pressure (see technical data, Sec. 8).

7.3 Error codes

The following indications will be given if the electronic system detects an error:

- ◆ The mixing operation is canceled.
- ◆ A warning signal will sound.
- ◆ The display will flash alternately „Err“ and an error code.

End of error display:

- ⇒ Press control button (except for Err 3, 5 and 6).
 ⇒ Err 3, 5 and 6: Switch device off, follow advice in section 5.5.2 “power failure”.

With Err 5 it may no longer be possible to remove the mixing bowl.

With an error code as shown in the following table, please proceed as advised.

Error code	Cause	Solution
Err. 2	<ul style="list-style-type: none"> • Minimum vacuum level not yet reached. • Vacuum level has dropped below 500 mbar. 	<ul style="list-style-type: none"> • Ventilate bowl and remove (see section 5.5.1, points no. 2 and 3). • If the error recurs: Have the device repaired.
Err. 3	<ul style="list-style-type: none"> • The vacuum is built-up too quickly due to the fact that the suction system is blocked. 	<ul style="list-style-type: none"> • Have filter cleaned or changed (Sec. 6.3).
Err. 4	<ul style="list-style-type: none"> • Motor fault. 	<ul style="list-style-type: none"> • Ventilate bowl and remove (see Sec. 5.5.1, points no. 2 and 3). • Have the device repaired.
Err. 10	<ul style="list-style-type: none"> • Too much material. • Mixing time too long, the material has already begun to set. 	<ul style="list-style-type: none"> • Ventilate bowl and remove (see section 5.5.1, points no. 2 and 3). • Only fill the bowl to the maximum fill level. This limit is valid for powders and liquids measured before mixing. • Ventilate bowl and remove (see section 5.5.1, points no. 2 and 3). • Select shorter mixing times.

For all other error codes:

- ⇒ Note error code.
 ⇒ Have the device repaired.
 ⇒ Quote the error code to the repair service.

8 Technical data

	Twister (No. 1826)	Twister venturi (No. 1827)
Mains voltage:	100 - 240 V, 50 / 60 Hz	
Power consumption:	180 VA	
Device fuse:	T4AL, 250 VAC	
Speed:	100 - 450 1/min	
Connection pressure:	--	4,5 - 6 bar [65 - 87 psi]
Air consumption, approx:	--	28 l/min
LpA *) (idling):	< 70 db(A)	
Dimensions (W x H x D):	152 x 285 x 235 mm [6.0 x 11.2 x 9.3 inch]	152 x 320 x 235 mm [6.0 x 12.6 x 9.3 inch]
Weight (without bowl) approx:	5,2 kg [11.5 lbs]	4,0 kg [8.8 lbs]

*) sound pressure level in accord. with EN ISO 11202

9 Guarantee

Renfert gives a **3-year guarantee** on all parts of the *Twister / Twister venturi* provided that the product is used in accordance with the operating instructions.

en

The original sales invoice of the specialist depot is required for a claim under guarantee.

Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the guarantee. These parts are marked in the spare part list.

The guarantee is rendered void in the case of incorrect use, non-adherence to the operating, cleaning, servicing or connection instructions, repairs carried out by the owner or repairs that are not carried out by a specialist firm, use of replacement parts from another manufacturer or unusual or unapproved uses not specified in the user guidelines.

Successful claims under guarantee do not extend the guarantee period.

10 Disposal instructions

The appliance should be disposed of by a specialist firm. The specialist firm should be informed of any hazardous residue in the appliance.

10.1 Disposal instructions for countries in the EU

To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic appliances for proper disposal or recycling.

Within the European Union appliances with this symbol should not therefore be disposed of in unsorted domestic waste.

Within the European Union appliances with this symbol should not therefore be disposed of in unsorted domestic waste.

For more information regarding proper disposal please apply at your local authorities.

Hochaktuell und ausführlich auf ...
Up to date and in detail at ...
Actualisé et détaillé sous ...
Aggiornato e dettagliato su ...
La máxima actualidad y detalle en ...
Актуально и подробно на ...

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declare under our sole responsibility, that the product(s)

Name	Part No.
Twister	1826 0000 / 1826 0000 03
Twister venturi	1827 0000 / 1827 0000 03
Twister evolution	1828 0000 / 1828 0000 03

is (are) in conformity with the relevant regulatory requirements by compliance with the UK designated standards.

UK legislation

Electromagnetic Compatibility Regulations 2016

Supply of Machinery (Safety) Regulations 2008

The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

Designated standards

EN 61326-1:2013

EN ISO 12100:2010

EN 61010-1:2010/A1:2019/AC:2019-04

EN IEC 61010-2-051:2021/A11:2021

EN IEC 63000:2018

Tilo Burgbacher,
Engineering Director

Hilzingen, 02.11.2022

The Technical Documentation has been retained by Renfert GmbH

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c/o Renfert GmbH
Renfert GmbH, Industriegebiet, 78247 Hilzingen / Germany