THE.R.MO. BRIDGE

In Ø 1.75 / 2.5 mm filament spool

de|en|fr|it|es|no|sv|sr|sl|cs

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HERSTELLER:
Pressing Dental s.r.l – Via Onesto Scavino, 4 – 47891 Falciano (Republik San Marino)

Monticelli Giovanni Alessandro – Via Emilia Romagna 233 – 47841 Cattolica (RN) Email: info@pressing-dental.com

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CONTENTS:
Acrylic resin in polymethyl methacrylate (PMMA) for prototyping with the FDM process, intend- ed for the manufacture of dental prostheses.

LOT:
The lot number is indicated on the packaging.

STORAGE:
Use within the use-by date indicated on the container.
Use within the use-by date indicated on the container. Store the product in its packaging. Do not expose the product to direct sunlight. Keep away from sources of heat, in a dry place. Follow the storage conditions indicated on the packaging.

INTENDED USE:
Material intended for printing of dentistry components, such as temporary crowns and bridg- es, using the FDM prototyping process. The transparent material can also be used to pro- duce surgical guides.

SPECIFIC INSTRUCTIONS FOR USE:
The.r.mo.Bridge® can be used with the specific 3D Tecno Master and SIMPLEX 2 SX equip- ment manufactured by Pressing Dental S.r.l. and Renfert GmbH.
Specifications for processing with 3D Tecno Master and SIMPLEX 2 SX equipment:
Processing program with relative temperature, pressure, cooling and flow rate specifications pre-set in the software.
Printing temperature: May vary depending on the type of product and color, ranging from 260 °C to 295 °C (Check the operating temperature on the label placed on the cellophane wrapping of each roll/package).
Use of the material on equipment not manufactured by Pressing or Renfert: Not contem- plated.
In order to learn how to use the device and the materials correctly, the user (dentist and/ or dental technician) must take a theoretical/practical course at the manufacturer or at au- thorized dealers.
During said training course, the user is taught all the currently known processing techniques for correct use and the benefits and limits of the product; instructions are also provided on specific use of the equipment that is essential for creating the prosthesis.
At the same time, the user is also provided with all the instructions needed to be able to use the product in combination with other devices.
If it is necessary to combine the The.r.mo. Bridge® resin with metal parts, it is essential to create mechanical retentions, as there is no adhesion between the materials.

RECOMMENDED COLOUR SCALE:
The.r.mo.Bridge® | VITA classical

FINISHING AND POLISHING:
Use non-aggressive polishing pastes e.g.: Renfert Polish 3D-printed materials.
Only use cotton brushes at low-rpm to avoid overheating of the material (e.g.: a handpiece brush with a diameter of 20 mm on a 20,000/40,000 rpm handpiece, polish at intervals of 5/10 seconds, exerting a slight pressure; also use the same procedure for brushes with a diameter of 80 mm and with 2800 rpm equipment).
The prosthesis must only be cleaned with running water and possibly with the use of ma- terials that are not too abrasive. Do not use ultrasound equipment with acids that could al- ter its characteristics.
Prior to use on patients, special disinfection procedures must be performed on surgical ga- ges.
If this material is used for the manufacture of bites, possible surface abrasion caused by con- tact with the antagonists must be considered.
The colour of the polymerised resin could change in the following cases:
• if, during the extrusion phase, the melting temperature indicated on the specific product is exceeded;
• if it is not polished properly.

WARNINGS:
The product is not resistant to chemical agents with an oxidizing effect and strong acids (pH <4). During finishing, use a suction unit and/or mask to avoid inhaling the dust.
During finishing and subsequent polishing, take extreme care not to overheat the product, as this is a thermoplastic material and there could therefore be an unintended alteration of the prosthesis and a release of irritant gases.
When manufacturing small devices or prostheses, ensure that the dimensions of the item are not such as to risk swallowing if it moves.
Attention: The material is not visible to X-rays.
Do not reuse the product.
Do not exceed the maximum melting temperature.
Do not change the pre-set parameters in the software
Do not exceed the melting time, as this could cause mechanical alterations of the material with probable consequences on the final prosthesis.
Do not contaminate the product during the prototyping phase.
Do not mix the material with other materials during processing.
In the case of allergic reactions to the material, suspend application immediately and seek medical advice.

H332	Harmful if inhaled	(intended for gases in case of thermal decom- position)
H317	May cause an allergic skin reaction	(intended for dust during finishing steps)
H314	Causes severe skin burns and eye damage	(intended for contact with the molten polymer)
H320	Causes eye irritation	(intended for dust during finishing steps)
H335	May cause respiratory irritation	(intended for dust during finishing steps)

Precautionary statements – General
P 101 If medical advice is needed, have product container or label at hand.

Precautionary statements – Prevention
P 202 Do not handle until all safety precautions have been read and understood.
P 210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources.
No smoking.
P235 Keep cool.
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P270 Do not eat, drink or smoke when using this product.
P281 Use personal protective equipment as required.
P235 + P410 Keep cool. Protect from sunlight.

Precautionary statements – Response
P312 Call a POISON CENTER or doctor/physician if you feel unwell.
P301 + P312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
P332 + P313 If skin irritation occurs: Get medical advice/attention.
P337 + P313 If eye irritation persists: Get medical advice/attention.

Precautionary statements – Storage
P410 + P412 Protect from sunlight. Do not expose to temperatures exceeding 45 °C/113 °F.
EUH401 To avoid risks to human health and the environment, comply with the instructions for use.

INFORMATION FOR THE FINAL USER:
When delivering the finished device, the dental technician or dentist must provide the user with at least the following instructions:
a) Do not clean the prosthesis with abrasive products or products intended for cleaning acrylic resin or metal prostheses; use only oral hygiene products.
b) Wash the prosthesis with cold water only and at a temperature lower than 42 °C in any case.
c) Make sure to perform normal oral hygiene procedures.
d) Inform the patient that the prosthesis is not radiopaque.
e) Teeth grinding and contact with abnormal antagonists cause wear of the material.
f) Inflammable device (synthetic resin).
g) In the case of allergic reactions to the material, suspend application immediately and seek medical advice.
h) Report to the dentist or dental technician or directly to Pressing Dental any serious inci- dent involving the resin of which the dental prosthesis is made

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FABRICANT :
Pressing Dental s.r.l – Via Onesto Scavino, 4 – 47891 Falciano (République de Saint-Marino)

Monticelli Giovanni Alessandro – Via Emilia Romagna, 233 – 47841 Cattolica (RI) e-mail : info@pressing-dental.com

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PRESSING DENTAL S.r.l.