making work easy

de
en
fr
it
es

SIMPLEX sliceware

Made in Germany

21-9051 04042025 - Vers. 1.1.4

making work easy

SIMPLEX sliceware

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
Contents

en

1	Introduction	3
1.1	Symbols used	3
2	Safety	3
3	Installation	3
4	Renfert mode	4
4.1	Starting SIMPLEX sliceware	4
4.2	Loading 3D objects	4
4.3	Editing 3D objects	5
4.3.1	View / zoom	5
4.3.2	Adding / removing 3D objects	6
4.3.3	Undoing / redoing actions	6
4.3.4	File functions	6
4.3.5	Moving 3D objects	7
4.3.6	Rotating 3D objects	8
4.3.7	Dropping 3D objects	9
4.3.8	Copying 3D objects	9
4.3.9	Scaling 3D objects	10
4.3.10	Optional functions	11
4.4	Slicing	11
4.5	USB	13
4.6	Options	13
5	Expert mode	13
6	Print	14

1 Introduction

SIMPLEX sliceware was developed with the goal of simplifying and perfecting the 3D printing process with filaments.

To print a 3D object, the SIMPLEX 3D filament printer requires tool path information (G-codes).

With SIMPLEX sliceware, this G-code is generated from 3D object files (e.g., *.stl or *.obj).

SIMPLEX sliceware features two modes, Renfert mode and Expert mode.

1.1 Symbols used

Symbols with the following meanings are used in these instructions or on the equipment:

Caution

Risk of damage to the equipment if the instruction is not followed.

Notice

Provides the operator with useful information that eases use.

- List, particular attention should be paid

- List
- Subordinate list

⇒ Instructions / required action / input / operational sequence:

You are asked to carry out the specified action in the order indicated.

- ◆ Result of an action / response by the equipment / response by the program:

The equipment or program responds as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur.

2 Safety

Before using SIMPLEX sliceware, read this user manual carefully and keep it for future reference. Renfert GmbH reserves the right to modify or revise this user manual at its own discretion at any time.

3 Download and Installation

In order to download and install the SIMPLEX sliceware program, you must register yourself and your SIMPLEX 3D filament printer with myRenfert.

[Login / Registration with myRenfert:](#)

The procedure is described in the Quick Start Guide supplied with your SIMPLEX 3D filament printer.

[You can also open the Quick Start Guide via this link:](#)

4 Renfert mode

4.1 Starting SIMPLEX sliceware

When you start SIMPLEX sliceware, Renfert mode opens automatically.

An empty virtual print bed that corresponds to the actual print bed of the SIMPLEX 3D filament printer is shown.


Fig. 1

4.2 Loading 3D objects

⇒ Load a 3D object using one of the following methods:

Icon	Shortcut	Name	Description
	Ins	Add	Loads a 3D object in the *.stl, *.ply, or *.obj file format using an open file dialog.
	----	Project	Opens the "Project" menu.
	CTRL-O	Open	Loads 3D objects from a previously saved project.

♦ The 3D object is loaded and displayed on the print bed.


Fig. 2

4.3 Editing 3D objects

The following actions are possible at any time using the mouse:

- Increasing / decreasing the size of the print bed view with or without a 3D object.
⇒ Scroll using the mouse wheel.
- Rotating the print bed with or without a 3D object.
⇒ Hold down the right mouse button and move the mouse.
- Moving the print bed with or without a 3D object.
⇒ Hold down the center mouse button (mouse wheel) and move the mouse.
- Selecting 3D objects
⇒ Click the 3D object with the left mouse button.
 - ♦ The 3D object is selected.
 - ♦ The 3D object is shown in color.
 - ♦ The space that the 3D object takes up is shown.

To select multiple 3D objects, click them with the Ctrl key pressed.

- Move the selected 3D objects in the X / Y direction.
⇒ While holding down the left mouse button, move the mouse to reposition the selected 3D objects.
- In doing so, if the 3D object is moved beyond the edge of the print bed, the color changes.

Additional actions are possible using the icons at the top of the screen.

4.3.1 View / zoom

In the "View" menu, you can choose preset perspectives and an enlarged display.

Icon	Shortcut	Name	Description
	End	View	Opens the menu for selecting different perspectives.
	F2	Home	Standard view
	F3	Top	Top view of the print bed and model
	F4	Maximum	Preset magnification
	F5	Left	Left view of the print bed and model
	F6	Front	Front view of the print bed and model
	F7	Right	Right view of the print bed and model
	F8	Bottom	Bottom view of the print bed and model
	F9	Back	Rear view of the print bed and model

If you no longer have a good overview when zooming or moving using the mouse, you can easily restore the standard view using the "View" menu and "Standard" or using the F2 function key.

If function keys are already reserved by other programs on the computer, some of the actions mentioned above may not be possible.

4.3.2 Adding / removing 3D objects

Icon	Shortcut	Name	Description
	Ins	Add	Loads a 3D object in the *.stl, *.ply, or *.obj file format using an open file dialog.
	Del	Remove	Removes all selected 3D objects.

You can place multiple 3D objects on the print bed.

All selected 3D objects must be printed using the same filament!

4.3.3 Undoing / redoing actions

Icon	Shortcut	Name	Description
	Ctrl-Z	Undo	The previous steps are reversed.
	Ctrl-Y	Redo	The previous steps are restored.

4.3.4 File functions

Icon	Shortcut	Name	Description
	---	Project - Menu	Opens the "Project" menu.
	Ctrl-N	New	Creates a new project.
	Ctrl-O	Open	Loads an existing SIMPLEX sliceware project (*.ssp file).
	Ctrl-S	Save	Saves the project.
	Ctrl-Shift-S	Save as	Saves the project under a new name.

4.3.5 Moving 3D objects

Icon	Shortcut	Name	Description
	M	Move	3D objects can be moved to any position.

⇒ Select one or more 3D objects

◆ The "Move" option can be selected.

⇒ Click "Move".

◆ An information window opens in which the move action is shown for the three spatial directions X, Y, and Z.


Fig. 3

◆ Three directional arrows are shown for the three possible directions of movement.


Fig. 4

⇒ Hover the mouse over the control point (sphere) of one of the directional arrows.

◆ The direction of movement is selected.

◆ The other directional arrows disappear.

⇒ Click the control point of the directional arrow with the left mouse button and move the mouse with the mouse button held down.

◆ The 3D object is moved in the direction of the arrow.


Fig. 5

If the 3D object is placed beyond or below the print bed during the movement, its color changes.

When moving the 3D object, the movement is shown in the different directions in the information window.

The 3D object can also be moved by entering a value or using the up/down keys.

If you have selected multiple 3D objects, they are all moved together.


Fig. 6

Additional actions are possible using the buttons in the information window.

Icon	Shortcut	Name	Description
		Reset	The 3D object is reset to the position it had before the movement was carried out.
		Drop	The 3D object is moved at the current position in the current orientation on to the print bed (Z plane).
		Realign	The 3D object is moved in the current orientation at the center of the print bed on to the print bed (Z plane).

4.3.6 Rotating 3D objects

Icon	Shortcut	Name	Description
	R	Rotate	Rotates 3D objects.

⇒ Select one or more 3D objects

- ◆ The "Rotate" option can be selected.

⇒ Click "Rotate".

- ◆ An information window opens in which the angle of rotation is shown around the three axes of rotation X, Y, and Z.


Fig. 7

- ◆ Three circles are shown with directional arrows for the three possible axes of rotation.


Fig. 8

⇒ Hover the mouse over the control point (sphere) of one of the directional arrows.

- ◆ The axis of rotation is selected and the circle is shown in white.

⇒ Click the control point of the directional arrow with the left mouse button and move the mouse with the mouse button held down.

- ◆ The 3D object is rotated along the circle.


Fig. 9

If the 3D object is placed beyond or below the print bed during the rotation, its color changes.

When rotating the 3D object, the angle of rotation around the different axes is shown in the information window.

The 3D object can also be rotated by entering a value or using the up/down keys.

If you have selected multiple 3D objects, they are all rotated together.



Fig. 10

Icon	Shortcut	Name	Description
		Reset	The 3D object is reset to the orientation it had before the rotation was carried out.
		Drop	The 3D object is moved at the current position in the current orientation on to the print bed (Z plane).
		angle snap	The orientation changes by this angular value when the mouse button is moved while holding down the Shift key.

4.3.7 Dropping 3D objects

Icon	Shortcut	Name	Description
	D	Drop	The 3D object is dropped with the selected surface on the print bed (Z plane).

⇒ Select a 3D object.

◆ The "Drop" option can be selected.

⇒ Click "Drop".

◆ The mouse cursor changes to a crosshair.

⇒ Using the left mouse button, click the surface of the 3D object with which the 3D object is to be placed on to the print bed.

To finish:

⇒ Click "Drop" again.


Fig. 11

To make it easier to select the surface, you can rotate the 3D object and print bed to a favorable orientation by holding down the right mouse button, or select an appropriate view from the "View" menu.


4.3.8 Copying 3D objects

Icon	Shortcut	Name	Description
  		More	Selects additional options.
	C	Copy	The selected 3D object is copied and placed on the print bed.

If you have selected several 3D objects, all the selected 3D objects are copied.

4.3.9 Scaling 3D objects

Icon	Shortcut	Name	Description
		More	Selects additional options.
	S	Scale	The selected 3D object can be scaled.

⇒ Select one or more 3D objects

- ♦ The "Scale" option can be selected.

⇒ Click "Scale".

- ♦ An information window opens in which the original size is shown in mm. This value corresponds to 100%.


Fig. 12

- ♦ A grid with control points (spheres) is shown for scaling in the X and Y directions (red/green) as well as for proportional scaling in all three directions (orange).


Fig. 13

⇒ Hover the mouse over the control point (sphere) of one of the scaling directions.

- ♦ The scaling direction is selected.
- ♦ The other scaling options disappear.

⇒ Click the control point with the left mouse button and move the mouse with the mouse button held down.

- ♦ The size of the 3D object is increased / decreased in accordance with the scaling direction.


Fig. 14

If the 3D object extends beyond the print bed when the size is increased, its color changes.

During scaling, the size of the 3D object and the percentage modification are shown in the information window.

Scaling in a particular direction can also be performed by entering a value or using the up/down keys.

If the "Uniform" checkbox is selected, the size will be changed in all directions even if a value was entered for only one direction.



Fig. 15

Additional actions are possible using the buttons in the information window.

Icon	Shortcut	Name	Description
		Reset	The 3D object is reset to the original size it had before scaling was carried out.
		Drop	The 3D object is moved at the current position in the current orientation on to the print bed (Z plane).
		Uniform	Scaling is performed simultaneously in all three directions.

4.3.10 Optional functions

The following functions are still being tested and developed.

en

Icon	Shortcut	Name	Description
  		More	Selects additional options.
	P	Add Plate	A support plate is placed on the print bed.
		Add Support	A support is placed between a selected location and the print bed.

The support plates can be selected, moved, and deleted again just like a 3D object.

Supports can be used to support overhangs during printing.

To work with supports:

⇒ Select a 3D object.

◆ Supports can be selected.

⇒ Click "More" - "Add Support".

◆ The mouse cursor changes to a crosshair.

⇒ Using the left mouse button, click the surface of the 3D object where a support should be placed linking to the print bed.

Several supports can be placed into position.

To finish:

⇒ Click "More" - "Add Support" again.

Once you have exited support mode, the supports can be selected, moved, and deleted again just like a 3D object.

4.4 Slicing

To print a 3D object, what is known as a G-code must be generated for the 3D object.

This process is known as slicing.

Icon	Shortcut	Name	Description
	Ctrl+Shift+L	Slice	Perform slicing, generate a G-code, and save it.

◆ There is at least one 3D object present on the print bed and placed on it.

◆ The "Slice" option can be selected.

⇒ Click "Slice".

◆ A window opens for filament/profile selection.

⇒ Select the filament and corresponding profile.

- SIMPLEX study model
- SIMPLEX working model
- SIMPLEX aligner model
- SIMPLEX multi-use model

◆ A window opens for selecting the valve and filament generation.

⇒ Select the valve that is to be used.

- Brass valve
- Steel valve

⇒ Select filament generation (Fig. 16).

- 1: Current filaments
- 2: Future filaments

When using multiple screens, these windows may open on a different screen.

⇒ Confirm your selection with OK.

- ♦ A "Save as" window opens.

⇒ Enter a file name to save the G-code. Select a different directory if required.

The procedure that follows can take up to 50 seconds depending on the size and number of 3D objects on the print bed.

- ♦ Slicing of the 3D objects is carried out.
- ♦ The G-code is generated for the 3D objects.
- ♦ The G-code is saved in the specified file.
- ♦ The 3D objects are shown with their print layers in the PrusaSlicer G-code viewer.


Fig. 16


Fig. 17

You can simulate the printing process using sliders 1 and 2 (Fig. 16).

⇒ Move slider 1 up/down

- ♦ Select the layer

Move slider 2 all the way to the left and then slide it from left to right

- ♦ The movement of the tool head on this layer is shown.
- ♦ The corresponding G-code is shown on the left of the screen.

Depending on the selected layer, slider 2 moves in discrete steps rather than smoothly.

4.5 USB

To print a model, the print data (the G-code) should preferably be transferred to the printer using a USB flash drive.

Icon	Shortcut	Name	Description
	Ctrl+F	Save to USB	Saves the G-code file to a USB flash drive.

- ◆ A USB flash drive is connected to the computer.
- ◆ Slicing was carried out, so a G-code file is available in the PrusaSlicer G-code viewer.
- ◆ USB can be selected

⇒ Click the USB icon

- ◆ A "Save as" window opens.
- ◆ The target directory is the USB flash drive
- ◆ The suggested filename is the name of the 3D object file

⇒ Click "Save"; change the filename beforehand if you wish.

4.6 Options

Icon	Shortcut	Name	Description
	Ctrl+Shift+O	Options	Select the directories and language.

In the "Options" window, you can specify directories that are preset for importing and exporting 3D objects, saving projects, and saving the G-code, as well as change the language of the sliceware.

5 Expert mode

Icon	Shortcut	Name	Description
		Expert Mode	Working with PrusaSlicer. The PrusaSlicer software is started. 3D objects are transferred from the sliceware to PrusaSlicer. Further processing is then only carried out in PrusaSlicer2:

6 Print

To print a model, save the G-code file of the 3D object on a USB flash drive (see Chap. 4.5 USB) and insert the USB flash drive into the 3D printer.

Before removing the USB flash drive from the PC, select the "Safely Remove Hardware and Eject Media" function in Windows.


Fig. 18

To run the G-code file on a 3D printer, follow the instructions provided with the 3D printer.

Hochaktuell und ausführlich auf ...
Up to date and in detail at ...
Actualisé et détaillé sous ...
Aggiornato e dettagliato su ...
La máxima actualidad y detalle en ...
Актуально и подробно на ...

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Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany