making work easy

de
en

SIMPLEX slice studio

Made in Germany

20-5199 03112025

making work easy

SIMPLEX slice studio

Inhalt

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1 Introduction

SIMPLEX slice studio was developed to simplify and perfect the 3D filament printing process. To print a 3D object, the SIMPLEX 3D filament printer requires tool path information (G-codes). Use SIMPLEX slice studio to generate this G-code from 3D object files (e. g., *.stl or *.obj).

1.1 Symbols used

In these Instructions For Use (IFU) or on the equipment you will find symbols with the following meanings:

Caution

Risk of damage to the equipment if the instruction is not observed.

Note

Provides useful information to make operation easier.

► List, special attention should be paid to

- List
- Subordinate list

⇒ Instructions / required action / input / operational sequence:

You are asked to perform the specified action in the given order.

◆ Result of an action / response by the equipment / response by the program:

The equipment or program responds as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur in the text.

1.2 Current status of the Instructions For Use

This edition 03112025 of the Instructions For Use (IFU) describes the range of functions of SIMPLEX slice studio, Build: 3.1.25204.1.

Software releases with lower version numbers may not include all of the functions described here.

Software releases with higher version numbers may contain functions that are not described in this version of the IFU.

2 Safety

Before using SIMPLEX slice studio, read these Instructions For Use (IFU) carefully and keep them for future reference.

Renfert GmbH reserves the right to modify or revise this user manual at its own discretion at any time.

3 Download and installation

To download and install SIMPLEX slice studio, you must register with myRenfert. This includes registering yourself and your SIMPLEX 3D filament printer.

[Login/registration with myRenfert:](#)

4 Operation

4.1 Overview


Fig. 1

The graphic user interface is divided into three areas:

A: Navigation area

	Zoom in/out navigation area
	Display project list
	Options (Settings)

B: Print bed

Display of the print bed with the positioned objects.

C: Work area

Work steps from selecting the printer and filament to loading an object, creating the G-code, and transferring it to a 3D printer.

1.1 Starting the software

When SIMPLEX slice studio is started, a list of existing projects opens automatically.

The following project-related actions can be selected from this list:

	Delete project
	Duplicate project
	Rename project
	Open parent folder
	Open project

In addition, a new project can be created using the “Create project” button.

5 Project

5.1 Create new project – project name

When a new project is opened, the software automatically suggests a name for the project.

You can replace this suggested name with one of your own.

Then, click “Create project”.


Fig. 2

5.2 Printer setup

- Select the printer (currently, only “SIMPLEX 2 SX” is available).
- Select the filament
- Select the nozzle to be used (if different nozzles are available for the filament).

Then, click “Next”.

5.3 Layout

5.3.1 Load objects

- ◆ An empty field with a large plus sign is displayed on screen.
- ◆ The possible actions are displayed in the work area (C, Fig. 1).

If the appropriate indication has already been selected in the work area (see „5.3.13 INDICATION“), the objects are automatically placed in the optimum position for this indication when they are loaded.

This positioning also takes into account the “Simple layout” or “Sequential layout” setting (see „5.3.8 LAYOUT“).

⇒ Drag one or more objects onto the plus sign and drop them, or click on the plus sign to open a dialog box and select objects to load.

- ◆ The object is placed on the print bed.
- ◆ If several objects are loaded at the same time, they are placed on the print bed. Placement depends on the selected indication and the “Simple layout” or “Sequential layout” setting.

5.3.2 View

You can change the view of the print bed and the object as follows:

- ⇒ Rotating the mouse wheel
 - ◆ Zooming in/out
- ⇒ Right-clicking the mouse and moving it
 - ◆ Swiveling the print bed
- ⇒ Clicking on one of the predefined views
 - ◆ The print bed swivels to the desired position

These actions can be performed at any time.

The views can also be selected using the numeric keypad:

	View from
Numeric keypad 5	Front
Numeric keypad 8	Top
Numeric keypad 4	Left
Numeric keypad 6	Right
Numeric keypad 2	Bottom
Ctrl + numeric pad 5	Back

Fig. 3

5.3.3 Select

A loaded object must first be selected for the following actions.

⇒ Click "Select"

⇒ Click the desired object

- ◆ The object is shown in color.

Hold down the Shift key and click on objects to select multiple objects for the following actions.

The following actions are possible with the selected object

- Move
- Rotate
- Scale
- Delete object(s)
- Duplicate object(s).

In the following actions, a different object can be selected directly by right-clicking.

5.3.4 Move:


- ◆ The selected object is highlighted in color.

⇒ Click "Move".

- ◆ Arrows indicating movement in the X, Y, and Z directions are displayed on the screen.

⇒ Click on the desired movement arrow and move the mouse in the desired direction.

- ◆ The object is moved in the corresponding direction.

5.3.5 Rotate:


- ◆ The selected object is highlighted in color.

⇒ Click "Rotate".

- ◆ Circles indicating the three spatial directions are displayed on the screen.

⇒ Click the desired movement circle and move the mouse in the desired direction.

- ◆ The object is rotated around the selected axis.

5.3.6 Scale:


- ◆ The selected object is highlighted in color.

⇒ Click "Scale".

- ◆ Arrows are displayed in the X, Y, and Z directions.

⇒ Click the arrow pointing in the desired scaling direction, then move the mouse.

- ◆ The object is compressed or stretched in the corresponding direction.

Scaling is only possible within the range of 95 % to 105 %.

You can also enter the parameters for moving, rotating, and scaling using the keyboard in the fields displayed for the corresponding action.

Movement parameters					
X	12,49	mm	Y	31,83	mm
Z	0	mm			
Rotation parameters					
X	-22,94	°	Y	21,38	°
Z	-49,26	°			
Scaling parameters					
X	105	%	Y	95	%
Z	105	%			

Fig. 4

5.3.7 Apply / reset

Click "Apply" to set the current position of the object as the starting point for all subsequent movements.

Clicking "Select" again has the same effect as clicking "Apply".

You can undo moving, rotating, and scaling by clicking "Reset".

This will restore the state to how it was before the first action was taken.

5.3.8 LAYOUT

The "Sequential printing" settings, when active or inactive, and the "Simple layout" or "Sequential layout" buttons allow to be place objects optimally on the print bed.

The "Enable sequential printing" button is used to toggle between normal printing and sequential printing.

Normal printing / simple layout

In the normal printing setting (where sequential printing is not active), a G-code is generated that prints all objects on the print bed simultaneously, layer by layer.

- ◆ This setting includes the "Simple layout" button.


⇒ Click "Simple layout"

Only one object on the print bed:

- ◆ The object is automatically moved to the center of the print bed.

Multiple objects on the print bed:

- ◆ The objects are automatically placed in an optimal position for printing on the print bed.

Sequential printing / sequential layout

In the "Sequential printing" mode, a G-code is generated in which the objects are printed one after the other.

This requires space for the print head of the 3D printer.

The objects on the print bed are moved apart accordingly.

Sequential printing is not always possible. If there are too many objects on the print bed, some objects are moved off the print bed and excluded from the G-code.

- ◆ This setting includes the "Sequential layout" button.


⇒ Click "Sequential layout".

Only one object on the print bed:

- ◆ The object is automatically moved to the center of the print bed.

Multiple objects on the print bed:

- ◆ The objects are automatically repositioned so they can be printed sequentially without colliding with the print head.


In the "Sequential layout" mode, objects may be removed from the print bed. In this case, it is not possible to print all present objects in a single print job.

5.3.9 Add new items


⇒ Click "Add new items"

- ◆ Load additional objects using the "File – Open File" dialog box.


5.3.10 Delete object(s)

- ◆ The selected object is highlighted in color.
- ⇒ Click “Delete object(s)”.
- ◆ After confirmation, the object is deleted.


5.3.11 Duplicate object(s)

- ◆ The selected object is highlighted in color.
- ⇒ Click “Duplicate object(s)”.
- ◆ The object is duplicated and placed on the print bed.

5.3.12 Drop

- Only one object can be dropped on the print bed at a time.
- ⇒ Click “Drop”.
- ⇒ Click the area where you want to place the object.
- ◆ The object is placed on the selected area on the print bed.

i ***If the correct indication is selected under “INDICATION”, the object will be placed in the optimal position for that indication. This step is necessary for successful printing.***

5.3.13 INDICATION

- ⇒ Click the selection field.
- ◆ A selection list of possible indications will open.
- ⇒ Click the desired indication to select it.

5.4 Finalize

After all changes and settings have been made in the “2. LAYOUT” window, click on “Next” to slice the object and proceed to the “3. Finalize” window.

In the “3. Finalize” window, you can decide how to proceed with the print data of the object.

5.4.1 G-code Viewer

⇒ Click “G-code Viewer”.

- ◆ The “PrusaSlicer G-code Viewer” software will start and display the object.

Use the sliders to display the individual layers of the object and the print history in one layer.

This software is intended for viewing purposes only and cannot be used to modify the generated G-code.


5.4.2 EXPORT

Printer

⇒ Click the button for the desired 3D printer.

- ◆ The print data for the object (G-code) is sent to the 3D printer via the network.


If multiple printers have been defined in the settings, a button will be displayed for each one (see „6.3 Printer network“).

File system

⇒ Click “File system”.

- ◆ A “Save as...” dialog box opens, and the print data for the object (G-code) can be saved to a file in any directory in the file system.

USB drive

If a USB drive or SD card is detected on a USB port on the PC, a “USB drive” button will appear for each USB drive / SD card, along with its drive letter.

⇒ Click the button corresponding to the desired USB drive.

- ◆ The object's print data (G-code) is saved directly to this USB drive.

5.4.3 Printing on SIMPLEX 2 SX

Once the G-code data has been successfully exported (either sent to a printer within the network or saved to an SD card), it can be printed on a 3D printer.

Export via network:

⇒ The print data is stored in the “SIMPLEX memory”.

Export to SD card

⇒ Insert the SD card into the SD card reader of the SIMPLEX 2 SX.

⇒ The print data can be found under “External SD card”.

6 Options (Settings)

6.1 Project

Set the root directory. Internal storage location of the software for internal projects.

This setting does not need to be changed.


Fig. 5

6.2 Customization

There are different settings for the graphical user interface and the language.

6.3 Printer network

Definition of the printer(s) within the network.

IP address

To find the IP address for the SIMPLEX 2 SX, go to the printer control panel and look under:

Options >> Settings >> Network settings

SIMPLEX 2 SX

IP address or domain name

172.16.3.132

Printer name

My Printer

Fig. 6

Printer name

If necessary, you can assign different names to printers on the network to distinguish between them.

The name of the printer is displayed when selecting the printer to which the data will be transmitted (see „Printer“).

6.4 Helpdesk

A direct internet link to the myRenfert service portal.

6.5 Feedback

Use this form to report suggestions or problems with the software directly to Renfert.

6.6 Information

Information about the software version and copyright.

6.7 License

Display of license information.

6.8 Additional packages

Displays installed additional packages and provides the option to purchase and install further additional packages.

There are currently no further additional packages available.

**Hochaktuell und ausführlich auf ...
Up to date and in detail at ...
Actualisé et détaillé sous ...
Aggiornato e dettagliato su ...
La máxima actualidad y detalle en ...
Актуально и подробно на ...**

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Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany