making work easy

FAQ - Deutsch

FAQ - English

FAQ - Français

FAQ - Italiano

FAQ - Español

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FAQ - Türkçe

FAQ - Polski

FAQ - Dansk

FAQ - Magyar

FAQ - Čeština

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FAQ - English

For chapter indications, refer to the complete instruction manual on the internet.

en

Twister (No. 1826)

Error	Cause	Solution
ON / OFF switch fails to work.	<ul style="list-style-type: none"> Power not connected. Device fuse blown. 	<ul style="list-style-type: none"> Check the power supply. Replace the fuse (Sec. 6.4).
Motor fails to start.	<ul style="list-style-type: none"> Motor fault. 	<ul style="list-style-type: none"> Have the device repaired.
No, or insufficient vacuum (not in the green range on the scale).	<ul style="list-style-type: none"> Suction filter system plugged. Seal surfaces dirty. Solenoid valve fault. Vacuum pump fault. 	<ul style="list-style-type: none"> Have filter cleaned or changed (Sec. 6.3). Clean the seal surfaces (Sec. 6.2). Have the device repaired. Have the device repaired.
Vacuum release / venting, too slow.	<ul style="list-style-type: none"> Suction filter system plugged. Solenoid valve fault. 	<ul style="list-style-type: none"> Have filter cleaned or changed (Sec. 6.3). Have the device repaired.
Vacuum pump doesn't start.	<ul style="list-style-type: none"> Bowl detection pin sticking. 	<ul style="list-style-type: none"> Clean area around pin. Have the device repaired.
The vacuum pump switches on and off intermittently during the mixing process.	<ul style="list-style-type: none"> At < 100 % vacuum, the vacuum pump will switch off once the set vacuum has been reached. 	<ul style="list-style-type: none"> None, If the vacuum should drop due to seepage, the vacuum pump will automatically switch back on until the set vacuum level has been reached.

Twister venturi (No. 1827)

Error	Cause	Solution
ON / OFF switch fails to work.	<ul style="list-style-type: none"> Power not connected. Device fuse blown. 	<ul style="list-style-type: none"> Check the power supply. Replace the fuse (Sec. 6.4).
Motor fails to start.	<ul style="list-style-type: none"> Motor fault. 	<ul style="list-style-type: none"> Have the device repaired.
No, or insufficient vacuum.	<ul style="list-style-type: none"> Compressor switched off. Operating pressure too low or too high. Compressed air hose not connected, leaky or bent. Air hose too long. Square are of air hose too small. Suction filter system plugged. Seal surfaces dirty. Inlet filter plugged. Silencer plugged. Solenoid valve fault. Venturi nozzle plugged. 	<ul style="list-style-type: none"> Switch the compressor on. Check dynamic operating pressure (see technical data, Sec. 8). Switch off other consumers. Check the hose. Maximum length 2 m. Minimum inner diameter 4 mm. Have filter cleaned or changed (Sec. 6.3). Clean the seal surfaces (Sec. 6.2). Change inlet filter (Sec. 6.5). Replace the silencer (Sec. 6.6). Connect filter pressure regulator in series. Have the device repaired. Have the device repaired.
Vacuum release / venting, too slow.	<ul style="list-style-type: none"> Suction filter system plugged. Silencer plugged. Solenoid valve fault. 	<ul style="list-style-type: none"> Have filter cleaned or changed (Sec. 6.3). Replace the silencer. (Sec. 6.6). Connect filter pressure regulator in series. Have the device repaired.
No vacuum being generated.	<ul style="list-style-type: none"> Bowl detection pin sticking. 	<ul style="list-style-type: none"> Clean area around pin. Have the device repaired.
Permanent venting noise.	<ul style="list-style-type: none"> Solenoid valve fault. Bowl detection pin sticking. 	<ul style="list-style-type: none"> Have the device repaired. Clean area around pin. Have the device repaired.
Fluttering noise from the silencer.	<ul style="list-style-type: none"> Working pressure too low. 	<ul style="list-style-type: none"> Check dynamic operating pressure (see technical data, Sec. 8).

Error codes

Error code	Cause	Solution
Err. 2	<ul style="list-style-type: none">• Minimum vacuum level not yet reached.• Vacuum level has dropped below 500 mbar.	<ul style="list-style-type: none">• Ventilate bowl and remove (see section 5.5.1, points no. 2 and 3).• If the error recurs: Have the device repaired.
Err. 3	<ul style="list-style-type: none">• The vacuum is built-up too quickly due to the fact that the suction system is blocked.	<ul style="list-style-type: none">• Have filter cleaned or changed (Sec. 6.3).
Err. 4	<ul style="list-style-type: none">• Motor fault.	<ul style="list-style-type: none">• Ventilate bowl and remove (see Sec. 5.5.1, points no. 2 and 3).• Have the device repaired.
Err. 10	<ul style="list-style-type: none">• Too much material.• Mixing time too long, the material has already begun to set.	<ul style="list-style-type: none">• Ventilate bowl and remove (see section 5.5.1, points no. 2 and 3).• Only fill the bowl to the maximum fill level. This limit is valid for powders and liquids measured before mixing.• Ventilate bowl and remove (see section 5.5.1, points no. 2 and 3).• Select shorter mixing times.

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