making work easy

FAQ - Deutsch
FAQ - English
FAQ - Français
FAQ - Italiano
FAQ - Español
FAQ - Português
FAQ - Русский
FAQ - Polski
FAQ - Svenska
FAQ - Magyar
FAQ - 中文
FAQ - 日本語
FAQ - 한국어
FAQ - Čeština
FAQ - Українська
FAQ - Suomi

MT plus

FAQ - frequently asked questions

Made in Germany

21-9873 31052016








Fig. 1

FAQ - Deutsch

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FAQ - English

i For information on operating elements, see [Fig. 1](#) on page 2 of this document.

i For chapter indications, refer to the complete instruction manual on the internet.

Fault	Possible cause	Corrective action
Device fails to start.	• Power supply not established. • Mains fuse faulty. • Front cover not properly closed. • Motor overheated.	• Check the mains power supply. • Check the mains fuse. • Close the front cover (Fig. 11). • Allow the motor to cool down, then restart at the main power switch.
Trimmer disc wobbles.	• Disc is not lying flush on the flange. • Disc not tightened sufficiently.	• Reinstall the disc according to the directions. Pay attention that the disc is mounted exactly on the lock pin (Fig. 30). • Remove all deposits from the back of the disc and from the stopping face of the flange. • Hand-tighten the disc (Fig. 29).
Excessive noise while the trimmer disc is running.	• Disc not adequately tightened.	• Hand-tighten the disc (Fig. 29), or reinstall, as required (Sec. 8.1).
Trimmer table wobbling.	• Trimmer table incorrectly installed.	• Reinstall the table according to the directions (Fig. 40).

Wet Trimming:

Fault	Possible cause	Corrective action
Disc buildup.	• Water volume incorrectly adjusted at the device. • Water supply not fully open. • Spray tube plugged. • Inlet strainer plugged. • Inlet water pressure too low.	• Increase the water volume at the trimmer (Fig. 9). • Open the water supply (Fig. 6). • Clean / decalcify the spray tube (Sec. 8.3). • Remove the strainer with pliers and clean / decalcify it. • The diameter of the inlet pipe is too small resp. too many users are connected to a small pipe >>> Check if water pressure is at least 1 bar (ask your plumber). If necessary the water pressure can artificially increased by a domestic water work (electric pump and pressure tank). • Ask your plumber.
Water sprays out of the grinding opening.	• Too much water. • Water drainage impeded.	• Reduce the water volume at the trimmer (Fig. 9). • Check/properly lay out drain tube >>> slope (Fig. 3a + 3b).
No rinse water.	• Water connection not properly established. • Spray tube soiled / calcified.	• Check the water connection (Sec. 7.2 steps 4 - 6). • Clean/decalcify the spray tube (Sec. 8.3).
Water is dripping out between the front door and the housing.	• Gypsum powder order plaster grains on the sealing surface.	• Clean the sealing surface and the sealing ring.

Dry Trimming:

Fault	Possible cause	Corrective action
Disc buildup.	• Models are too wet.	• Clean the disc with a brush under running water or with plaster solvent (Renfert GO-2011 Art. No. 2011-0000).
Excessive dust being blown out of the grinding opening.	• Extraction not switched on. • Insufficient suction. • Extractor or dust bag full. • Extractor hose blocked.	• Switch the extraction on. • Increase the suction. • Replace the dust bag. • Check / clean the hose.

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