
making work easy

de
en
it

SIMPLEX 2 SX

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
Made in Germany

21-5159 15072026

making work easy

SIMPLEX 2 SX

Made in Germany

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
SIMPLEX 2 SX

Made in Germany

1	Online help and information	4
1.1	Icons and symbols used	4
2	Safety	5
2.1	Intended use	5
2.2	Meaning of these Instructions For Use (IFU)	5
2.3	Requirements for operators	5
2.4	Safe condition of the unit	5
2.5	Hazards as a result of modifications and repairs	5
2.6	Risk as a result of electric shock and fire	5
2.7	Risk of skin injuries	6
2.8	Mechanical hazards	6
3	Product description	6
3.1	General description	6
3.2	Assemblies and functional elements	7
3.3	Scope of delivery	7
3.4	Accessories	8
4	Startup	8
4.1	Unpacking	8
4.2	Setting up the unit	8
4.3	Removing the protective films	8
4.4	Removing the transport locks	8
4.4.1	Securing the print head	8
4.4.2	Transport brackets	9
4.5	Connecting the display	9
4.6	Electrical connection	9
4.7	Switching on	9
4.8	Installation Wizard	9
4.9	Registering the printer / Downloading SIMPLEX slice studio	10
4.10	Operating the display	10
4.11	Door/lid monitoring	11
5	Operation	11
5.1	Switching on	11
5.2	Filament	11
5.2.1	Position of the print head for loading/unloading the filament	11
5.2.2	Filament clamping lever	11
5.2.3	Loading the filament	12
5.2.4	Unload the filament	12
5.2.5	Filament sensor	13
5.3	File menu	13
5.4	Printing an object	13
5.4.1	Loading the print file	13
5.4.2	Clean and free print bed	14
5.4.3	Start printing	14
5.4.4	Indications during a print job	14
5.4.5	Actions during a print job	14
5.4.6	Opening the door during printing	15
5.4.7	Messages	15
5.4.8	Completing a print job	16
5.4.9	Removal of the printed object	16

6	Options	16
6.1	Move axes	17
6.1.1	Move axes manually	17
6.2	Extruder	17
6.3	Print bed leveling – meshbed measurement	18
6.4	Information	18
6.5	Service / Support	19
6.5.1	myRenfert Service Page	19
6.5.2	Print bed conditioning	19
6.5.3	Z-offset	19
6.5.4	Maintenance	19
6.5.5	Fans	19
6.5.6	Filament list – import/export	19
6.5.7	Firmware update	20
6.5.8	Log file	20
6.5.9	Installation Wizard	20
6.5.10	Factory reset	20
6.5.11	Move to positions	20
6.6	Settings	21
6.6.1	Screen lock	21
6.6.2	Lighting	21
6.6.3	Status lighting	21
6.6.4	Display	21
6.6.5	Energy saving mode	22
6.6.6	Customized information	22
6.6.7	Units	22
6.6.8	Language	22
6.6.9	Custom filament settings	22
6.6.10	Network settings	22
7	Cleaning / Maintenance	22
7.1	Maintenance intervals	22
7.2	Cleaning	23
7.3	Cleaning print bed	23
7.4	Changing the nozzle	23
7.5	Spare parts	24
8	Technical data	24
9	Warranty	25
10	Obligations regarding information	25
10.1	Disposal instructions	25
10.1.1	Disposal information for countries in the EU	25
10.1.2	Special instructions for customers in Germany	25

1 Online help and information

<https://www.renfert.com/simplex-2-sx>

The following documents and information are available for download here:

- Instructions For Use (IFU)
- Quick Start Guide
- Spare parts lists and spare parts drawings
- Repair instructions
- FAQs, cause and rectification of faults
- How-to videos
- Declarations of conformity

Note: The documents and information listed are not available for all products.

1.1 Icons and symbols used

In these Instructions For Use (IFU) or on the equipment you will find icons and symbols with the following meanings:

Danger
Immediate risk of injury. Consult accompanying documents!

Electrical voltage
Danger due to electrical voltage.

Hot surface
Danger of burns!

Caution
Risk of damage to the equipment if the instruction is not observed.

Note
Provides useful information to make operation easier.

The equipment complies with applicable EU directives.

This product complies with the relevant UK legislation.
See the UKCA declaration of conformity on the Internet at www.renfert.com.

The equipment is subject to the provisions of the EU directive 2002/96/EC (WEEE Directive).

- List, special attention should be paid to
 - List
 - Subordinate list

⇒ Instructions / required action / input / operational sequence:

You are asked to perform the specified action in the given order.

- ◆ Result of an action / response by the equipment / response by the program:

The unit or program responds as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur in the text.

2 Safety

en

2.1 *Intended use*

The unit is intended for commercial use in dental offices, dental labs or in-office labs in accordance with the IFU (Instructions For Use).

The unit is designed for manufacturing three-dimensional objects from molten filament.

These are, for example:

- Orthodontic models (models for aligner production and thermoforming, planning and diagnostic models, working models)
- Production of impression trays
- Temporary restorations

2.2 *Meaning of these Instructions For Use (IFU)*

These IFU will help you use the unit safely.

If these IFU are not observed, accidents may lead to life-threatening electric shock, injury, or material damage.

- ⇒ Only use the unit in accordance with these IFU.
- ⇒ Keep these IFU ready with the unit.
- ⇒ Provide these IFU to all subsequent users of the unit.

2.3 *Requirements for operators*

- ⇒ This product is intended for use by persons 14 years of age or older who are familiar with the procedures and regulations used in dental offices, dental labs or in-office labs.

2.4 *Safe condition of the unit*

Defective, broken or leaking components may cause injuries.

- ⇒ Check the power cord, housing, and other components such as connection cables, tubes, and the operating panel pad for damage, such as kinks, cracks, porosity, and signs of wear and tear.
- ⇒ Defective equipment must be taken out of service immediately. Disconnect the power plug. Send the unit in for repair or dispose of it safely. Until then, secure the unit to prevent it from being switched on again so that it is not used inadvertently.
- ⇒ Do not expose the unit to mechanical impact. Do not drop it.

2.5 *Hazards as a result of modifications and repairs*

Improper modifications and repairs to the unit may cause accidents, fire, and electric shock.

- ⇒ Only have the unit repaired by qualified specialist dealers or return it.

The use of unsuitable accessories or spare parts may result in unpredictable accidents and damage.

- ⇒ Only use accessories and spare parts supplied or approved by Renfert GmbH.

2.6 *Risk as a result of electric shock and fire*

Risk of electric shock and of fire as a result of excessive operating voltage or if there is no protective conductor.

- ⇒ Ensure that the power plug is easily accessible.
- ⇒ Only operate the unit using a power cord equipped with a plug appropriate to the local power supply. Only have the power cords changed by a qualified electrician.
- ⇒ Only operate the unit if the information on the type plate matches the power supply.
- ⇒ Only connect the unit to electrical sockets that are connected to the protective conductor system.

2.7 *Risk of skin injuries*

Danger of burns!

Print table, print bed and nozzle get hot during operation.

- ⇒ Do not touch the print table, print bed, or nozzle while the printer is heating up or during printing.
- ⇒ Wait for the print bed, print table, and nozzle to cool down before working on them.
- ⇒ Do not touch these components during printing!
- ⇒ After printing, allow the components to cool before removing them!

2.8 *Mechanical hazards*

Risk of crushing!

- ⇒ Only operate the unit with the door and lid closed.
- ⇒ Movement commands are still executed until they are completed, even after the door or cover has been opened.
- ⇒ Do not open the door or lid while the print head or console is in motion.
- ⇒ During operation, never reach under the console.

In the event of a crushing hazard:

- ▶ Switch the unit off.
- ▶ Push the crushing component aside by hand and turn the spindle to move the console.

3 **Product description**

3.1 **General description**

The SIMPLEX 2 SX 3D filament printer is used to create three-dimensional objects from molten filament.

The SIMPLEX 3D filament printing system was developed by Renfert to simplify and perfect the 3D filament printing process. It offers an excellent manufacturing experience without prior programming or 3D printing experience.

To print 3D objects you need:

- the SIMPLEX 2 SX 3D filament printer
- The SIMPLEX slice studio software (can be downloaded from the Internet after registering the unit at myRenfert).
- a SIMPLEX filament of your choice.

3.2 Assemblies and functional elements

- 1 Display / touchscreen
- 2 ON/OFF switch / connection to mains power supply
- 3 Filament sensor
- 4 SD card slot
- 5 Print table (heated) with print bed (removable metal plate)
- 6 Print head (nozzle, extruder, and print bed sensor)
- 7 Adjustment dials for print bed leveling
- 8 Bowden (filament guide tube)
- 9 Console


Fig. 1

3.3 Scope of delivery

- 1 3D filament printer SIMPLEX 2 SX
- 1 Power cord
- 1 Quick Start Guide
- 1 SIMPLEX slice studio
(can be downloaded from the Internet after registering the machine at myRenfert)
- 1 Offset sensor
- 1 Micro SD card
- 1 USB adapter for using the micro SD card on a PC
- 1 Brass brush
- 1 Wire cutter
- 1 TORX20 screwdriver
- 1 TORX10 screwdriver
- 1 SW2 Allen key

3.4 Accessories

1735 0110	SIMPLEX study model 2 – 1000 g
1735 0210	SIMPLEX working model 2 – 1000 g
1735 0320	SIMPLEX aligner model 2 – 1000 g
1735 0610	SIMPLEX multi-use model 2 – 1000 g
1734 0060	Tool set
1734 0061	Spare parts package 0.4 mm nozzle
1734 0062	Spare parts package 0.25 mm nozzle

4 Startup

You can also find a video on the Renfert website that shows how to install and start-up the unit.

Before startup, the unit must be stored at room temperature for at least 3 hours.

Do not apply force to the unit during unpacking, setup, operation, or maintenance.

4.1 Unpacking

Pay attention to the weight of the printer (see technical data) and ask a second person to help if necessary.

⇒ Remove the unit and accessories from the packaging.

⇒ Check the delivery for completeness (compare with the scope of delivery).

4.2 Setting up the unit

Select a setup location where the printer is not exposed to direct sunlight or drafts (e.g., from air conditioning).

⇒ Set up the unit on a firm, stable, level workbench or worktop.

⇒ Position the unit so that the display does not protrude beyond the work surface.

4.3 Removing the protective films

Remove the protective films from the trim parts, the door, the lid, and the display.

4.4 Removing the transport locks

4.4.1 Securing the print head

Two cable ties secure the crossbar and print head for transport.

⇒ Use the provided side cutter to cut and remove the cable ties.


Fig. 2

4.4.2 Transport brackets

For transport, the console is attached to the housing with two brackets.

- ⇒ Loosen the screws with the provided TORX20 screwdriver.
- ⇒ Remove the brackets and screws and store them in a safe place.

Do not dispose of the brackets and screws! Keep them for later use!

If the printer is to be transported and the transport locks need to be reattached, the console and print head will automatically move into the necessary positions.

Refer to chapter „6.5.11 Move to positions“


Fig. 3

4.5 Connecting the display

- ⇒ Remove the adhesive tape that was used to secure the connection cable during transport.
- ⇒ Remove the display from the accessory packaging.

Do not touch the connection pins on the display! The display can be damaged by electrostatic discharge.

- ⇒ Insert the green connector into the port on the back of the display until it clicks into place.
- ⇒ Insert the display vertically into the support from above.
- ⇒ Carefully swing the display backward until it clicks into place.


Fig. 4

4.6 Electrical connection

Observe the instructions in chapter 2.6 “Risk as a result of electric shock and fire!”

Make sure that the ON/OFF switch (2, Fig. 1) is set to “Off” (O).

- ⇒ Plug the power cord into the connection for the mains power supply (2, Fig. 1) on the rear of the unit.
- ⇒ Connect the power plug to an electrical socket of the building.

4.7 Switching on

- ⇒ Switch the unit on using the ON/OFF switch (2, Fig. 1) on the rear of the unit.
 - ◆ The build chamber fan starts up briefly.
 - ◆ The Renfert logo appears on the screen, and the display begins to initialize.
 - ◆ The measuring tip of the print bed sensor is extended twice briefly. After that, the print bed sensor lights up red.

If the print bed sensor flashes red, the print table was positioned too high when the printer was switched on. In this case, no action is required.

The error status will be deleted during the next reference run.

4.8 Installation Wizard

The Installation Wizard appears on the display after the unit is switched on for the first time.

This wizard will guide you through the initial steps of starting up the unit.

The following actions are performed:



- ◆ Setting the language
- ◆ Setting the units
- ◆ Print bed leveling (see also chapter „6.3 Print bed leveling – meshbed measurement“)
- ◆ Inserting a filament (see also chapter „5.2.3 Loading the filament“)
- ◆ Printing a test file

4.9 Registering the printer / Downloading SIMPLEX slice studio

The SIMPLEX slice studio software is required to create print files.

You can download this software from the Renfert website.

To do this, you must first register the printer.

Unit registration and software download via myRenfert	SIMPLEX Tutorials
 myRenfert Login	 Video

4.10 Operating the display

The display is a touchscreen.

To select a menu item, make a selection from a list. Alternatively, you can perform an action by tapping the corresponding, visually highlighted buttons or by tapping a list entry.

For long lists, you can scroll forward and backward by swiping up and down.

The main menu with selection options is located on the left side of the display:

 Home screen	 File menu	 Options
--	--	--

The following button icons are used for operation and control

	OK / Confirm / Continue / Next
	Cancel / Cancel print <i>Even if the "Cancel" button has been pressed, actions that have already been started will still be completed.</i>
	Back
	Decrease/increase value
	Start printing
	Pause printing
	Resume printing
	Deleting the print file (G-code)

4.11 Door/lid monitoring

The door and lid must be closed during any movement of the print head or print table.

Before the 3D printer performs any movement, it checks that the door and lid are closed.

If this is not the case, you will be prompted to close the door/lid. You will then need to restart the action.

Keep the door and lid closed during all movements until the movement is complete.

If the door or lid is opened during printing, the print table moves down slightly and the print head moves to a pause position. Printing can only be resumed with the door/lid closed.

5 Operation

5.1 Switching on

The unit is switched on and off using the ON/OFF switch (2, Fig. 1) on the back of the device.

⇒ Switch on the unit.

- ◆ The build chamber fan starts up briefly.
- ◆ The Renfert logo appears on the screen, and the display begins to initialize.
- ◆ The measuring tip of the print bed sensor is extended twice briefly. After that, the print bed sensor lights up red.
- ◆ The display will return to the “Home” screen.

5.2 Filament

Danger of burns!

If you change the filament during or shortly after printing, the print bed may still be hot. Do not touch the print table and/or print bed!

The 3D printer does not contain a device for detecting the type of filament loaded.

It is the operator's responsibility to ensure that the filament intended for printing is loaded. Only a plausibility check is performed to verify that the loaded filament matches the one defined in the print job.

If this information is different, a warning message is indicated.

5.2.1 Position of the print head for loading/unloading the filament

When loading and unloading the filament, the print head is preferably positioned on the left side so that it does not interfere with the loading/unloading of the filament due to tight radii in the filament guide tube.

The print head automatically moves to this position at the end of a print job (regular end and abort) and after the print bed has been levelled.

If, for some reason, the print head is too far to the right and a reference run has not yet been performed after switching on the printer (see chapter 5.4.3 “Moving axes / Reference run”), the print head can be manually moved slightly to the side so that it does not interfere with filament loading/unloading.

5.2.2 Filament clamping lever

The filament clamping lever adjusts the extruder's filament feed clamping mechanism.

This lever should always be in its middle position and does not need adjusting.

The correct position is essential for good printing results! Different lever positions may cause printing problems and interruptions in the print flow.

In the left position, the clamping mechanism is open and the filament is not fed.

The right position is only intended for extremely hard filaments.



Fig. 5

5.2.3 Loading the filament

- ⇒ Open the door.
 - ⇒ Check that the two support rollers for the filament roll are properly positioned in the guides.
 - ⇒ Place the filament roll on the support rollers so that the filament on the top of the roll is facing into the unit.
 - ⇒ Insert the filament into the filament sensor (3) and push it in until you feel resistance when it reaches the print head.
- Do not push further with force!

- ⇒ Tap the "Filament" button on the "Home" screen (Fig. 6).
 - ◆ A list of stored filaments will appear.
- ⇒ Select the filament you want to load by tapping it.
 - ◆ The nozzle will heat up to the temperature assigned to the filament.
 - ◆ When the correct temperature is reached, the extruder will start and an animation will show how to insert the filament.
- ⇒ Push the filament in with light pressure so that it can be grasped by the extruder.
- ⇒ Once the extruder grabs the filament, stop pushing and let the extruder feed it.
- ⇒ Extrude a small quantity of filament.


Fig. 7

If the filament is being conveyed smoothly:

- ⇒ Tap the "Okay" button.
 - ◆ The window will close and the "Home" screen will reappear.
 - ◆ The extruder will stop.
 - ◆ The nozzle heater will be turned off again.


Fig. 8

If the "Confirm" button is not pressed, the extruder stops automatically after 1 minute.

5.2.4 Unload the filament

- ⇒ Open the door.
- ⇒ Tap the "Filament" button on the "Home" screen (Fig. 8).
- ⇒ Tap the loaded filament in the selection.
- ⇒ Confirm unloading in the following window.
 - ◆ The nozzle is heated to the required unloading temperature.
 - ◆ Once the temperature is reached, an animation will appear on the display to guide you through the next steps.
 - ◆ First, a small amount of filament is extruded.
 - ◆ The extruder then pushes the filament back approximately 30 mm.
- ⇒ Gently pull the filament out of the filament sensor.
- ⇒ Confirm unloading.


Fig. 9

If the filament cannot be pulled out, you can tap the "Repeat inloading" button to make the extruder move back another 30 mm.

If there is no filament piece protruding from the sensor that can be pulled out, proceed as follows:

- ⇒ Press the clamping ring of the Bowden clamp (A) downwards on the print head.
- ⇒ Pull the Bowden (8) out upwards (B).
- ⇒ First, remove the filament residue from the print head. Then, remove the filament from the Bowden tube.
- ⇒ Reinsert the Bowden into the clamping mechanism.


Fig. 10

5.2.5 Filament sensor

The filament sensor monitors the filament flow during printing.

If the filament flow is too low—for example, if the nozzle is clogged or the filament has run out—the filament sensor will generate a signal that pauses printing.

At the same time, the following message is displayed on the screen:

Filament error on extruder 0: tooLittleMovement

5.3 File menu

In the "File" menu, you will see the print files stored on the printer. These can be selected for printing or deleted.

The print data is located in the printer at the following locations:

- SIMPLEX memory
 - Sample files
 - Print data that has been transmitted to the printer via SIMPLEX slice studio.
- External SD card
 - Print data created on a PC and transferred to the printer via an SD card.

Fig. 11

5.4 Printing an object

The print file for the object to be printed is a so-called G-code file.

You created this print file using the SIMPLEX slice studio software.

The print file can be transmitted to the printer via an SD card or loaded into the printer using SIMPLEX slice studio.

5.4.1 Loading the print file

- ⇒ Remove the SD card from the USB adapter on the PC.
- ⇒ Insert the SD card carefully into the SD card slot (4) on the 3D printer until it clicks into place.
- ⇒ Select the "File" menu.
- ⇒ Select "External SD card".

- ◆ The maps and print files on the SD card will be displayed.

Alternatively, if the print file was transmitted to the printer via SIMPLEX slice studio:

- ⇒ Select "SIMPLEX memory."

- ◆ The unit displays the print files transferred with SIMPLEX slice studio and stored in the printer's internal memory.

- ⇒ Navigate to the desired print file and
- ⇒ select it by tapping on it.


Fig. 12

- ♦ The display shows the following information:
 - Name of the print file.
 - The (estimated) printing time.
 - The filament defined in the print file.
 - The estimated material consumption in g [oz.] (estimated value).
 - The currently loaded filament.

If the loaded filament matches the filament defined in the print file, both are displayed in white.

If they do not match, both are shown in yellow, and a warning triangle indicates the mismatch.

⇒ If a wrong filament is loaded, cancel the operation by tapping the "Back" button.

- ♦ The display will return to the "Home" screen.

⇒ Load the required filament (see chapters „5.2.4 Unload the filament“ and „5.2.3 Loading the filament“) and reload the print file.

You can also delete the print file by tapping the Trash Can icon.

CAUTION: Tapping on the Trash Can icon immediately deletes the file without further confirmation.


Fig. 13

5.4.2 Clean and free print bed

Before printing, make sure that the print bed is clean, especially free of grease.

See chapter „7.3 Cleaning print bed“ for instructions on cleaning the print bed.

Make sure that there are no foreign objects under the print table to avoid interfering with the printing process and damaging the print table.

5.4.3 Start printing

The desired print data has been loaded and the correct filament has been inserted.

⇒ Close the lid and door.

⇒ Start printing by tapping the "Start" button.

⇒ Confirm that the print bed is clear.

- ♦ The print bed will heat to the required temperature.
- ♦ Once the print bed has reached the desired temperature, the print head and print table perform a reference run (see chapter 6.1, "Moving Axes / Reference Run") and the nozzle heats up to the required temperature.
- ♦ Once the nozzle has reached the target temperature, the printing process begins automatically.

5.4.4 Indications during a print job

During printing, the display indicates the following information:

- The remaining printing time (hh:mm).
- How much of the print job has been printed (numerical value in % and graphical animation).
- The current nozzle temperature (target and actual value).
- The current print bed temperature (target and actual value).
- The filament defined in the print data.
- The file name of the print data.

The remaining printing time is estimated by the controller based on the print commands that have been read in so far. This estimate is constantly updated and may cause "jumps" in the displayed progression.


Fig. 14

5.4.5 Actions during a print job

The following actions can be performed during printing:

Pause printing:

⇒ Tap the "Pause" button

- ♦ The print job is interrupted.
- ♦ The print head moves to a pause position.


Fig. 15

Resume printing

- ⇒ Tap the "Play" button.
- ◆ The print job is resumed.


Fig. 16

Cancel print

**To cancel a print job, it must first be paused.
A canceled print job cannot be resumed!**

- ⇒ Tap the "Pause" button
 - ◆ The print job is interrupted.
 - ◆ The print head moves to a pause position.
- ⇒ Tap the "Cancel" button
 - ◆ Confirmation prompt asking if you really want to cancel the print.
- ⇒ Confirm to cancel the print.
 - ◆ The print job is canceled.
 - ◆ The display returns to the "Home" screen.


Fig. 17


Fig. 18

Changing the nozzle or print bed temperature:

- ⇒ Tap the corresponding button.
- ⇒ Set the temperature using the adjustment dial.
- ⇒ Confirm the setting.

A small orange triangle in the temperature display window indicates a change in temperature.

The triangle will remain displayed even if the temperature is changed back to its original setting.


Fig. 19

Changing the nozzle or print bed temperature may negatively affect the print result!

Changing the filament

If the filament needs to be changed during printing as a precautionary measure, this can be done during a pause in printing.

- ⇒ Tap the "Pause" button.
 - ◆ The print job is interrupted.
 - ◆ The print head moves to a pause position.
- ⇒ Tap the "Filament" button.
 - ◆ The display toggles to the "Unload Filament" menu.
- ⇒ Unload the inserted filament and load the new filament.
 - Refer to chapters „5.2.4 Unload the filament“ and „5.2.3 Loading the filament“.
- After loading the filament, the print process display will reappear, allowing you to continue printing.
- ⇒ Tap the "Play" button.
 - ◆ The print job is resumed.

Fig. 20

If a loaded filament does not match the one defined in the print data, the name of the filament is displayed in yellow.

5.4.6 Opening the door during printing

Opening the door acts as a pause command.

- ⇒ The door is opened.
 - ◆ Message: **Print paused by external trigger**
 - ◆ The print job is interrupted.
 - ◆ The print head moves to a pause position.

To continue the print job, close the door and press the "Play" button.

5.4.7 Messages

The printer's status is continuously monitored by the printer's control system.

Various events trigger specific responses. These events are indicated on the display.

Messages can be closed by tapping the "Okay" button.

Event	Message	Response
• Filament flow too low or filament empty	• Filament error on extruder 0: tooLittleMovement	• Print job paused
• Door/lid opened during printing	• Print paused by external trigger	• Print job paused

5.4.8 Completing a print job

When the print job is finished properly, the print bed moves to its lowest position to allow easy removal of the print object.

The display shows:

- The duration of the print.
- The amount of filament used.

⇒ Confirm the information.

- ◆ Displaying the "Home" screen

5.4.9 Removal of the printed object

Danger of burns!

The nozzle and print bed are still hot after printing!

Allow the print bed and nozzle to cool down before removing the printed object.

⇒ Observe the temperature of the nozzle and print bed on the "Home" screen.

⇒ Wait until both components have cooled sufficiently (below 40 °C).

⇒ Open the door.

⇒ Remove the print bed (removable metal plate).

⇒ Remove the printed object by gently bending the print bed.

⇒ Reinstall the print bed, ensuring that it is in the correct position.

The recess at the rear edge must be against the two screws (Fig. 18).

Be sure to remove the print bed to remove the printed object. Otherwise, applying force during removal may damage the guides, which will negatively impact print quality over time.

If printed objects stick too strongly or not strongly enough to the print bed, it should be reconditioned. Refer to chapter „6.5.2 Print bed conditioning“ for more information.


Fig. 21

6 Options

⇒ Tap the "Options" button.

- ◆ The "Options" menu is displayed.

In this menu, various functions can be started directly, or additional functions can be selected from lists.


Fig. 22

6.1 Move axes

In this menu, the axes and the print table can be moved manually. This is only possible after a reference run has been completed. Until then, the movement buttons are grayed out and inactive.


Fig. 23


Fig. 24

EXCEPTION:

The “Print table all the way down” button is always active and can be used to move the print table down to its lowest position.

This may be necessary, for example, if printing has been canceled and the nozzle is still on the model.

The movement of the print table is limited by a limit switch.

6.1.1 Move axes manually

During the reference run, make sure that there are no print objects or other objects on the print bed to avoid collisions and damage to the print head.

The reference run, as well as the movement of the axes and the print table, are only possible when the door is closed.

In the "Axes" menu:

⇒ Tap the "Reference Run" button.

⇒ If the print bed is free, tap the "Okay" button.

◆ The 3D printer will perform a reference run.

◆ After that, the arrow keys can be used to move print head and print table.

⇒ Use the arrow keys to move the print head and print table to the best position for the desired action.

6.2 Extruder

This menu allows you to extrude or retract the filament manually.

If a filament is loaded, the nozzle is heated to the filament temperature before this display is indicated.

The filament can be extruded or retracted using the up/down arrow buttons.

Tapping on the temperature will bring up an adjustment dial that can be used to change the temperature of the nozzle.

If a filament is loaded or unloaded via this menu item, the display will not be updated. This means that the display of loaded filaments in the “Home” screen will not change.

To prevent this, you can jump to the “Filament loading/unloading” function by pressing the respective button (see chapters „5.2.4 Unload the filament“ and „5.2.3 Loading the filament“).


Fig. 25


Fig. 26

6.3 Print bed leveling – meshbed measurement

Print bed leveling is important because it ensures that the first layer of filament adheres evenly to the bed. This ensures a stable start to printing, prevents warping, and improves overall print quality.

After installation, check the print bed leveling and adjust if necessary.

The display will guide you through the steps necessary to level the print bed.

- ⇒ Switch on the printer.
 - ⇒ Make sure the print bed is clear.
 - ⇒ Close the door and lid.
 - ⇒ Tap the “Options Menu” button (Fig. „22“) and select “Print Bed Leveling.”
 - ⇒ Confirm that you want to start print bed leveling.
 - ◆ The printer will perform a reference run.
 - ◆ The print head will move to the center rear of the print bed and measures the height.
 - ◆ The print head moves to the front left and measures the height.
 - ◆ The display shows how far the height at this point is from the target value and the direction to turn the adjustment dial to adjust the height.
 - ⇒ Open the door.
 - ⇒ Turn the adjustment dial in the indicated direction until the display turns green and the "Okay" button is enabled.
 - ⇒ Close the door.
 - ⇒ Tap the "Okay" button.
 - ◆ Next, the print head moves to the front right and measures the height.
 - ◆ The display shows how far the height at this point is from the target value and the direction to turn the adjustment dial to adjust the height.
 - ⇒ Open the door.
 - ⇒ Turn the adjustment dial in the indicated direction until the display turns green and the "Okay" button is enabled.
 - ⇒ Close the door.
 - ⇒ Tap the "Okay" button.
 - ⇒ The printer performs a meshbed measurement. This takes approximately 5 minutes.
 - ⇒ After the measurement is complete, the print head moves to a parking position.
- The meshbed measurement allows the printer to adjust the nozzle height during printing to compensate for unevenness in the print bed.


Fig. 27

The Meshbed measurement determines and stores the height of 100 dots. These measurements are used when printing objects.

6.4 Information

Display of printer information:

- Name of printer model
- Software versions
- Operating hours
- Printing hours
- Printing time, bed (can be reset after changing the print bed)
- Printing time, nozzle (can be reset after changing the nozzle)


Fig. 28

6.5 Service / Support

This menu contains the submenus described below.

6.5.1 myRenfert Service Page

⇒ Options – Service/Support – myRenfert Service Page

Indication of various QR codes and internet links that lead to various myRenfert pages:

- Service videos
- Device Registration
- myRenfert Service Portal
- Service requests


Fig. 29

6.5.2 Print bed conditioning

⇒ Options – Service/Support – Print bed conditioning

Print bed conditioning can improve the adhesion of the models to the print bed.

To do this, first place a damp cellulose cloth (e.g., paper towel) on the print bed. Then, heat the print bed. Finally, wait until the cellulose cloth has dried.

For better accessibility, move the console to its lowest position. The door must be closed for this.

The display will guide you through the steps necessary to condition the print bed.

For safety reasons, the print bed heater automatically turns off after 10 minutes.

6.5.3 Z-offset

⇒ Options – Service / Support – Z-offset

The Z-offset is a value specific to each unit that describes the difference between the switching point of the print bed sensor and the nozzle tip.

This value must be redetermined after modifications to the print head, e.g. a nozzle change.

The Z-offset must be determined and adjusted to ensure optimal adhesion of the first filament layer to the print bed. An incorrectly set Z-offset can result in adhesion issues, as well as damage to the print bed and print head.

To correctly determine the Z-offset, the filament must be removed prior to the measurement and the underside of the nozzle must be cleaned.

The Z-offset probe is required for the determination process (see scope of delivery).

The display will guide you through the process of determining the Z-offset step by step.

Before the measurement starts, the nozzle is heated to 200 °C. This results in an approximate wait-time of 60 seconds before the actual measurement is performed.

Upon completion of the measurement, the result is checked to see if it falls within a plausible range. If this is not the case, the Z-offset is set to zero.

Check the displayed Z-offset. If the value is zero, exit the menu item and repeat the measurement!

Reinsert the short-circuit bridge after the measurement! Otherwise, this menu item cannot be exited.


Fig. 30

6.5.4 Maintenance

⇒ Options – Service/Support – Maintenance

Displays the operating hours until various maintenance tasks are required, such as checking the belt tension or lubricating the Z-axis spindle.

The hours displayed in orange show how long the corresponding maintenance task has been overdue.

After performing the maintenance task, reset the hour counter by tapping the “Okay” button.

Refer to chapter „7.1 Maintenance intervals“.

A factory reset does not change the hour counters for maintenance work.

6.5.5 Fans

⇒ Options – Service / Support – Fans

The individual fans can be switched on and off for test purposes and their speed can be adjusted.

6.5.6 Filament list – import/export

⇒ Options – Service / Support – Filament list import / export

The default filament data stored in the display can be exported to a USB stick and new default filament data can be imported.

Insert the USB stick into the USB port on the display (see Fig. 28).

Import is only possible when no filament is loaded.

6.5.7 Firmware update

⇒ Options – Service / Support – Firmware update

The firmware update installs new software for the display.

The update is available from the Renfert website or directly from a USB stick.

If you obtain the update from the Renfert website, you must first copy the file to a USB stick.

Proceed as follows to perform the update:

⇒ Remove the cover from the USB port.

⇒ Insert the USB stick into the USB port.

⇒ Tap on Options – Service/Support – Firmware update.

⇒ Follow the instructions on the display.

- ◆ The update will be carried out.

- ◆ At the end of the update, the “Home” screen will appear.

⇒ Remove the USB stick.

⇒ Replace the cover on the USB port.

Fig. 31

DO NOT remove the USB stick during the update and DO NOT switch off the unit!

Replace the cover immediately after the update. This will prevent it from being lost and will prevent water or detergent from entering and damaging the USB port during cleaning.

This update only affects the software displayed on the screen.

The printer firmware is updated via Renfert Slice studio. The procedure is described there.

6.5.8 Log file

⇒ Options – Service / Support – Log file

The display records various actions in a log file, which can help identify and analyze errors more quickly.

This allows problems to be solved more efficiently, thereby minimizing downtime.

The list includes the last 50 actions.

The list can be exported to a USB stick so that it can be sent to Renfert Service if necessary.

⇒ Remove the cover from the USB port (Fig. 28).

⇒ Insert the USB stick into the USB port.

⇒ Tap the “Export” button.

- ◆ The log file is copied to the USB stick.

- ◆ The name of the log file is “logs.txt”.

Once the export is complete:

⇒ Remove the USB stick.

⇒ Replace the cover on the USB port.

6.5.9 Installation Wizard

⇒ Options – Service / Support – Installation Wizard

You can run the Installation Wizard again.

See chapter „4.8 Installation Wizard“.

6.5.10 Factory reset

⇒ Options – Service / Support – Factory reset

All settings are reset to the factory settings.

The language resets to “English” and the Installation Wizard opens again when the printer is turned on.

6.5.11 Move to positions

⇒ Options – Service / Support – Move to positions

The print head and print table can be moved to two defined positions.

Transport position

If you need to transport the printer, move the console to this position so that you can reattach the transport brackets and secure the print head with cable ties.

See chapter „4.4 Removing the transport locks“.

6.6 Settings

This menu contains the submenus described below.


Fig. 32

6.6.1 Screen lock

⇒ Options – Settings – Screen lock

It is possible to lock the operation of the display.

You can activate and deactivate this option to lock the screen.

The screen lock can be unlocked with a 4-digit code. The default code is “0000”.

Display:

	The screen can be locked, but is not currently locked
	The screen is locked.

Locking the screen:

⇒ Press the white lock for 5 seconds.

- ◆ A yellow rectangle appears, which disappears after 5 seconds
- ◆ The lock turns yellow and the screen is locked.

Unlocking the screen:

⇒ Press the yellow lock for 5 seconds

- ◆ A yellow rectangle appears and disappears after 5 seconds.
- ◆ A window for entering the unlock code appears when this option is activated.
- ◆ The lock turns white and the screen is unlocked again.

Turning the printer off and on again unlocks the screen.

6.6.2 Lighting

⇒ Options – Settings – Lighting

The interior lighting can be turned on and off.

Use the slider or the arrow keys to set the desired brightness.

6.6.3 Status lighting

⇒ Options – Settings – Status lighting

You can set the color of the status lighting individually for the following modes:

- Stand-by
- The printer is active

⇒ Select status lighting

- ⇒ Tap the color field for “Standby” or “Printer active.”
 - ◆ A color palette displaying the available colors is shown.
- ⇒ Tap the desired color on the color palette.
- ⇒ Confirm your selection with the "Okay" button.

6.6.4 Display

- ⇒ Options – Settings – Display
- The brightness of the display can be set within a range of 10 % – 100 %.
- Use the slider or the arrow keys to set the desired brightness.

6.6.5 Energy saving mode

- ⇒ Options – Settings – Energy saving mode
- The display and interior lighting will turn off after 1 to 60 minutes if the display is not touched during that period.
- A setting of 0 min. switches off the energy-saving mode, i. e. the display and/or interior lighting remain permanently on.
- The energy-saving mode is terminated by touching the display.

6.6.6 Customized information

- ⇒ Options – Settings – Personalized information
- The printer can be given a name, which is shown on the display.

6.6.7 Units

- ⇒ Options – Settings – Units
- Toggles the units for:
 - Temperature: °C / °F
 - Weight: g / oz

6.6.8 Language

- ⇒ Options – Settings – Language
- Toggles the display language.

6.6.9 Custom filament settings

- ⇒ Options – Settings – Custom Filament Settings
- Two custom filaments can be activated.
- The following can be defined for each filament:
 - Name
 - Loading/unloading temperature
 - Printing temperatures
 - Print bed temperature.
- Activated filaments are displayed in the filament list for loading and unloading.

6.6.10 Network settings

- ⇒ Options – Settings – Network settings
- The printer name used on a network and the current IP address are displayed.
- Only network specialists should make changes or modify any settings.**

7 Cleaning / Maintenance

The unit does not contain any internal parts that require maintenance. Opening the unit is not permitted other than as described below!

Before cleaning and maintenance work:

- Switch off the unit and disconnect the power plug!

7.1 Maintenance intervals

The software includes operating hour counters for the following maintenance tasks.

- Check belts
- Perform print bed leveling
- Check the fan function, especially the object fan and hot-end fan
- Lubricate Z-spindle
- Clean the X and Y rails (dry cellulose cloth, e.g., kitchen paper; do not use cleaning agents)
- Check the print bed/nozzle, replace if necessary

Once the specified operating hours have been reached, a message appears on the display indicating that maintenance is required.

Refer to chapter „6.5.4 Maintenance“.

Regardless of this, these maintenance tasks should be performed at least once a year.

Helpful service videos can be found at [myRenfert Website](#), which can also be accessed via the QR code below.

7.2 Cleaning

Do not steam clean the unit under any circumstances.

To clean the unit, only wipe the outside with a damp cloth.

Do not use any solvent-based or abrasive cleaning agents.

7.3 Cleaning print bed

⇒ Keep the print bed free of grease (fingers, oil, etc.).

Regular cleaning with unsuitable cleaning agents can damage the print bed and render it unusable!

Do NOT use the following cleaning agents:

- ▶ Solvent-based cleaning agents
- ▶ Alcohol
- ▶ Isopropanol or methylated spirit
- ▶ Denatured alcohol
- ▶ Soap
- ▶ Disinfectant

For regular cleaning:

⇒ Remove the print bed.

⇒ Rinse the print bed under running water.

⇒ Wipe thoroughly with a cellulose cloth (e.g., kitchen paper) and rub dry.

In exceptional cases, cleaning solvent can be used to remove heavy grease stains. Then rinse thoroughly with water.

Cleaning solvent is not suitable for regular cleaning because it damages the print bed surface over time!

7.4 Changing the nozzle

**The filament must be unloaded before changing the nozzle.
See chapter „5.2.4 Unload the filament“.**

Caution – risk of burns!

The nozzle may still be hot from unloading the filament.

Allow the nozzle to cool down!

To change the nozzle:

⇒ Move the print table to its lowest position, see chapter „6.1 Move axes“ – Exception.

⇒ Switch off the unit so that the print head can be moved freely.

⇒ Loosen and unscrew the old nozzle (20).

**Only turn the nozzle (20).
Do not turn the heater (21)!**

⇒ Screw the new nozzle in as far as possible, then tighten it by hand.

If the nozzle is not tightened properly, it may result in filament leaks, blockages, and uneven printing.

Therefore, it is essential to screw the nozzle in completely. This avoids mechanical issues, ensures tightness, guarantees optimum heat transfer and maintains the mechanical stability of the nozzle.

After replacing the nozzle, or after removing and reinserting the installed nozzle, the difference between the print bed sensor's switching point and the nozzle tip must be recalculated.

This requires performing the Z-offset. See Chapter „6.5.3 Z-offset“.


Fig. 33

7.5 Spare parts

You can find wear parts and spare parts in the spare parts list available online at, see chap. „1 Online help and information“.

The components excluded from the warranty (wear parts, consumables) are marked in the spare parts list.

The serial number, date of manufacturing and equipment version are shown on the nameplate of the equipment.

8 Technical data

Item number:	1734 0000
Nominal voltage / V:	230
Permissible mains voltage / V:	100 - 240
Mains frequency / Hz:	50/60
Power consumption / W *1), max.:	350
Mains input fuse:	T5A
Sound pressure level LpA *2) / dB(A):	41
Dimensions (width x height x depth) / mm [inch]:	523 x 450 x 577 [20.6 x 17.7 x 22.7]
Weight / kg [lbs.]:	30 [66.1]
This product contains a light source of energy efficiency class:	G
Ambient conditions for operation *3)	
Ambient temperature / °C [°F]:	10 - 30 [50 - 86]
Ambient conditions for storage and transport	
Ambient temperature:	- 20 to + 60 °C [- 4 to + 140 °F]
Maximum relative humidity:	80 %
Ambient conditions according to EN 61010-1 *3)	
Ambient conditions:	Indoors
Installation height, max.:	2,000 m above sea level
Ambient temperature:	5 - 40 °C [41 - 104 °F]
Relative humidity:	Max. 80 % at 31 °C [87.8 °F], decreasing linearly to 50 % at 40 °C [104 °F]
Voltage fluctuations in the mains power supply from the nominal value, max.:	10 %
Contamination level:	2
Over-voltage category:	II

- *1) Power consumption at nominal voltage
- *2) Noise level in accordance with EN ISO 11202
- *3) The product performs as specified in the technical data when operated within ambient conditions. If the unit is operated outside these ambient conditions, safety devices contained in the product may stop operation, the optimum performance specified in the technical data may not be achieved, and the service life of components may be reduced. The product fulfills the safety requirements of standard EN 61010-1 when operated within ambient conditions. Operation outside of these ambient conditions is not permitted, as electrical safety may be compromised.

9 Warranty

Renfert offers a **3-year warranty** on all parts of the equipment, provided that the equipment is properly used.

Warranty claims may only be made upon presentation of the original sales receipt from the authorized dealer.

Parts that are subject to natural wear and tear (wear parts) and consumables are excluded from the warranty. These parts are marked in the spare parts list.

The warranty is voided in the case of improper use; failure to observe the operating, cleaning, maintenance, and connection instructions; in the case of independent repairs or repairs that are not carried out by an authorized dealer; if spare parts from other manufacturers are used, or in the case of unusual influences or influences not in compliance with the utilization instructions.

Warranty claims shall not extend the warranty period.

10 Obligations regarding information

- Information on REACH and SVHC can be found on our website in the service area (www.renfert.com/service).

10.1 Disposal instructions

The equipment must be disposed of by a specialist facility. The specialist facility must be informed of any residue in the equipment that is harmful to health.

10.1.1 Disposal information for countries in the EU

To conserve and protect the environment, prevent environmental pollution, and improve the recycling of raw materials, the European Commission has adopted a directive that requires the manufacturer to accept the return of electrical and electronic equipment for proper disposal or recycling.

For this reason, equipment with this symbol may not be disposed of in unsorted domestic waste within the European Union.

Please contact your local authorities for more information on proper disposal.

10.1.2 Special instructions for customers in Germany

Renfert electrical equipment is equipment for commercial use.

This type of equipment may not be disposed of at municipal collection points for electrical equipment but must instead be returned directly to Renfert.

You can find out about the current options for return online at www.renfert.com

**Hochaktuell und ausführlich auf ...
Up to date and in detail at ...
Actualisé et détaillé sous ...
Aggiornato e dettagliato su ...
La máxima actualidad y detalle en ...
Актуально и подробно на ...**

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Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany