
making work easy

FAQ - Deutsch

FAQ - English

Renfert CONNECT

FAQ - frequently asked questions

Made in Germany

[21-XXXXX] [TTTTMMJJJJ] [X]

- 1 App Allgemein
- 2 App Simplex
- 3 App EASY view+

1 App Allgemein

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FAQ - English

- 1 App – General Questions
- 2 App – Simplex
- 3 App – EASY view+

1 App – General Questions

Question / Error / Problem	Cause	Answer / Troubleshooting / Solution
Where can I download the app?		<ul style="list-style-type: none">• Android: from Google Play Store:• iPhone iOS: from the Apple App Store:  
Do I have to pay for the app?		<ul style="list-style-type: none">• No, the app is available free of charge.
What equipment can I operate with the app?		<ul style="list-style-type: none">• We are constantly expanding the app. Generally, it will be possible to operate all Renfert smart equipment using the app. You can check if your unit is already included either directly in the app or on our web page: https://www.renfer.com/connect
What information / data is provided in the app?		<ul style="list-style-type: none">• In the app, you can use certain operating functions for the various units. What's more, the app also ensures that your unit always has the latest software installed. You will also find service videos, FAQs, and operating instructions.
Where is the customer data from the Renfert CONNECT app stored?		<ul style="list-style-type: none">• The data is stored in accordance with GDPR requirements in a database reserved specifically for this purpose.
Who should I ask if I have additional questions about the app?		<ul style="list-style-type: none">• If you have any problems, please contact our service hotline.
Where can I find contact information for the Renfert Service team?		<ul style="list-style-type: none">• Click the two menu strips on the top left of the screen. An overlapping menu is shown. Under Support/Service you can see who to contact.
What operating system version must be installed on my mobile device so that I can use the app?		<ul style="list-style-type: none">• Android operating system: 8.0 or newer• iPhone iOS operating system: 14.0 or newer.
Will the app work on my tablet too?		<ul style="list-style-type: none">• The app is not optimized for use on tablets, but it will work on most tablets with the latest operating system.
What do I need in order to control my Renfert smart equipment using the app?		<ul style="list-style-type: none">• You will need a Renfert CONNECT stick. You can check if the Renfert CONNECT stick is available for your region on our web page: https://www.renfer.com/connect
Is voice control available?		<ul style="list-style-type: none">• No, voice control is not planned at this time.

Question / Error / Problem	Cause	Answer / Troubleshooting / Solution
In what countries is Renfert CONNECT available?		<ul style="list-style-type: none"> The app is available worldwide. You can check if the Renfert CONNECT stick is available for your region on our web page: https://www.renfert.com/connect
Can I upgrade my existing Renfert equipment to use Renfert CONNECT?		<ul style="list-style-type: none"> This requires that the equipment is intended for use with the Renfert CONNECT stick. Ask our Service team for details.
How do I delete my Renfert CONNECT account?		<ul style="list-style-type: none"> You can delete your account in the app under "My Profile". Your data will be anonymized and can no longer be restored.
What operating systems does Renfert CONNECT support?		<ul style="list-style-type: none"> Android operating system: 8.0 or newer iPhone iOS operating system: 14.0 or newer.
What are the minimum requirements for my Wi-Fi router?		<ul style="list-style-type: none"> The Renfert CONNECT system should work with all commercially available Wi-Fi routers that are not more than 5 years old (approximately). However, Wi-Fi networks must be encrypted for security reasons (WPA or WPA2 encryption). If your Wi-Fi network is not encrypted, it will not be possible to log on your unit.
How are data privacy and data security handled?		<ul style="list-style-type: none"> The data is stored in accordance with GDPR requirements in an online platform database reserved specifically for this purpose.
Are you also planning a LAN interface for the equipment?		<ul style="list-style-type: none"> No, a LAN interface is not planned for the equipment as Renfert CONNECT uses Wi-Fi to exchange data. Routing a LAN cable for every unit is a very time-consuming task. We chose Wi-Fi because it is now almost universally available.
Do I need other equipment in addition to a unit that supports Renfert CONNECT functionality and a Wi-Fi router?		<ul style="list-style-type: none"> Yes. You will need the Renfert CONNECT stick so that your unit can communicate with your router. You can check if the Renfert CONNECT stick is available for your region on our web page: https://www.renfert.com/connect
The signal on my unit / stick is weak because my Wi-Fi router is quite far away. Does this mean that I cannot use Renfert CONNECT after all?		<ul style="list-style-type: none"> It is possible to strengthen the Wi-Fi from your router using a Wi-Fi repeater. The signal should be strong enough if the Wi-Fi repeater is installed between the unit and the router.
Can I connect my equipment using WPS?		<ul style="list-style-type: none"> No, WPS functionality is not planned at this time.
How can I find out if the Wi-Fi signal is strong enough where my unit is installed?		<ul style="list-style-type: none"> You can see how strong the signal is directly at the installation location based on the Wi-Fi coverage of your smartphone or tablet.
What routers (Wi-Fi & LTE) provide particularly good connection quality for equipment with Renfert CONNECT?		<ul style="list-style-type: none"> Generally, most routers are suitable for connecting equipment with the app using the Renfert CONNECT stick. Experience has shown that routers such as Fritz!Box from AVM, Speedport from Deutsche Telekom, and AirPort Express from Apple provide a very strong connection.

Question / Error / Problem	Cause	Answer / Troubleshooting / Solution
<p>What are the Wi-Fi requirements?</p>		<ul style="list-style-type: none"> Renfert CONNECT usually operates with the standard Wi-Fi settings of commercially available routers. For a stable connection, the Renfert CONNECT stick connects with the 2.4 GHz Wi-Fi frequency band of your router. Make sure that this is activated. With newer routers, both frequency bands (2.4 and 5 GHz) are enabled as standard. <p>In some cases, certain settings must be adjusted manually. If necessary, check the following settings on your router. You can find information on configuration settings in the installation manual for your router.</p> <p>If required, e.g., during manual configuration, please note the following tips and settings: Encrypt your Wi-Fi, use security level WPA2 wherever possible or at least WPA. Open networks are not suitable for use with Renfert CONNECT. Renfert CONNECT has been tested using network mode IEEE 802.11 g. Please make sure that your router is configured accordingly.</p> <p>Renfert CONNECT does not work with Wi-Fi networks that require browser registration, for example by entering a user name and password in a web browser.</p> <p>You should have a stable Wi-Fi network with a strong signal and Internet access. Deactivate MAC filtering in your router to avoid saving the MAC addresses of your household appliances and devices with communication enabled. The use of proxies is not recommended.</p>
<p>After I log on, my unit that I have already registered is not visible.</p>	<ul style="list-style-type: none"> Unstable Internet connection. 	<ul style="list-style-type: none"> Make sure that your Internet connection is stable. Sometimes it's sufficient to wait 30 seconds. You may need to log off and on again.
<p>Can I see on the stick if it is connected?</p>		<ul style="list-style-type: none"> Yes. If the stick's open Wi-Fi is being transmitted, the LED flashes green. If it has connected with a router, the LED lights up in blue. If the router is out of range, then the LED flashes blue.
<p>What does the "mysimplelink" Wi-Fi mean?</p>		<ul style="list-style-type: none"> This is an open Wi-Fi transmitted by the stick. It does not enable an Internet connection but is required initially so that the equipment can be located in order to then connect with a router via the app. Once the stick has connected with the router, this open Wi-Fi disappears. Remember that no more than 4 mobile devices can connect with this open Wi-Fi.
<p>The app is not connected to the unit.</p>	<ul style="list-style-type: none"> The unit or the stick have been updated. 	<ul style="list-style-type: none"> Switch the unit off and then on again.
<p>The password reset email has not arrived.</p>	<ul style="list-style-type: none"> The incorrect email address has been entered. 	<ul style="list-style-type: none"> Addresses ending in .com and .de are considered different email addresses by the system. The system does not indicate if the email address has already been registered.
<p>An update to the app was installed and now the app keeps crashing.</p>	<ul style="list-style-type: none"> Junk data in the memory can cause problems. 	<ul style="list-style-type: none"> After uninstalling and reinstalling the app, start the mobile device again.

2 App – Simplex

en

Question / Error / Problem	Cause	Answer / Troubleshooting / Solution
Can I delete my registered Simplex from the app?	<ul style="list-style-type: none">• Currently it is not possible to remove the unit yourself from the app.	<ul style="list-style-type: none">• You have the following options:<ul style="list-style-type: none">- You can delete your account. Your data will then be anonymized, and you can create a new account if required.- You can contact our support team. They can remove the unit from the database and so from your app the next time synchronization occurs.
Where can I find the serial number of my unit in the app?		<ul style="list-style-type: none">• Open the "My Equipment" menu. Then select the corresponding unit. You will find the serial number that you stored during registration on the start page under the name of the unit.
After I log on, my unit that I have already registered is not visible.	<ul style="list-style-type: none">• Unstable Internet connection.	<ul style="list-style-type: none">• Make sure that your Internet connection is stable. Sometimes it's sufficient to wait 30 seconds. You may need to log off and on again.
I don't have a mobile device that can run the app. Is there another way to access Simplex software?		<ul style="list-style-type: none">• Please contact our support team. They will need your email address, the name and address of your dental office or laboratory, and the serial number of your Simplex unit. We will then store your information in the database so that you can access it at any time in the future once you have a suitable mobile phone. Once this record is saved in the database, you will not only receive a link to your software via email, every software update will be sent to you automatically and conveniently via email too.

3 App – EASY view+

Question / Error / Problem	Cause	Answer / Troubleshooting / Solution
Can the unit be accessed multiple times simultaneously?		<ul style="list-style-type: none"> You can use the same account to access the same EASY view+ from different devices. Access from two devices at the same time is not possible. Downloaded images are saved locally, however they remain on the stick and can also be downloaded to a different device at any time.
After I log on, my unit that I have already registered is not visible.	<ul style="list-style-type: none"> Unstable Internet connection. 	<ul style="list-style-type: none"> Make sure that your Internet connection is stable. Sometimes it's sufficient to wait 30 seconds. You may need to log off and on again.
The sequence of the images is random / mixed-up.		
Our dental office / our laboratory has purchased a second EASY view+. It was not possible to add it to the app.		<ul style="list-style-type: none"> Currently it is only possible to register one unit for each account. Create a second account for the second unit.
I have registered my EASY view+ but no images were downloaded during initial synchronization.	<ul style="list-style-type: none"> Please check the color of the LED on the Renfert CONNECT stick. If it flashes green, the unit was not switched off and on again. 	<ul style="list-style-type: none"> Switch the EASY view+ off and on again. Contact our support team. You will find the number in the app.
I moved an image to a folder and renamed it. In the gallery, it is still showing the old name.	<ul style="list-style-type: none"> Compare for synchronization. 	<ul style="list-style-type: none"> Currently EASY view+ needs to compare using the gallery so that images that have already been synchronized are not downloaded again. A solution to this problem will be provided automatically via an update as development of the app continues.
It was not possible to upload an image to EASY view+.		<ul style="list-style-type: none"> Check the format of the image. Due to the amount of data, it is currently only possible to upload JPG and PNG files. If you convert your image to this type of format, uploading should be possible without problems.
I uploaded an image but can't find it under this name in the splitscreen mode of my EASY view+.		<ul style="list-style-type: none"> To ensure that images are not uploaded twice, the app renames the image when uploading. Look carefully right at the start or end of the image. There should be a long sequence of digits (e.g., 6276B04A-DE7E-430F-BB3E).
Where can I find the serial number of my unit in the app?		<ul style="list-style-type: none"> Open the "My Equipment" menu. Then select the corresponding unit. You will find the serial number of your unit on the start page under the name of the unit.
I deleted an image in the app but it's still showing on the EASY view+ stick and was downloaded again at the next synchronization.		<ul style="list-style-type: none"> To prevent images from being accidentally deleted from the stick it is currently not possible to delete images from the stick using the app. You can delete them at any time in the usual way via your PC.
In the device setup in the app, setup stops at step 3 (network settings).	<ul style="list-style-type: none"> Unstable connection to the network. 	<ul style="list-style-type: none"> Go back to the network settings and set your own network. Then repeat the step.
Can I connect my unit to another router in another room?		<ul style="list-style-type: none"> Yes, this is possible. To do so you must press the right mouse button on EASY view+ and start the "RESET Wi-Fi" function. You must also disconnect the unit accordingly in the app. You can then connect the unit again in the app and via the other router using "Add device".