making work easy

de
en
fr
it
es
pt
tr
ru
pl
zh
ja
ko
ar

POWER steamer water softener

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making work easy

POWER steamer water softener

Made in Germany

Inhalt

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Contents

en

1	Introduction	3
1.1	Symbols used	3
2	Safety	3
2.1	Intended use	3
2.2	Improper use	3
2.3	Ambient conditions for storage and transport	3
2.4	Hazard and warning information	3
2.4.1	General information	3
2.4.2	Specific information	4
2.5	Authorized persons	4
2.6	Disclaimer	4
3	General description	4
3.1	Assemblies and functional elements	5
3.2	Scope of delivery	5
3.3	Accessories	5
4	Startup	6
4.1	Unpacking	6
4.2	Installation	6
4.3	Blend level / capacity of the filter cartridge	6
4.3.1	Using the test strips	6
4.3.2	Setting the blend level	6
4.4	Assembling the filter head / flowmeter / filter cartridge	7
4.5	Programming and display unit	8
4.5.1	Entering the filter capacity	9
5	Operation	9
5.1	Replacing a filter cartridge	9
5.1.1	Resetting the counter / checking the filter capacity	10
5.2	Programming and display unit	10
5.2.1	Selecting a menu item	10
5.2.2	Filter status display	10
5.2.3	Reset the filter capacity (when replacing the filter)	10
5.2.4	Entering the filter capacity	11
5.2.5	Counter reading (when inserting the cartridge)	11
5.2.6	Alarm	11
5.2.7	Setting the unit (liters or gallons)	11
5.2.8	Setting the pulse per liter	12
6	Cleaning / Maintenance	12
6.1	Flushing the system	12
6.2	Spare parts	12
7	Troubleshooting	13
8	Technical data	13
8.1	Filter unit	13
8.2	Flowmeter	13
9	Guarantee	14
10	Disposal information	14
10.1	Disposing of consumables	14
10.2	Disposing of the equipment	14
10.2.1	Disposal information for countries in the EU	14

1 Introduction

1.1 Symbols used

Symbols with the following meanings are used in these instructions or on the equipment:

Danger
Immediate danger of injury. Consult accompanying documents!

Electrical current
Danger due to electrical current.

Caution
Risk of damage to the equipment if the instruction is not followed.

Notice
Provides the operator with useful information that eases use.

The equipment complies with applicable EU directives.

The equipment is subject to the provisions of the EU directive 2002/96/EC (WEEE directive).

- List, particular attention should be paid

- List
 - Subordinate list

⇒ Instructions / appropriate action / input / operational sequence:

You are asked to carry out the specified action in the order indicated.

◆ Result of an action / response by the equipment / response by the program:

The equipment or program responds as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur.

2 Safety

2.1 Intended use

The POWER steamer water softener is looped in between the water connection and the POWER steamer 2.

2.2 Improper use

Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the equipment, increase the risk of serious injury and lead to damage to the environment or the product itself.

2.3 Ambient conditions for storage and transport

The following ambient conditions must be observed during storage and transport:

- Ambient temperature of -15 – +45 °C [-4 – +140 °F]
- Maximum relative humidity 80 %

2.4 Hazard and warning information

2.4.1 General information

- If the unit is not used in accordance with the supplied operating instructions, the intended safety can no longer be guaranteed.
- Check connection cables (such as power cords), tubes and housing (such as the key pad) regularly for damage (such as kinks, cracks and porosity) or signs of ageing.
Do not operate equipment with damaged connection cables, tubes or housing parts or other defects!
- Only operate the equipment under supervision.
- Please observe the national accident prevention regulations!
- Information on REACH and SVHC is available on our website at www.renfert.com, in the Support area.

2.4.2 Specific information

- Only cold water of drinking water quality may be used to supply the system.
- When supplying with softened water, only use filters for drinking water with > 4° dH total hardness.
- Store components in a dry location.
- The system must not come into contact with chemicals, solvents, or any fumes.
- Before the filter system is put into operation, the consumer unit supplied must be free from lime-scale.
- No copper pipes or galvanized / nickel-plated pipes or connecting parts may be installed between the filter system and the consumer unit.
- The filter cartridge must not be opened or mechanically damaged.
- Water that is already microbiologically contaminated or of unknown microbiological quality may not be used under any circumstances without sufficient disinfection before or after the filter.
- The filter head and the wall bracket must be replaced after 5 years of use / at the latest 6 years after the production date. This also applies to hoses and seals. Pay attention to the date of the production stamp.
- All components must be installed in accordance with country-specific guidelines. DIN 1988 must be observed for the assembly and operation of the system.
- If the system pressure is greater than 8 bar, a pressure reducer must be installed upstream of the filter system.
- The sensor unit must not be exposed to any mechanical loads. Pay particular attention to the lever effect of hoses and bend radii. If necessary, provide support for supply lines!
- Only use commercially available AA batteries (2 x) for the sensor unit.
- Ensure that the battery cover of the sensor unit is properly closed. Otherwise, foreign objects or moisture could penetrate.

2.5 Authorized persons

Installation and maintenance of the filter systems may only be carried out by trained and authorized personnel.

Operation and maintenance of the equipment may only be performed by trained personnel.

Any repairs not described in these operating instructions may only be performed by a qualified electrician.

2.6 Disclaimer

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- The product is employed for any purposes other than those specified in the operating instructions.
- The product is altered in any way other than those alterations described in the operating instructions.
- The product is not repaired by an authorized dealer or if original Renfert spare parts are not used.
- The product continues to be used despite obvious safety defects or damage.
- The product is subjected to mechanical impact or is dropped.
- The product is installed incorrectly or not used as specified in the operating instructions.

3 General description

With the POWER steamer water softener, the carbonate hardness of the water for the POWER steamer 2 is reduced in the flow process using ion-selective filter media. The reduction of hardness can be adapted to the quality of the local drinking water via the blending unit in the filter head.

In addition, the filter material binds heavy metal ions such as copper and cadmium.

The integrated activated carbon block reduces unwanted turbidity, organic impurities, aromas, and flavorings as well as chlorine residue in the filtrate and in the blended water.

3.1 Assemblies and functional elements

- | | |
|--|--|
| 1 Programming and display unit | 7 Blending levels |
| 2 Filter head with wall bracket | 8 Filter cartridge |
| 3 Inlet | 9 Key for blending adjustment |
| 4 Outlet | 10 Sensor unit |
| 5 Purge/relief valve | 11 Holder for programming and display unit |
| 6 Flushing connection with flushing hose | 12 Wall bracket for programming and display unit |


Fig. 1

3.2 Scope of delivery

- 1 Filter head with flushing hose
- 1 Flushing hose
- 1 Connection hose
- 2 Adapters 3/4" - 3/8" incl. gaskets
- 1 Filter cartridge M for POWER steamer water softener
- 1 Flowmeter (items 1, 10, 11, 12, Fig. 1) for POWER steamer water softener
- 1 Operating instructions
- Test strips for determining the water hardness

3.3 Accessories

- 1846 0101 Flowmeter for POWER steamer water softener
- 1846 0102 Filter head for POWER steamer water softener
- 1846 0103 Filter cartridge M for POWER steamer water softener
- 1846 0104 Filter cartridge L for POWER steamer water softener

4 Startup

4.1 Unpacking

- ⇒ Remove the equipment and accessories from the package.
- ⇒ Check the delivery for completeness (compare with the scope of delivery).

Necessary batteries 2x AA are not included in the scope of delivery.

4.2 Installation

Before you begin installation:

A shut-off valve, e.g., a water faucet, must be installed upstream of the filter system.

- ⇒ Determine a suitable location for installing the filter system:
The location of the system must be frost-proof and protected against direct sunlight.
- ⇒ Determine the carbonate hardness and set the corresponding blend level (see Chap. 4.3 and 4.4).
- ⇒ Before the filter system is put into operation, the consumer unit supplied must be free from limescale.
- ⇒ Check the filter system and accessories for damage before assembly. This applies in particular to the O-rings and seals.
- ⇒ If stored below 0° C, the filter cartridge must be stored for at least 24 hours at ambient temperature at the installation site before installation.

4.3 Blend level / capacity of the filter cartridge

The required blend level and the capacity of the filter cartridge depends on the carbonate hardness of the local tap water.

Use the test strips provided for determining the carbonate hardness.

4.3.1 Using the test strips

- ⇒ Briefly immerse the test strips (1 sec.).
- ⇒ Gently shake off excess liquid from the test strip.
- ⇒ Read the result after 1 minute:

◆ 4 green	
◆ 1 red	
◆ 2 red	
◆ 3 red	
◆ 4 red	

4.3.2 Setting the blend level

The blend valve in the filter head of the system enables precise adjustment of the carbonate hardness in the filtered water. As a result, the required blend level for the POWER steamer 2 can be set individually and the performance of the filter cartridges can be optimally utilized.

To set the blend level, proceed as follows:

- ⇒ From the table, take the blend level value that corresponds to the carbonate hardness that has been determined.
- ⇒ Press the key (9) for blending adjustment and turn the blending adjustment dial to the corresponding level.
- ⇒ Remove the blending adjustment key after setting has been completed and store the key in a safe place for possible future adjustments to the settings.


Fig. 2

Test strips					
German degrees ($^{\circ}\text{dH}$) / degree of General Hardness (dGH)	< 3	> 4	> 7	> 14	> 21
English degrees (Clark degrees) ($^{\circ}\text{e}$)	< 3.75	> 5	> 8.75	> 17.5	> 26.25
French degrees ($^{\circ}\text{fH}$)	< 5.4	> 7.2	> 12.6	> 25.2	> 37.8
Parts per million (ppm CaCO_3)	< 53.4	> 71.2	> 125	> 249	> 374
Blend level	3	3	2	2	1
Capacity Filter cartridge M in liters	3570	3060	1440	990	510
Capacity Filter cartridge L in liters	6670	5710	2690	1840	940

4.4 Assembling the filter head / flowmeter / filter cartridge

The flowmeter for the POWER steamer water softener consists of the following components:

- Sensor unit
- Programming and display unit

Assembly of the product comprises the following steps:

- Assembly of the filter head.
- Assembly of the sensor unit.
- Connection to the water supply and to the POWER steamer 2.
- Assembly of the filter cartridge.
- Assembly and connection of the programming and display unit (Chap. 4.5).

⇒ Disconnect the POWER steamer 2 from the power source.

⇒ Turn off the water supply to the POWER steamer 2 and disconnect the water hose to the POWER steamer 2 from the water supply.

⇒ Securely attach the filter head to a wall in a vertical position using appropriate screws (not included in the scope of delivery).

i *The filter system can also be operated either free-standing or on its side without being secured. Vertical wall mounting is recommended.*

! *When assembling the sensor unit, the adapter and the water hoses, make sure that you use the corresponding gaskets.*

- ⇒ Screw one adapter (3/4" - 3/8") onto the connection of the sensor unit (10, Fig. 1).
- ⇒ Screw one adapter (3/4" - 3/8") onto the inlet of the filter head (3).
- ⇒ Screw the sensor unit to the outlet of the filter head (4) using the union nut.
- ⇒ Connect the water hose for the POWER steamer 2 to the adapter of the sensor unit.


Fig. 3

! *The sensor unit must not be exposed to any mechanical loads.*

- ⇒ Connect the water hose provided to the adapter at the inlet of the filter head and to the water supply.
- ⇒ Note the flow direction – markings on the filter head and sensor unit!
- ⇒ Max. torque 10 Nm.

Fig. 4: Typical installation example


Fig. 4

- ! Only use the hoses provided.
- Do not use hoses or adapters with conical screw connections. They will damage the connections and cause the warranty to be voided.
- ! When using screw adapters:
- Only use matching screw adapters of a suitable length. When screwing on, do not place them axially on the filter head. Unsuitable adapters can damage the connections and cause the warranty to be voided.

- ⇒ Run the flushing hose into a suitable container (e.g., a bucket) or into a drain.
- ⇒ Open the purge valve (5) (OPEN).
- ⇒ Turn on the water supply.
- ⇒ Insert the filter cartridge into the filter head.


Fig. 5

- i** Insert the filter cartridge with its marking at insertion mark A and screw it in as far as it will go, operating mark B.

- ◆ Inserting the filter cartridge vents the system and flushes the filter cartridge:
 - Filter cartridge M: Flush with min. 5 liters of water
 - Filter cartridge L: Flush with min. 10 liters of water


Fig. 6

- ⇒ Close the purge valve (CLOSE).
- ⇒ After the initial installation of the filter system, vent and flush the water drain hose and the POWER steamer 2 with at least 2 liters of water. To do this, open the service opening on the POWER steamer 2 and run the rinsing program 3 times (see operating instructions for the POWER steamer 2).

- ! After installing the system and inserting or replacing a filter cartridge, check all components for leak-tightness. Water must not leak out at any location.

Fig. 7

4.5 Programming and display unit

The programming and display unit can either be attached to the head of the filter unit using the holder (11, Fig. 1) or to a wall using the wall bracket provided (12, Fig. 1).

Attachment to the head using the holder (11, Fig. 1):

- ⇒ Push the holder onto the head unit until it snaps into place.
- ⇒ Open the battery compartment and insert batteries, 2 x type AA.
- ⇒ Pay attention to correct polarity.
- ⇒ Hang the programming and display unit in the holder.
- ⇒ Insert the connection cable of the sensor unit into the port on the programming and display unit (20).


Fig. 8

Attachment using wall bracket (12, Fig. 1):

⇒ Screw the wall bracket securely to a wall using suitable flat-head screws, max. shank diameter 4 mm.

Alternatively:

⇒ Attach the supplied adhesive pad to the back of the wall bracket.

⇒ Completely remove the protective film from the adhesive pad.

⇒ Stick the wall bracket to the required surface in the required position.

⇒ Open the battery compartment and insert batteries. 2 x type AA, pay attention to correct polarity.

⇒ Hang the programming and display unit in the holder.

⇒ Insert the connection cable of the sensor unit into the port on the programming and display unit (20, Fig. 8).

⇒ Press the connection cable into the cable guide (21) to prevent the connection plug from being unintentionally pulled out.


Fig. 9

4.5.1 Entering the filter capacity

The capacity of the filter cartridge must be entered on the flowmeter so that the flowmeter can provide a reminder in good time to replace a used filter cartridge.

⇒ Determine the capacity of the filter cartridge used (M / L) for the local carbonate hardness (see table in Chap. 4.3.2.).

⇒ Enter the capacity in the programming and display unit, see Chap. 5.2.4.

The system is now ready for operation.

5 Operation

5.1 Replacing a filter cartridge

⇒ Run the flushing hose into a suitable container (e.g., a bucket) or into a drain.

⇒ Slowly unscrew the used filter cartridge anti-clockwise. This releases it from the filter head and it can be removed.

◆ During this procedure, the consumer unit and the water supply line are hydraulically disconnected from the filter head and the system is decompressed. In this case, it can happen that a minimal quantity of decompression water can escape from the flushing hose due to pressure peaks. Take this into consideration and place a suitable container underneath.

⇒ Open the purge valve (5, Fig. 5).

⇒ Remove the new filter cartridge from the packaging and check for damage.

⇒ Insert the filter cartridge into the filter head.

⇒ In doing so, pay attention to the insertion mark on the filter head and the marking on the filter cartridge. Then turn the filter cartridge as far as it will go. The mark on the filter cartridge must form a line with the operating mark (Fig. 7) on the filter head.

◆ This vents the system and flushes the filter cartridge:

- Filter cartridge M: Flush with min. 5 liters of water

- Filter cartridge L: Flush with min. 10 liters of water

⇒ Close the purge valve.

⇒ Enter the date of installation and the time of replacement in the corresponding fields on the sticker on the filter cartridge (a ballpoint pen can be used).

After installing the system and inserting or replacing a filter cartridge, check all components for leak-tightness. Water must not leak out at any location.

NOTICE: The rinsing water is initially milky and cloudy. This is due to the dispersed air and clears after a short time.

Check the correct position of the filter cartridge in the filter head using the markings on the filter head and filter cartridge (7, Fig. 1). When inserting the filter cartridge, pay attention to the position of the sticker. This should face forward so that all the necessary information is visible.

If the filter cartridge is removed from the filter head without being replaced by a new one, the water supply to the filter head must be shut off.

5.1.1 Resetting the counter / checking the filter capacity

If a filter cartridge of the same size was used when replacing the filter cartridge, only the counter needs to be reset.

To reset the filter capacity after replacement, see Chap. 5.2.3.

If a different size is used, the filter capacity must be determined and entered again (see Chap. 4.5.1).

5.2 Programming and display unit

The display of the display unit is switched off 30 seconds after the last key is pressed.

Pressing any key switches the display on again and the filter status is shown.

5.2.1 Selecting a menu item

By pressing the menu key, the following menu items are displayed one after another.

Key	Display
	The display is switched off. ---
	Press the "MENU" key: - Display: Filter status. 1004.0 L
	Press the "MENU" key: - Display: Flow. 14.9 l/min
	Press the "MENU" key: - Display: Total counter. Σ 3560 L
	Press the "MENU" key: - Display: Set filter capacity. CAP SET 1004 L
	Press the "MENU" key: - Display: Counter reading, 0 = Current cartridge was inserted at a counter reading of 287 l. 0 Σ 287 L
	Press the "MENU" key: - Display: Alarm. ON
	Press the "MENU" key: - Display: Settings. L

5.2.2 Filter status display

Key	Display
	The display is switched off. ---
	Press the "MENU" key: - Display: Filter status. OR - Display: Filter status, filter capacity specified has been reached or the maximum period of use of 12 months exceeded. Replace filter! 1004.0 L 72.1 L

5.2.3 Reset the filter capacity (when replacing the filter)

Key	Display
	The display is switched off. ---
	Press the "MENU" key: - Display: Filter status. 72.1 L
	Press and hold the "RESET" key for approx. 5 seconds. - The display shows "reset 5s", the countdown runs down. reset 5s
	The filter capacity is now reset and the last value entered is shown. 1004.0 L

After the reset, a timer is started that indicates replacement of the filter after 360 days.

5.2.4 Entering the filter capacity

Key		Display
	Press the "MENU" key until CAP SET is shown. - Display: Set filter capacity.	CAP SET 1004 L
	Press and hold the "MENU" key for approx. 5 seconds. - The active digit flashes.	CAP SET 01004 L
	Press the "MODIFY" key: - Set the active digit.	CAP SET 01005 L
	Press the "MENU" key: - Confirm the active digit and move to the next digit. - The active digit flashes (repeat the procedure until the last digit).	CAP SET 01005 L

Based on the local carbonate hardness of the tap water, use the table in Chapter 4.3 to determine the resulting capacity of the filter cartridge.

5.2.5 Counter reading (when inserting the cartridge)

It is possible to call up the total counter reading when replacing the last 5 filter cartridges.

Key		Display
	Press the "MENU" key until the counter reading is shown. - Display of the counter reading at which the current cartridge was inserted.	0 $\Sigma 287 \text{ L}$
	Press and hold the "MENU" key for approx. 5 seconds. - The digit "-1" is shown = Total counter reading at the last but one replacement of the cartridge.	-1 $\Sigma 5 \text{ L}$
	Press the "MENU" key: - The digit "-2" is shown, it is possible to continue to press until "-4".	-2 $\Sigma 5 \text{ L}$

Output of the counter reading (total number of liters) of the last 5 cartridge replacements.

5.2.6 Alarm

Key		Display
	Press the "MENU" key until the alarm is shown. - Display: Alarm.	ON
	Press and hold the "MENU" key for approx. 5 seconds. - "ON" or "OFF" is shown.	ON
	Press the "MODIFY" key: - Switch between "ON" and "OFF".	ON
	Press the "MENU" key: - The current setting is applied.	OFF

When switched on (beep on), an acoustic signal sounds when the remaining capacity is 10% or after a maximum period of use of 12 months when water is drawn off.

5.2.7 Setting the unit (liters or gallons)

Key		Display
	Press the "MENU" key until the settings are shown. - Display: Settings.	L
	Press and hold the "MENU" key for approx. 5 seconds. - "L" or "Gal" flashes.	L
	Press the "MODIFY" key: - Switch between "L" and "GAL".	L
	Press the "MENU" key: - The current setting is applied.	Gal

5.2.8 Setting the pulse per liter

For service and expansion purposes only!
Only change the setting if instructed to do so by Renfert!

Key	Display
	Press the "MENU" key until the settings are shown. - Display: Settings.
	Press and hold the "MENU" key for approx. 5 seconds. - "L" or "Gal" flashes.
	Press and hold the "MODIFY" key for approx. 10 seconds. - Pulse per liter is shown, the active digit flashes.
	Press the "MODIFY" key: - Set the active digit.
	Press the "MENU" key: - Confirm the active digit and move to the next digit. - The active digit flashes (repeat the procedure until the last digit).

6 Cleaning / Maintenance

Safe functioning of the system is only possible if the filter cartridge is replaced regularly. The replacement cycles depend on the carbonate hardness of the drinking water, usage, and the blend level set.

We recommend replacing the filter cartridge after 6 months, but no later than after 12 months.

The operator is required to check the system for leaks on a daily basis.

When replacing the filter cartridge, all parts must be checked for contamination and damage. Damaged parts must be replaced and contamination removed.

If the system is not used for an extended period, the filter cartridge must be replaced.

6.1 Flushing the system

The system must be flushed thoroughly after maintenance work or if it has not been used for an extended period.

- ⇒ Run the flushing hose into a suitable container (e.g., a bucket) or into a drain.
- ⇒ Open the purge valve (5, Fig. 5).
- ⇒ Turn on the water supply.
- ⇒ Allow the quantity of water specified in the table below to flow through.
- ⇒ Close the purge valve (CLOSE).

	Flush volume after 1 week of stagnation:	Flush volume after 4 weeks of stagnation:
Filter cartridge M for POWER steamer water softener	3 liters	15 liters
Filter cartridge L for POWER steamer water softener	5 liters	25 liters

- ⇒ Then carry out 3 rinsing cycles using the POWER steamer.

6.2 Spare parts

You can find wear parts and spare parts in the spare parts list available online at www.renfert.com/p918.

Enter the following item number: 18460100.

The components excluded from the warranty (wear parts, consumables) are marked in the spare parts list.

7 Troubleshooting

en

Errors	Cause	Corrective action
POWER steamer is not filled.	• Water connections assembled against the direction of flow on the filter head.	• Check the assembly of the water connections. The direction of flow is indicated on the filter head and the sensor unit.
	• Water supply is turned off.	• Turn on the water supply.
A negative value is shown when the "Filter status" (Chap. 5.2.2) is displayed.	• The filter capacity was not reset (see Chap. 5.2.3) the last time the filter was replaced.	• Reset the filter capacity (see Chap. 5.2.3). Notes: If the filter cartridge then remains in operation until the filter replacement display is shown, this may be considerably longer than intended for the filter cartridge. If necessary, replace the filter cartridge with a new one before the filter replacement display is shown.

8 Technical data

8.1 Filter unit

Item no.	1846 0100 + 1846 0103	1846 0100 + 1846 0104
Dimensions		
Height of the entire system / mm [inch]	475 [18.7]	410 [16.1]
Height of the filter cartridge / mm [inch]	425 [16.7]	360 [14.2]
Diameter of the filter cartridges / mm [inch]	95 [3.7]	136 [5.6]
Minimum distance to the floor / mm [inch]	40 [1.6]	40 [1.6]
Weight of the filter cartridge / kg [lbs]	1.8 [4.0]	3.2 [7.1]
Operating data		
System pressure (without pressure surges) / hPa (bar) [psi]	2000 - 8000 (2 - 8) [29 - 116]	
Water temperature / ambient temperature / °C [°F]	4 - 30 [39 - 86]	

8.2 Flowmeter

Item no.	1846 0101
Sensor unit	
Connection thread	Inlet 3/8" BSP IG / Outlet 3/8" BSP AG
Inlet pressure (min/max) / hPa (bar) [psi]	1000 - 8000 (1 - 8) [14.5 - 116]
Flow range / l/h [gal/l]	15 - 500 [3.96 - 132]
Pressure loss / hPa (bar) [psi]	≤ 500 (0.5) [7.25]
Measurement (horizontal installation) / %	+/- 5
Water temperature (min/max) / °C [°F]	4 - 30 [37.4 - 86]
Ambient temperature (min/max) / °C [°F]	4 - 40 [37.4 - 104]
Dimensions (W x H x D) / mm [inch]	80 x 39 x 42 [3.15 x 1.54 x 1.65]
Installation position	Horizontal recommended
Programming and display unit	
Graphic display	5 digits in liters and gallons
Counter	Decreasing from 99999 to -9999
Dimensions (W x H x D) / mm [inch]	80 x 54 x 29 [3.15 x 2.13 x 1.14]
Power supply, battery type	2 x AA

This unit corresponds to Part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) This unit may not cause harmful interference and
- (2) This unit must accept any interference received, including interference that may cause undesired operation.

The label with the FCC ID and logo is located inside the battery compartment.

9 Guarantee

Renfert offers a **24-month warranty** on all parts of the equipment, provided that the equipment is properly used.

Warranty claims may only be made upon presentation of the original sales receipt from the authorized dealer.

Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the warranty. These parts are marked in the spare parts list.

The warranty is voided in the case of improper use; failure to observe the operating, cleaning, maintenance, and connection instructions; in the case of independent repairs or repairs that are not carried out by an authorized dealer; if spare parts from other manufacturers are used, or in the case of unusual influences or influences not in compliance with the utilization instructions.

Warranty claims shall not extend the warranty period.

10 Disposal information

10.1 Disposing of consumables

10.2 Disposing of the equipment

The equipment must be disposed of by a specialist facility. The specialist facility must be informed of any residue in the equipment that is harmful to health.

10.2.1 Disposal information for countries in the EU

To conserve and protect the environment, prevent environmental pollution, and improve the recycling of raw materials, the European Commission has adopted a directive that requires the manufacturer to accept the return of electrical and electronic equipment for proper disposal or recycling.

For this reason, equipment with this symbol may not be disposed of in unsorted domestic waste within the European Union.

Please contact your local authorities for more information on proper disposal.

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Renfert GmbH, Industriegebiet, 78247 Hilzingen / Germany
c/o Renfert GmbH
Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany