making work easy

- de
- en
- fr
- it
- es
- pt
- tr
- ru
- pl
- zh
- ja
- ko
- ar

fi

sv

cs









uk








SILENT TS

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














Made in Germany

21-9191 01032023

making work easy

SILENT TS

Made in Germany

Inhalt

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Contents

en

1	Introduction	3
1.1	Symbols used	3
2	Safety	3
2.1	Intended use	3
2.2	Improper use	4
2.3	Ambient conditions for safe operation	4
2.4	Ambient conditions for storage and transport	4
2.5	Hazard and warning information	4
2.5.1	General information	4
2.5.2	Specific information	5
2.6	Authorized persons	5
2.7	Exclusion of liability	5
3	Product description	6
3.1	General description	6
3.2	Assemblies and functional elements	6
3.3	Scope of delivery	6
3.4	Accessories	6
4	Startup	7
4.1	Setup	7
4.2	Connection to the suction point	7
4.3	Electrical connection	7
4.4	Connecting the Electrical Devices	7
4.5	External air duct	8
5	Operation	8
5.1	Switching on	8
5.2	Suction Power	8
5.3	Select operating mode: Automatic operation / continuous operation	8
5.4	Automatic switch-on for automatic operation	9
5.5	"Change filter" - display	9
5.5.1	Adjust the sensitivity of the „Replace filter“ indicator	9
5.5.2	Activating / deactivating the „Replace filter“ indicator tone	10
5.6	Operation with the SILENT flow sensor	10
6	Cleaning / Maintenance	10
6.1	Seals	10
6.2	Changing the dust bag	11
6.3	Filters	11
6.3.1	Replacing the fine filter	11
6.3.2	Replacing the Exhaust Filter	11
6.3.3	Replacing the Electronics Filter	12
6.4	Fuses	12
6.5	Self-diagnosis	12
6.6	Replace the suction turbine	12
6.7	Replacing the Exhaust Air Labyrinth	13
6.8	Spare parts	13
7	Troubleshooting	13
8	Technical data	15
9	Warranty	16
10	Disposal information	16
10.1	Disposing of consumables	16
10.2	Disposing of the equipment	16
10.2.1	Disposal information for countries in the EU	16

1 Introduction

1.1 Symbols used

The following symbols are employed in these instructions and on the unit itself:

Danger

Immediate risk of injury.

Electrical voltage

Danger due to electrical voltage.

Caution

Risk of damage to the equipment if the instruction is not followed.

Notice

Provides the operator with useful information that eases use.

The equipment complies with applicable EU directives.

This product complies with the relevant UK legislation.

See the UKCA declaration of conformity available online at www.renfert.com.

The equipment is subject to the provisions of the EU directive 2002/96/EC (WEEE directive).

► List, particular attention should be paid

- List
- List

⇒ Instructions / appropriate action / input / operational sequence:

You are asked to carry out the specified action in the order indicated.

- ◆ Result of an action / response by the equipment / response by the program:

The equipment or program responds as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur.

2 Safety

2.1 Intended use

This device is designed to extract dry, non-explosive dust.

The unit is intended solely for use in a commercial dental laboratory and dental practice.

The intended use also includes compliance with the conditions specified by the manufacturer concerning operation and maintenance.

The extraction unit can be used in combination with a SILENT flow sensor within the scope of measures to protect against exposure to dust *) as defined by the German ordinance on hazardous substances ("GefStoffV") / respective employers' liability insurance association.

For this purpose, the information on low-dust disposal of the extracted material (chapter "Emptying the dust drawer") and on low-dust replacement of the fine filter (chapter "Replacing the fine filter") in particular must be observed.

In doing so, both the instructions for use of the extraction unit as well as those of the SILENT flow sensor must be observed.

*) Note for customers in Germany: A validated and recognized dust protection measure in accordance with testing principle GS-IFA-M20 (issued by the German Institute for Work and Health IFA) shall only be deemed to have been provided when all components, i.e. the flow rate monitoring unit such as the SILENT flow sensor, the extraction unit, and the detection unit, are validated and recognized.

In the case of multi-workbench extraction units and the use of Y-junctions or extractor switches, every suction channel / every suction point must be fitted with a flow rate monitoring unit.

2.2 Improper use

! The unit is not suitable for operation on dental CAM systems! Visit www.renfert.com for information on Silent series extractors suitable for CAM systems.

Fire-promoting, easily flammable, red-hot, burning or explosive materials must not be suctioned into the device.

It is not permitted to suction liquids.

The unit is not intended for private, household use.

Any use other than specified in these instructions is deemed improper use.

The manufacturer shall not be liable for resulting damages.

Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the equipment, increase the risk of serious injury and lead to damage to the environment or the product itself.

2.3 Ambient conditions for safe operation

The unit may only be operated:

- Indoors
- Up to an altitude of 2,000 m above sea level,
- At an ambient temperature of between 5 - 40 °C [41 - 104 °F] *),
- At a maximum relative humidity of 80 % at 31 °C [87.8 °F], decreasing linearly to 50 % relative humidity at 40 °C [104 °F] *),
- With mains power where the voltage fluctuations do not exceed 10 % of the nominal value,
- Under contamination level 2 conditions
- Under over-voltage category II conditions

*) At between 5 - 30 °C [41 - 86 °F], the unit can be operated at a humidity of up to 80 %. At temperatures of 31 °C - 40 °C [87.8 °F - 104 °F] the humidity must reduce proportionately to ensure that the unit can be operated (e.g. at 35 °C [95 °F] = 65 % humidity, at 40 °C [104 °F] = 50 % humidity). The unit may not be operated at temperatures above 40 °C [104 °F].

2.4 Ambient conditions for storage and transport

The following ambient conditions must be observed during storage and transport:

- Ambient temperature - 20 – + 60 °C [- 4 – + 140 °F].
- Maximum relative humidity 80 %.

2.5 Hazard and warning information

2.5.1 General information

- ▶ If the unit is not operated in accordance with the supplied instructions for use, the intended safety can no longer be guaranteed.
- ▶ The unit may only be operated using a power cord equipped with a plug appropriate to the local power supply. Any necessary modifications may only be carried out by a qualified electrician.
- ▶ Only operate the unit if the information on the nameplate corresponds to the specifications of your local mains power supply.
- ▶ The unit may only be plugged into electrical sockets that are connected to the protective conductor system.
- ▶ The power plug must be easily accessible.
- ▶ Disconnect the unit from the mains power supply before carrying out work on the electrical parts.
- ▶ Check connection cables (such as power cords), tubes and housing (such as the key pad) regularly for damage (such as kinks, cracks, and porosity) or signs of aging. Do not operate equipment with damaged connection cables, hoses or housing parts or other defects!
- ▶ Defective equipment must be taken out of service immediately. Disconnect the power plug and secure the equipment against further use. Send the equipment in for repair!
- ▶ Only operate the equipment under supervision.
- ▶ Please observe the national accident prevention regulations!
- ▶ It is the responsibility of the operator to ensure that national regulations are complied with during operation and regarding recurring safety inspection of electrical equipment.
For Germany these are the regulation 3 by DGUV (German Statutory Accident Insurance) in relation with VDE 0701-0702 (Association for Electrical, Electronic and Information Technology).
- ▶ Information on REACH and SVHC is available on our website at www.renfert.com, in the Support area.

2.5.2 Specific information

- ▶ The German ordinance on hazardous substances ("GefStoffV") or equivalent national regulations must be observed during operation, cleaning, and maintenance.
- ▶ The mains socket on the device is only designed for purposes specified in the operating instructions. Connecting other devices may cause material damage.
- ▶ Switch off the electrical device before connecting an electrical appliance to the appliance outlet.
- ▶ Read the operating instructions for the electrical appliance which is to be connected and observe the stipulated safety instructions.
- ▶ Observe national regulations and permitted dust exposure levels in the working environment. Refer to your national association for occupational safety and health or the responsible authority.
- ▶ Please observe the Safety Data Sheets of the materials to be extracted.
- ▶ Wear personal protective equipment when suctioning hazardous materials.
- ▶ Depending on the type of extracted material, wear appropriate personal protective equipment when emptying the dust drawer or cleaning.
- ▶ When disposing of the extracted material or used filters, observe local provisions and accident prevention regulations!
- ▶ Suctioning may only be performed with the dust drawer closed.
- ▶ Do not operate without a suction hose.
- ▶ Do not extract flammable or explosive gasses, fumes or dust.
- ▶ The following applications are associated with increased risks and are therefore not permitted: When cleaning milling chambers and SLM printing systems by suction, a sufficiently high concentration and purity (i.e., unmixed with other dental dusts such as plaster, resin) can result in spontaneous combustion or the milling dust or printing powder exploding because of a potential exothermic reaction (e.g., due to oxidation).

An alternative cleaning method must be selected, particularly for the following milling or printing media (e.g., manual cleaning):

- Wood
- Titanium / titanium-aluminum
- Light metals and light metal alloys (e.g., aluminum, magnesium)
- Cobalt-chrome powder (e.g., for use in SLM systems)

If large quantities of light metals such as titanium alloy are processed (e.g., with sandpaper), producing very fine grinding dust, a sufficiently high concentration and purity may result in spontaneous combustion because of a potential exothermic reaction.

- ▶ Do not suction hot materials.
- ▶ Do not suction liquids.
- ▶ If the dust extractor is employed to extract hazardous materials, appropriate personal protective gear must be worn and steps must be taken to ensure that the exhaust air is properly ventilated. Refer to the safety data sheets for the corresponding requirements.
- ▶ Dispose of extracted material in accordance with statutory regulations.

2.6 Authorized persons

Operation and maintenance of the equipment may only be performed by trained personnel.

Minors and pregnant persons may only operate and service the extraction unit if they are wearing appropriate protective equipment, in particular when suctioning hazardous materials.

Any repairs not described in this user information may only be performed by a qualified electrician.

2.7 Exclusion of liability

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- ▶ The product is employed for any purposes other than those specified in the instructions for use.
- ▶ The product is altered in any way other than those alterations described in the instructions for use.
- ▶ The product is not repaired by an authorized dealer or if original Renfert spare parts are not used.
- ▶ The product continues to be used despite obvious safety defects or damage.
- ▶ The product is subjected to mechanical impact or is dropped.

3 Product description

en

3.1 General description

This is a workplace and appliance suction unit for extracting dust in dental laboratories.

The extraction unit can be operated both manually and automatically, depending on the connected, electrical, dust-generating equipment.

3.2 Assemblies and functional elements

See Fig. 4, 5

-  (A) OFF switch
-  (B) ON switch
-  (C) Display LED Automatic operation
-  (D) Operating mode button, continuous / automatic mode
-  (E) Display LED Continuous operation
-  (F) Extraction force / activation sensitivity adjustment switch
-  (G) Extraction power / activation sensitivity display
-  (H) Extraction force / activation sensitivity adjustment switch
-  (K) „Change Filter“ LED indicator
- (L) Mains fuse (Fig. 4)
- (M) Motor protection switch (Fig. 4)
- (N) Device socket (Fig. 4)

3.3 Scope of delivery

- 1 Laboratory extraction Silent TS
- 1 Operating manual
- 1 Suction hose
- 1 Dust bag (inserted ready for use in dust drawer)
- 1 Fine filter (built-in)
- 1 Adapter for hose connection
- 1 Schuko plug (only for 2921-0050)

3.4 Accessories

- 2921-0002 Dust bag kit (5 pieces).
- 2921-0001 Caster set
- 2925-0000 Extractor clamp (without glass pane)
- 2925-1000 Glass pane with attachment (to the extractor clamp)
- 2926-0000 Extractor Switch incl. 4 mufflers (without hose)
- 2921-0003 Set of bushings, 2 pieces
- 2921-0004 External air duct
- 2934-0007 90°-Angled connector Silent
- 15-0823 Suction hose (Yard goods)
- 90003-4305 Adapter for hose connection
- 90003-4240 Suction hose (3 m), cpl.
- 90003-4826 Suction hose antistatic, 3 m, incl. 2 end bushings
- 90003-4314 Y-junction
- 90003-4410 H- + Hepa Filter Silent TS

4 Startup

4.1 Setup

Remove the unit and accessories from the shipping box.

Check the delivery for completeness (compare with the "Scope of delivery").

The extraction unit is a standing device. Do not operate it on its side.

In particular, after mounting the roller set (see accessories), the suction unit may only be operated standing on the floor.

Place the extraction unit so that:

- The exhaust air outlet is not blocked.
- The front is easily accessible for removing the dust drawer.

When inserting the new fine filter, pay attention to correct seating, as otherwise leaks will occur.

- Circle: at least 120 mm diameter
- Rectangle: at least 170 x 65 mm

If the warm exhaust air is transported out of the cabinet to the rear through openings, it must be ensured that the warm exhaust air can escape from there unhindered.

4.2 Connection to the suction point

Caution: Risk of injury!

When shortening the suction hose, make sure to cut the integrated wire as straight as possible.

Long suction hoses, tight bends and kinks greatly reduce the suction power at the suction point.

⇒ Shorten the suction hose if necessary.

⇒ Plug the suction tube into the suction socket (Fig. 1).

⇒ Connect the suction tube to the correct suction point.

If the diameter is not suitable, use an adapter (see Accessories) to avoid loss of suction.

Avoid steep inclines and "sagging" along the hose.

Never operate the Silent TS without an extraction hose.

4.3 Electrical connection

Prior to electrical connection, check that the voltage indicated on the nameplate matches the local voltage supply.

Ensure that the layout of live parts (socket, plug and coupling) and extension cables are arranged so that the protection class is maintained.

- Unroll the power cord.
- Press the OFF switch (Fig. 2a).
- Plug the power cord into the wall outlet (Fig. 3).

4.4 Connecting the Electrical Devices

For the connection of dust-generating devices, a socket for two-pole plugs with protective contact (according to DIN 49441 or NEMA) (N, Fig. 4) is located on the rear of the device.

For devices with 230 V nominal voltage and socket according to DIN 49441, a suitable plug is included (Fig. 32). With the help of this plug, an adapter can be made to the country-specific connector system.

This adapter may only be manufactured by an electrotechnical specialist!

The protective conductor system must not be interrupted by the adapter!

The device socket carries voltage as soon as the Silent TS is connected to the power supply, even if it is switched off. This means that short jobs can be carried out even without operating the extraction system.

When connecting electrical appliances to the suction unit, ensure that the maximum permissible power for connected devices is not exceeded (see chapter 8, Technical data).

4.5 External air duct

The exhaust air can be transported out of the laboratory via an external air duct (see Accessories).

The corresponding installation instructions are included with the external air duct.

When using extraction units in conjunction with an external air duct, a significant quantity of air is extracted from the room per hour.

This can lead to vacuum pressure, which can cause poisonous gases (e.g., carbon monoxide) to be sucked into the (working) area when using gas, liquid, or solid fuels with furnaces dependent on ambient air.

Depending on the structural situation, an additional air supply or vacuum pressure monitoring must therefore be ensured and checked where appropriate by the responsible institutions (e.g., chimney sweep).

5 Operation

Before beginning to work, check that a dust bag has been inserted.

If operated without a dust bag, the intended protection against exposure to dust can no longer be guaranteed.

5.1 Switching on

The extraction unit is switched on and off via the on / off switch (A / B).

Switching on (b), Fig. 2.

Afterwards, the extraction unit is in the last operating mode set.

5.2 Suction Power

The suction power for the suction unit can be set in 4 levels.

The current level is shown in the display (G).

Change the suction power:

Increase suction

Reduce suction

The suction always starts with the last used suction level.

5.3 Select operating mode: Automatic operation / continuous operation

The extraction unit has two operating modes.

The set operating mode is indicated by the displays  / 

- Automatic mode  :

The suction is dependent on the electrical device connected to the appliance outlet (8, Fig. 1).

- Continuous operation  :
Suction running continuously.

⇒ Press operating mode key (D).

- ♦ Switch the operating mode.

To start or stop the unit during permanent operation, use the operating mode key (D). Do not use the on/off switch for this.

If the suction to a connected electrical appliance (e.g. a handpiece) does not react as required (e.g. does not activate, even though a connected device is in operation) then the switch-on threshold for automatic operation must be adjusted.

5.4 Automatic switch-on for automatic operation

Only delivery, the switch-on automatic is set in such a way that operation with the most handpieces is possible.

If the suction to a connected electrical appliance (e.g. a handpiece) does not react as required (e.g. does not activate, even though a connected device is in operation) then the switch-on threshold for automatic operation must be adjusted.

To change:

⇒ Switch suction on

⇒  Press for 3 sec.

♦ LED  and  blinken

♦ in the display a „0“ blinks

⇒ Switch off electrical devices without standby mode.

⇒ Electrical device with a standby mode (e.g. handpieces), switch to standby mode (e.g. with a handpiece, switch only the control unit on without activating the handpiece).

⇒  Press key

♦ in the display „1“ blinks

⇒ An electrical device without standby mode must be switched on.

⇒ Operate the handpiece at the speed at which the suction is to be started and run for approx. 3 - 5 seconds (so that overcurrent during start-up is not detected).

⇒  Press key

♦ An acoustic signal acknowledges the information.

For devices with standby mode, these should be switched on approx. 5 minutes before setting to avoid incorrect measurements.

5.5 "Change filter" - display

♦ Indicator LED (K) lights up, 3 beeps sound.

⇒ Change dust bag (see chap. 6.2)

or

⇒ Replace fine filter (see Chap. 6.3)

If the suction unit continues to be operated despite the "Change filter" display, it stops automatically when the air flow falls below a fixed set value.

At the same time, the LED (K) flashes.

By switching OFF and ON, this detection is deleted and the suction starts again.

Operation with a clogged filter system may cause hazards and the suction unit may be damaged. If the suction stops automatically when the filter is full, be sure to change the filter.

5.5.1 Adjust the sensitivity of the „Replace filter“ indicator

It may be necessary to reduce responsivity if the suction support of the connected appliance has a smaller diameter than the Silent TS. (<< 32.5 mm).

The „Replace filter“ sensitivity can be adjusted within limits.

To do this:

- Switch the Silent TS off;

- Switch the Silent TS on, and in doing so keep the  button firmly pressed until the display begins to blink. (5 = factory set parameter)

-  Press button,

- Reduce level of sensitivity,

- „Change filter“ signal will appear later.

- ==> Dust bag will get very full

-  Press button,

- Increase level of sensitivity,

- „Change filter“ signal appears more quickly.

- ==> Dust bag will get less full

-  Press button,

- Acoustic signal confirms successful adjustment.

5.5.2 Activating / deactivating the „Replace filter“ indicator tone

The „Replace filter“ indicator tone can be activated or deactivated.

To do this:

- Switch the Silent off;
- While holding down the  and  keys, switch the Silent back on and wait for the indicator tone
 - Brief indicator tone: acoustic signal deactivated.
 - Longer indicator tone: acoustic signal active.

5.6 Operation with the SILENT flow sensor

If the extraction unit is operated in combination with a flow rate monitoring unit such as the SILENT flow sensor, the following must be observed:

- **If the flow rate is less than the defined threshold value for a certain period, this is detected by the SILENT flow sensor, the corresponding indicator begins to flash, and a warning tone is output every 30 seconds.**

If the flow rate is too low at the suction point being monitored, safe operation is no longer guaranteed and the level of exposure to dust can be harmful to health!

In this case, the following measures must be taken:

- ⇒ Stop working at all suction points of the affected extraction unit. This applies to the individual suction channels of any multi-workbench extraction unit as well as to all suction points operated via a Y-junction or via extractor switches.
- ⇒ Determine why the flow rate is too low and resolve the issue.

Cause	Measure
Suction level is too low.	• Choose a higher suction level.
Dust bag is full.	• Change the dust bag.
Changing the dust bag does not result in the expected improvement.	• Replace the fine filter. • For low-dust replacement of the fine filter, observe Chap. 6.3.1.
Blockage in the suction point (suction mouth).	• Disconnect the suction hose at the suction point and check whether there is a blockage in the suction point; remove where appropriate.
Blockage in the suction line.	• Check the suction line at all disconnection points / service openings section by section for blockages.
Leakage in the section of the suction line between the SILENT flow sensor and the extraction unit.	• Check the installation of the SILENT flow sensor in accordance with its instructions.

Determining the cause and measures can be carried out in the order shown.

After every measure, you should switch on the extraction unit again and check whether the minimum flow rate is again exceeded so that correct operation is possible again.

- ⇒ Work should only resume after the fault has been rectified and the minimum flow rate is exceeded again.

6 Cleaning / Maintenance

**The unit does not contain any internal parts that require maintenance.
Opening the unit is not permitted other than as described below!**

6.1 Seals

To ensure the proper function of the dust extractor, it is vital that the following three seals:

- Profile seal on the dust drawer (Fig. 12a);
- V-ring seal on the dust drawer (Fig. 12b);
- Seal on the upper front panel (Fig. 20);

are not damaged.

Inspect these seals when replacing the associated filters and replace them if they are found to be damaged (refer to the spare parts list).

6.2 Changing the dust bag

The dust extractor may only be operated with a complete filtration system. The dust bag must be immediately replaced with a new one as soon as the "Change Filter" indicator goes on.

- Pull the front panel forward and off (Fig. 8).
- Pull the dust drawer out towards the front (Fig. 9).
- Leave the dust bag in the dust drawer and take it to the disposal location.
- Remove the dust bag from the drawer.
- Pull off the protective film, seal the filter opening dust-tight (Fig. 10).
- Dispose the dust bag accordingly.

Always comply with all local ordinances governing proper disposal and accident prevention! Depending on the filter contents, protective clothing may be required.

- Insert the new dust bag into the dust drawer, being careful to ensure that the dust bag is properly aligned in the guide grooves and the glue lap faces up (Fig. 11).
- Inspect the dust drawer seal for damage and replace it, as required (Fig. 12a, b).
- Insert the dust drawer, making sure it is straight, and push it fully against its stops. Make sure the dust drawer is seated correctly on the guide (Fig. 13).
- Install the front panel on the bottom hooks first (Fig. 14), then push it into the upper lock bolts, pressing firmly until it locks in place (Fig. 15).

6.3 Filters

NEVER operate the dust extractor without the complete filtration system.

6.3.1 Replacing the fine filter

Examine the fine particle filter regularly, but at least annually, and replace it. It must always be replaced if the "Change Filter" indicator (K) continues to go on, despite the dust bag having been replaced.

There is a risk of increased exposure to dust when replacing the fine filter.

That's why the following measures must be taken:

- Keep the plastic bag ready that contains the new filter.
- Keep an additional dust extraction device or suction tube from another workbench extraction unit ready.
- Wear at least an FFP2 mask for respiratory protection.
- Pull the front panel forward and off (Fig. 16).
- Turn the fine particle filter anticlockwise to release it (Fig. 17).
- Pull the fine particle filter straight forwards and out of the unit (Fig. 18).
- Place the fine filters immediately in the plastic bag and seal the bag tightly to prevent dust escaping.
- Ensure proper disposal.

Always comply with all local ordinances governing proper disposal and accident prevention! Depending on the filter contents, protective clothing may be required.

- Insert the new fine filter and push it in fully via the retaining bracket (Fig. 19).
- Turn the fine particle filter clockwise and hand-tighten it in place (Fig. 17).
- Inspect the seal on the front panel for damage and replace it, as required (Fig. 20).
- Snap the front panel in place (Fig. 21, 22).
- Use the additional dust extraction device to suction up dust that has been released.

6.3.2 Replacing the Exhaust Filter

The exhaust filter primarily traps particles rubbed off the suction turbine graphite brushes. The filter should be replaced annually:

- Press on the four lock tabs on the filter cassette and pull the cassette down and off the unit (Fig. 23).
- Properly dispose of the filter and the filter cassette.

Always comply with all local ordinances governing proper disposal and accident prevention! Depending on the filter contents, protective clothing may be required.

- Insert a new filter mat in the filter cassette so that the smooth, compressed side of the filter mat faces down or outwards when the filter is installed.
- Install the new filter cassette on the exhaust air vent.
- Make sure the filter cassette is properly seated and locked in place on all four tabs.

6.3.3 Replacing the Electronics Filter

Cooling air for the unit's electronics passes through the electronics filter. The type of material trapped by the filter and the degree of contamination depend on the ambient laboratory conditions.

The filter should be examined annually and replaced, as necessary:

- Press on the two lock tabs on the filter cassette and pull the cassette off towards the rear of the unit (Fig. 24).
- Properly dispose of the filter and the filter cassette.

Always comply with all local ordinances governing proper disposal and accident prevention! Depending on the filter contents, protective clothing may be required.

- Install a new filter cassette and filter on the opening.
- Make sure the filter cassette is properly seated and locked in place on both tabs.

6.4 Fuses

Instead, two device overload switches, (L, Fig. 4) are provided to protect these units.

Press the button to reset a switch which has tripped.

6.5 Self-diagnosis

The dust extractor is equipped with a self-diagnosis, which checks various functions and displays errors on the display (G).

Before the diagnosis takes place, insert a new dust bag and ensure that the fine filter is clean and the suction pipe is not blocked.

Activate the self-diagnosis:

- Switch Silent TS off;
- Press and hold the operating mode selector switch (D) and switch Silent TS on (Fig. 7a).
- For approx. 3 sec.:
 - all LED displays are illuminated,
 - An „8“ is indicated in the display
 - The signalling device gives an acoustic signal.
- A “d” appears on the display during the unit's self-test (Fig. 7b). During this time, the unit's electronics check various internal elements and functions. The suction turbines will briefly run at Level 1 during this test phase.
- If no errors are detected, the Silent TS automatically goes into the previously selected operating mode at the end of the self-test (after approx. 10 sec.).
- If an error is detected during the self-test, a flashing “E” will appear on the display, together with a number (1 to 3).
- These indicate:
 - E1: Airflow sensor fault; have the unit repaired.
 - E2: Plug on grey motor wire out of its socket (Fig. 26) or suction turbine fault. Replace suction turbine if necessary (refer to the spare parts list).
 - E3: Error in the electronics; have the unit repaired.
- In this case, the Silent TS remains in the self-diagnosis mode until it is switched off at the OFF switch (A).

6.6 Replace the suction turbine

The suction turbine is encased in an encapsulating housing to form a single unit, thus allowing it to be easily replaced without the need for any tools.

Before opening the unit, disconnect it from the mains power supply by unplugging the power cord from the wall outlet.

The motor may be hot. Allow the motor to cool off before proceeding.

The suction turbine may only be operated when it is installed in the unit. Only an authorized electrician may perform a functional inspection of or repair the unit.

Please observe the documents included with the spare part!

- Turn the lock knob 90° anticlockwise (Fig. 25).
- Remove the motor compartment cover.
- Open the fixing clip and remove the motor plug (Fig. 26).
- Loosen the two internal, grey knurled nuts (Fig. 27) and unscrew them approx. 1 cm. These knurled nuts do not need to be completely removed.
- Release the suction turbine by turning it anticlockwise, then pull it straight out towards the back (Fig. 28).

- Install the new suction turbine (Fig. 28) by pushing it in straight and turning it clockwise to lock it in place (Fig. 29). The connecting wire should be on top (Fig. 28a).
- Tighten the internal knurled nuts (Fig. 27).
- Insert the motor plug and close the fixing clip (Fig. 26).
- Check that the lock knob is in the open position (vertical).
- Place the motor compartment cover on the motor compartment. Two guide pins have been provided on the cover to assist in positioning it correctly.
- Turn the lock knob 90° clockwise (Fig. 31).

Always comply with all local ordinances governing proper disposal and accident prevention!

6.7 Replacing the Exhaust Air Labyrinth

Dust from the suction turbine graphite brushes also settles on the walls of the exhaust air labyrinth. Over time, the labyrinth will become discoloured.

Disconnect the power plug before opening the device.

Replacing the exhaust air labyrinth:

- Turn the lock knob 90° anticlockwise (Fig. 25).
- Remove the motor compartment cover.
- Pull the exhaust air labyrinth out towards the back.
- Install the new exhaust air labyrinth by pushing it in straight, making sure it is properly seated. The openings in the exhaust air labyrinth must face towards the rear (Fig. 30) and the metal grating downwards.
- Check that the lock knob is in the open position (vertical).
- Place the motor compartment cover on the motor compartment. Two guide pins have been provided on the cover to assist in positioning it correctly.
- Turn the lock knob 90° clockwise (Fig. 31).

Always comply with all local ordinances governing proper disposal and accident prevention!

6.8 Spare parts

You can find components subject to wear and the spare parts on the spare part list in the internet at www.renfert.com/p918.

The components excluded from the warranty (wear parts, consumables) are marked in the spare parts list. The serial number, date of manufacturing and equipment version are shown on the nameplate of the equipment.

7 Troubleshooting

Error	Possible cause	Corrective action
Dust extractor stops suddenly, for no apparent reason.	• Dust extractor continued to be used, despite "Change Filter" indicator.	• Check and replace the dust bag and fine particle filter (Chap. 6.2 and 6.3).
	• Unit fuse blown.	• Check and reinstall the fuse.
	• Motor protection switch (M, Fig. 4) has been activated due to overload (temperature, electricity).	• Check the motor protection switch, if necessary reset. • Change the motor if repeatedly activated.
	• Motor over-temperature protector switch activated.	• Check dust bag and filter, replace as required. (Chap. 6.2 and 6.3). • Check exhaust air filter, replace as required. (Chap. 6.3.2). • Check whether suction hose is blocked. Remove obstruction as required. • Dust extractor can be restarted after being allowed to cool off for approx. 1 hour.
	• Plug on grey motor wire out of its socket or loose.	• Check plug (Fig. 26).
	• Suction turbine fault.	• Carry out self-diagnosis (chapt. 6.5), at fault „E2“ exchange suction turbine (see spare parts list).
	• In-house fuse blown.	• Replace the in-house fuse. • Check the total power consumption of the dust extractor plus dust-generating devices.

Error	Possible cause	Corrective action
Dust extractor stops suddenly, for no apparent reason. E4 is indicated on the display.	The electronics may be overheated due to excessively high ambient temperature, e.g., if the unit is installed in a cabinet.	- Allow the extractor to cool down. - Attempt to improve the ventilation conditions. Refer also to Chap. 4.1.
After switching on the unit, the display immediately indicates E4.	Temperature sensor fault.	Have the device repaired.
	After overheating, the device has not yet cooled down sufficiently.	Disconnect the device and let it continue to cool down.
The self-test returns a suction turbine error, even though the turbine is running.	Filter full.	Check dust bag and fine particle filter, replace as required. (Chap. 6.2 and 6.3.1).
	Suction hose blocked.	Remove obstruction from the suction hose.
	Suction turbine fault.	Replace the suction turbine (Chap.6.6).
Insufficient extraction force.	Extraction force level too low.	Increase the extraction level (Chap. 5.3).
	Filter blocked.	Check dust bag and fine particle filter, replace as required.(Chap. 6.2 and 6.3.1).
	Suction hose blocked.	Remove obstruction from the suction hose.
	Dust drawer seal not tight.	- Check the proper seat of the dust drawer (Chap. 6.2). - Check the dust drawer seal for damage and replace as required (Chap. 6.1).
	Upper front panel fails to close properly.	- Check the proper seat of the upper front panel (Chap. 6.3.1). - Check the upper front panel seal for damage and replace as required (Chap. 6.1).
	Suction hose damaged.	Inspect the suction hose and replace as required.
“Change Filter” indicator goes on shortly after filters were replaced.	Fine particle filter blocked.	Replace the fine particle filter (Chap. 6.3.1).
	Suction hose blocked.	Remove obstruction from the suction hose.
	Cross-section of the dust-generating device's hose fitting is too small.	Adjust the sensitivity of the „Replace filter“ indicator (Chap. 5.6.1.).
	Dust is too fine.	Adjust the sensitivity of the „Replace filter“ indicator (Chap. 5.6.1.)
	Flow sensor defect.	Carry out self-diagnosis (chapt. 6.5), at fault „E1“ equipment requires professional repair.
Dust bag ruptures.	The “Change Filter” indicator's sensitivity has not been adjusted for the material being extracted.	Reduce the sensitivity of the „Replace filter“ indicator (Chap. 5.6.1.).
Dust extractor fails to start when a dustgenerating device is being used.	Dust extractor not switched on.	Switch the dust extractor on (Chap. 5.1).
	Dust-generating device not plugged into the unit power outlet on the Silent TS.	Connect the dust-generating device to the Silent TS (Chap. 4.4).
	Wrong operating mode selected.	Select the automatic mode (Chap. 5.4).
	Automatic start feature's threshold value too high.	Adjust the automatic start feature (Chap. 5.5).
	Device fuse has blown.	Check and reinstall the fuse.(Chap. 6.4).
	Automatic start feature fails to detect the attached device's operation.	Some older handheld instruments do not generate a sufficiently strong signal for the automatic start feature. Switch the dust extractor to the continuous mode.
Dust extractor starts, although a connected device is not being used.	Automatic start feature's threshold value too low.	Adjust the automatic start feature (Chap. 5.5).

Error	Possible cause	Corrective action
Dust extractor fails to stop when a connected device is switched off.	• Dust-generating device not plugged into the unit power outlet on the Silent TS.	• Connect the dust-generating device to the Silent TS (Chap. 4.4).
	• Wrong operating mode selected.	• Select the automatic mode (Chap. 5.4).
	• Activation threshold too low.	• Adjust the automatic start feature (Chap. 5.5).
No satisfactory setting of the automatic start feature can be adjusted.	• Some older handheld instruments do not generate a sufficiently strong signal for the automatic start feature.	• Switch the dust extractor to the continuous mode.
Extraction briefly starts up when the unit is switched on.	• Activation sensitivity too low.	• Adjust the automatic start feature (Chap. 5.5).
Panels covering the fine particle filter and the dust drawer are difficult to open.	• Lock bolts dirty.	• Clean the lock bolts and lubricate them slightly, if required.

8 Technical data

	2921 0050	2921 1050
Nominal voltage	230 V	120 V
Permissible mains voltage:	220 - 240 V	100 - 120 V
Mains frequency:	50 / 60 Hz	50 / 60 Hz
Suction turbine power: *)	1400 W	1400 W (120 V) 1000 W (100 V)
max. connecting power: *)	2000 W	360 W (120 V) 500 W (100 V)
Total connected power: *)	3400 W	1800 W (120 V) 1500 W (100 V)
Mains input fuse (L, Fig. 4):	2 x 15 A (T)	
LpA **) (at max. volume flow):	56.3 dB(A)	
Ø suction support: - inner - outer	32 m [1.26 inch] 40 m [1.57 inch]	
Volume flow, max.: ***)	3300 l/min (120 V / 230 V) 2900 l/min (100 V)	
Vacuum pressure, max.: ****)	200 hPa [2.9 psi]	
Fine filter: - Filter surface, approx. - Filter quality	0.8 m² [1240 sq inch] M	
Fill volume, dust bag, approx.:	7.5 l [2 US gal]	
Dimensions (Width x Height x Depth):	224 x 591 x 564 mm [8.8 x 23.3 x 22.2 inch]	
Weight (empty), approx.:	26 kg [57 lbs]	

*) Performance values at nominal voltage

**) Noise level according to EN ISO 11202

***) Free blowing turbine at nominal turbine voltage

****) At nominal turbine voltage

9 Warranty

Renfert offers a 3-year warranty on all parts of the equipment, provided that the equipment is properly used.

A 3-year warranty is provided for the suction motor, however with a maximum of 1000 operating hours (motor running time).

Warranty claims may only be made upon presentation of the original sales receipt from the authorized dealer.

Parts that are subject to natural wear and tear (wear parts) and consumables are excluded from the warranty. These parts are marked in the spare parts list.

The warranty is voided in the case of improper use; failure to observe the operating, cleaning, maintenance, and connection instructions; in the case of independent repairs or repairs that are not carried out by an authorized dealer; if spare parts from other manufacturers are used, or in the case of unusual influences or influences not in compliance with the utilization instructions.

Warranty claims shall not extend the warranty period.

10 Disposal information

10.1 Disposing of consumables

Full dust bags and filters, including motor and electronic filters, must be disposed of in accordance with country-specific regulations.

When doing so, personal protective equipment must be worn, depending on the material trapped in the filter.

10.2 Disposing of the equipment

The equipment must be disposed of by a specialist facility. The specialist facility must be informed of any residue in the equipment that is harmful to health.

10.2.1 Disposal information for countries in the EU

To conserve and protect the environment, prevent environmental pollution, and improve the recycling of raw materials, the European Commission has adopted a directive that requires the manufacturer to accept the return of electrical and electronic equipment for proper disposal or recycling.

For this reason, equipment with this symbol may not be disposed of in unsorted domestic waste within the European Union.

Please contact your local authorities for more information on proper disposal.

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declare under our sole responsibility, that the product(s)

Name	Part No.
SILENT TS	2921 0050
SILENT TS2	2930 0050
SILENT TC	2935 0000 / 2935 0000 03
SILENT TC2	2936 0000 / 2936 0000 03
SILENT EC2	2937 0000 / 2937 0000 03
SILENT compact	2934 0000 / 2934 0000 03
SILENT compactCAM	2934 2000 / 2934 2000 03 2934 2500 / 2934 2500 03
SILENT powerCAM EC	2939 0000
SILENT powerCAM TC	2938 0000
iVAC silent	36 2934 2000

is (are) in conformity with the relevant regulatory requirements by compliance with the UK designated standards.

UK legislation

Electromagnetic Compatibility Regulations 2016

Supply of Machinery (Safety) Regulations 2008

The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

Designated standards

EN 61326-1:2013

EN ISO12100:2010

EN 61010-1:2010/A1:2019/AC:2019-04

EN IEC 63000:2018

Tilo Burgbacher,
Engineering Director

Hilzingen, 02.11.2022

The Technical Documentation has been retained by Renfert GmbH

Hochaktuell und ausführlich auf ...
Up to date and in detail at ...
Actualisé et détaillé sous ...
Aggiornato e dettagliato su ...
La máxima actualidad y detalle en ...
Актуально и подробно на ...

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Renfert GmbH, Industriegebiet, 78247 Hilzingen / Germany
Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany