Top spin

No. 1835

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Instruction manual • Mode d'emploi

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ENGLISH

1. Introduction

We are pleased that you have decided to purchase the *Top spin* pin hole drilling unit. This unit sets new standards in function, performance, safety and design.

To ensure that you enjoy using this product for many years please take note of the following instructions.

2. Description / Applications

Top spin is designed to accurately drill holes for the secure fixation of pins and various pin systems within dental technology.

2.1 Proper use

Top spin is a device for precisely drilling positioned fixture holes for pins in plaster dental models or epoxy resin for a range of pin systems.

It is exclusively for use in dental facilities.

2.2 Improper use

The product must not be used for drilling into metals.

The product must not be used for milling. Breakage of the drill presents a risk of injury.

Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the device, increase the risk of serious injury and lead to damage to the environment or the device itself.

The protection provided is no longer assured in the case of improper use and there is a danger of injury.

2.3 Ambient Conditions

(in accordance with DIN EN 61010-1)

The *Top spin* may only be operated:

- Indoors;
- Up to an altitude of 2,000 m [6,500 ft] above sea level;
- At an ambient temperature of 5 - 40°C [41 - 104°F] *);
- At a maximum relative humidity of 80% at 31°C [87.8°F], decreasing linearly to 50% relative humidity at 40°C [104°F] *);
- With a mains power supply where current fluctuations do not exceed 10% of the nominal value;
- Under level 2 contamination conditions;
- Under over voltage category II conditions.

*) Between 5 - 30°C [41 - 86°F], the unit can be operated at a relative humidity of up to 80%. At temperatures between 31 - 40°C [87.8 - 104°F], the humidity must decrease proportionally in order to ensure operational readiness (e.g., at 35°C [95°F] = 65% humidity; at 40°C [104°F] = 50% humidity). The unit may not be operated at temperatures above 40°C [104°F].

3. Areas of Caution

3.1 Employed Symbols

The following symbols are used in these instructions and on the unit:

Danger

This indicates an imminent risk of injury.

Electrical hazard

There is a risk of electrical shock.

Attention

Failure to observe the associated information can result in damage to the unit.

Note

The operating instructions contain useful information to make handling easier.

Laser

Do not look into the laser beam!

For indoor use only.

Note the operating instructions.

3.2 Areas of Caution

- **Top spin** is an electrical appliance with a danger potential. This appliance may be connected / used only by authorized and qualified persons after checking its conformity with respect to the national norms.
- Observe the regulations for prevention of accidents of the Employer's Liability Insurance Association and the following safety precautions when drilling:
- Before switching on always check the voltage on the unit matches the voltage supply in your laboratory.
- Units with faulty connections or other defects should no longer be put into operation.
- Pull the power plug before working on the appliance.
- Never reach into the drill when the unit is connected to the electricity supply. There is a risk of injury if the unit is accidentally activated!
- Do not work with the **Top spin** with long, untied hair or loose / baggy clothing. If anything gets tangled up or drawn in there is a risk of injury.
- The unit is solely designed for drilling plaster or epoxy-resin models.
- Do not touch the drill when in motion.
- The unit should only be activated when the model is placed in position.
- When drilling do not hold the model directly over the laser beam - risk of injury.

Attention!

Laser beam!

Laser class 2

Do not look into the beam!

- It is the responsibility of the operator that national regulations during operation and regarding a repeated safety inspection of electrical equipment are complied with. For Germany these are BGV A3 in relation with VDE 0701-0702.

3.3 Liability Exclusion

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- The product is employed for any purposes other than those cited in the operating instructions;
- The product is altered in any way other than those alterations described in the operating instructions;
- The product is repaired by other than an authorized facility or if any but Renfert OEM parts are employed;
- The product continues to be employed, despite obvious safety faults.

4. Installation / Operation

Place the **Top spin** on a level surface. Two operating positions are possible:

1. Operation from a standing position = 0° setting (Fig. 1).
2. To change position: turn the housing base by (Fig. 2).
3. Operation from a sitting position = 16° setting (Fig. 3).

Check that the voltage in your laboratory matches the voltage printed on the identification plate.

4. Connect the plug - **Top spin** (Fig. 4).
5. Connect the plug - mains (Fig. 5).
6. Set the power switch to I (Fig. 6).
7. Check adjustment of drilling depth. Make a trial boring. Drilling depth can progressively be adjusted (Fig. 7).

The **Top spin** is now ready to use.

5. Operation / Drilling

1. Focus the laser beam on the model die (Fig. 8).
2. Hold the die firm with both hands (Fig. 8).
3. Motor starts automatically when the model table is depressed (Fig. 9).
4. For greater pressure use the *Top spin*'s counter-balances (Fig. 10).
5. Switch off the *Top spin* (Fig. 11).

Attention

Do not use *Top spin* to drill metal!

Do not use lubricants or cutting abrasives!

5.1 Changing the drill bit

1. Switch off the *Top spin* (Fig. 11) and pull the power plug.
2. Lift the work table and click it into position (Fig. 12).
3. Remove the dust bucket (Fig. 13).
4. Set the groove of the drill spindle to locking pin (Fig. 14).
5. Press the locking pin into the groove of the drill spindle (Fig. 15) and hold.
6. Loosen the screw cap by turning anti-clockwise (Fig. 16).
7. Remove the drill bit (Fig. 17) and replace it with a new one.
8. Hand-tighten the screw cap in a clockwise direction (Fig. 18).
9. Replace the dust bucket (Fig. 13).
10. Lower the work table (Fig. 19).

Attention

Use only original Renfert drill bits!

5.2 Changing of the fuse

1. Switch off the *Top spin* (Fig. 11) and pull the power plug.
2. Detach clips at the top and at the bottom (Fig. 22a).
3. Pull out holder together with fuses (Fig. 22a).
4. Change defect fuse (Fig. 23), for fuse type see „Technical Data“.
5. Put back the holder together with fuses (Fig. 22b).

6. Cleaning / Maintenance

Only clean *Top spin* in a dry state.

1. Switch off the *Top spin* (Fig. 11).
2. Lift the work table and click it into position (Fig. 12).
3. Remove the dust bucket (Fig. 13).
4. Empty the dust bucket and replace it.
5. Lower the work table (Fig. 19).

6.1 Changing or cleaning the drill spindle

1. Remove the drill bit (see 5.1, points 1 to 7).
2. Loosen the drill spindle by turning clockwise (Fig. 20).
3. Clean or change the drill spindle.
4. Hand-tighten the drill spindle in an anti-clockwise direction (Fig. 21).
5. Insert the drill bit (see 5.1, points 8 to 10).

7. Spare Parts

For the order numbers of consumables and spare parts please refer to the spare parts list at the end of this manual.

8. Guarantee

When used properly Renfert guarantees all parts of the *Top spin* for 3 years. In case of any claims for warranty original dealers' invoice is required. Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the guarantee. These parts are marked in the spare part list.

The guarantee is void for damage caused by misuse, disregard of the instructions for use, cleaning, maintenance and connection, for improper self-repaired units or repair by untrained personnel, for use of spare parts from other manufacturers and on influences outside the normal practice or unauthorized in the instruction manual. Servicing made under guarantee does not extend the guarantee.

9. Disposal instructions for countries in the EU

To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic units for proper disposal or recycling. Within the European Union units with this symbol should not therefore be disposed of in unsorted domestic waste:

For more information regarding proper disposal please apply at your local authorities.

10. Technical Data

Voltage:	230 V, 50 Hz 115 V, 60 Hz
Power consumption:	150 W
Fuse:	
1835-0000 (230V/50Hz):	M1A (1A, medium slow)
1835-4000 (115V/60Hz):	M2A (2A, medium slow)
Revolutions:	16000 1/min.
Dimensions (BxHxD):	166 x 322 x 149 mm [6.54 x 12.69 x 5.87 inch]
Weight:	5 kg
Drill type:	HM-drill with shaft-Ø 3 mm [0,12 inch]
Table hub:	0 - 12 mm [0 - 0.47 inch], infinitely variable
Level of noise pressure according to DIN EN ISO 11202:	LpA < 70dB(A)
Laser:	
Laser-Class:	2
Wavelength:	630-680nm
Power output P0:	max. < 1mW

11. Parts Supplied

- 1 *Top spin*
- 1 stepped drill bit (see drawing 1)
- 1 instruction manual with spare parts list
- 1 mains cable
- 1 set of hardware

12. Accessories


5010-0001 1 stepped drill bit

(see drawing 1)


367-0000 3 Smart-Pin drill bits

(see drawing 2)

Drawing 1

Drawing 2

13. Troubleshooting Guide

Fault	Cause	Remedy
<ul style="list-style-type: none">The unit does not switch on.	<ul style="list-style-type: none">The electricity supply is not properly connected.The electrical safety device is faulty.	<ul style="list-style-type: none">Check the electricity supply.Check the electrical safety device.
<ul style="list-style-type: none">The drill holes are too big or not round.	<ul style="list-style-type: none">The drill bit is damaged or worn.The drill jig is not tight enough.Chuck dirty.	<ul style="list-style-type: none">Replace the drill bits, clean the drill jig.Tighten the jig by hand.Cleaning chuck.
<ul style="list-style-type: none">The drilling is not effective.	<ul style="list-style-type: none">The drill bit is worn.	<ul style="list-style-type: none">Replace the drill bit.

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Renfert GmbH • Industriegebiet • 78247 Hilzingen/Germany