making work easy

FAQ - Deutsch

FAQ - English

FAQ - Français

FAQ - Italiano

FAQ - Español

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FAQ - Türkçe

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FAQ - English

For information on operating elements, see [Fig. 1](#) and [Fig. 2](#) on page 2 of this document.

For chapter indications, refer to the complete instruction manual on the internet.

EN

Problem	Possible cause	Corrective action
No air and no abrasive.	- Blasting nozzle blocked at the hand-piece. - Foot switch connecting lines kinked. - Tank selection switch in an intermediate position.	- Unscrew the nozzle from the hand-piece and blow it out. - With the nozzle removed, blow out the blasting hose by activating the foot pedal. - Inspect / correct the connection lines. - Turn the switch until you feel it click into place.
Only air, no or very little abrasive.	The dosage nozzle in the mixing chamber is blocked due to contaminated or moist sand.	Cleaning the dosing nozzle, see chapter 6.4.
Only air, no or very little abrasive.	- Mixing chamber filter combination blocked. - Abrasive contains impurities or is moist.	- Change filter cartridge (22 / Fig. 1) (wrench SW14). - Replace the abrasive.
Unit blasts continuously.	- Foot switch connected incorrectly. - Valve in foot switch is defective.	- Inspect / correct the foot switch connection. - Replace foot switch.
Lamp is not alight.	- Power cord not plugged in. - Illumination defect.	- Plug the power cord into the wall outlet. - Change light using the service or see spare part list.
Compressed air leaking from the foot switch.	Pneumatic hoses reversed.	Reverse the connections.
Foot switch fails to switch on.	- No pressure at the compressed air connection. - Pneumatic hoses are kinked.	- Inspect the compressed air lines and / or the compressor. - Reposition the hoses.
Foot switch hoses cannot be disconnected from the unit.	Push-in connection is not released (Fig. 2).	Press the push-in connection ring on the hose side firmly while pulling on the hoses. Please use the release tools supplied (24).
Quick-stop not functioning (the nozzle blasts continuously).	- Pneumatic connections are loose. - Sandblasting hose worn. - The connecting pressure at the external air connection (15, Fig. 1) is too small. - Outer pressure regulator on the maintenance unit is altered (fixing screw is loosened).	- Inspect the pneumatic connections for correct seating. - Replace the sandblasting hose. - Increase the external connecting pressure at a min. of 5 bar. The external connecting pressure is displayed on the manometer (8, Fig. 1). - Turn the inner pressure regulator clockwise until stop. Activate foot pedal and set the pressure regulator on the maintenance unit at 5 bar (see manometer).
Sound of rushing air when the unit is operating.	- A small amount of leakage is always present at the control valves, thus allowing (a minimum) of air to escape. - Continue to turn the tank selection switch whilst activating the foot switch at the same time.	- This is not a fault, but instead is the result of a pressure release opening which must not be blocked. - Only turn the tank selection switch when the foot switch is not activated.

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