making work easy

de  
en

# SIMPLEX sliceware

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Made in Germany

21-7140 24082022

making work easy

# SIMPLEX sliceware

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Made in Germany

# Inhalt

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# 1 Introduction

SIMPLEX sliceware was developed with the goal of simplifying and perfecting the 3D printing process with filaments.

To print a 3D object, the SIMPLEX 3D-filament-printer requires tool path information (G-Codes).

With SIMPLEX sliceware, this G-Code is generated from 3D object files (e.g., \*.stl or \*.obj).

SIMPLEX sliceware features two modes, Renfert mode and Expert mode.

## 1.1 Symbols used

Symbols with the following meanings are used in these instructions or on the equipment:

### Caution

**Risk of damage to the equipment if the instruction is not followed.**

### Notice

*Provides the operator with useful information that eases use.*

#### ► List, particular attention should be paid

- List
- Subordinate list

⇒ Instructions / appropriate action / input / operational sequence:

You are asked to carry out the specified action in the order indicated.

- ◆ Result of an action / response by the equipment / response by the program:

The equipment or program responds as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur.

# 2 Safety

Before using SIMPLEX sliceware, read this user manual carefully and keep it for future reference. Renfert GmbH reserves the right to modify or revise this user manual at its own discretion at any time.

## 2.1 Disclaimer

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- **The product is employed for any purposes other than those specified in the operating instructions.**
- **The product is altered in any way other than those alterations described in the operating instructions.**



## 3 Installation

With the help of the Renfert CONNECT app you can download the programs SIMPLEX sliceware and SIMPLEX print and install them on your computer as well as receive firmware updates for your SIMPLEX 3D filament printer







**If you do not have a mobile device on which you can install the Renfert CONNECT app, or if you have problems installing the software, please contact Renfert Service at [Support@renfert.de](mailto:Support@renfert.de).**

- ⇒ Install the Renfert CONNECT app from the App Store (Apple) or the Google Play Store (Android) on your mobile device.
- ⇒ Open the app and log in with an existing account, or register to create a Renfert account for your laboratory / dental office.
- ⇒ Go to “My Equipment”.

- ⇒ Select “Add device”.
- ⇒ Select “SIMPLEX” as the device type.
- ⇒ Enter the serial number of your printer.  
The serial number is provided on the nameplate on the rear of the unit.
- ⇒ Select “Add”.

◆ You can now find all the important information about your SIMPLEX 3D filament printer under “My SIMPLEX”.

### 3.1 Installing SIMPLEX sliceware / SIMPLEX print

Under “My SIMPLEX”, you can request the software package to use your SIMPLEX 3D filament printer. Go to “My SIMPLEX”, select “Software Downloads”.

A link is sent to the e-mail address that you stored during registration.

- ⇒ Open the e-mail using an e-mail client on your computer.
- ⇒ Click on the link in the e-mail.
- ⇒ Download the installation program to your computer.
- ⇒ Select Windows 7 or a higher version.
- ⇒ Open the installation program and follow the instructions for installing the software.

## 4 Renfert mode

### 4.1 Starting SIMPLEX sliceware


When you start SIMPLEX sliceware, Renfert mode opens automatically.

A virtual print bed that corresponds to the actual print bed of the SIMPLEX 3D-filament-printer is shown together with the last object that you have edited.

There are two views in Renfert mode:


- Objects view
- G-Code view

### 4.2 Objects view

In the Objects view, you can use the icons on the right, top, and left of the screen to perform different actions.


#### 4.2.1 View / zoom

Icon	Name	Description
	Camera	Changes the view of the print bed.

A window opens where you can select different views of the print bed.

You can also change the view using the NUM block keys:

- ⇒ Num 2: Bottom
- ⇒ Num 3: Back
- ⇒ Num 4: Left
- ⇒ Num 5: Front
- ⇒ Num 6: Right
- ⇒ Num 7: Home
- ⇒ Num 8: Top
- ⇒ Num 9: Extends

In addition, you can use the mouse to change the view and size (zoom).

- ⇒ Scrolling using the mouse wheel: Increases / decreases size
- ⇒ Right mouse button pressed: Rotates the view
- ⇒ Center mouse button (mouse wheel) pressed: Moves



**If you no longer have a good overview when zooming or moving the view using the mouse, you can easily restore the standard view by selecting the Camera icon and "Home" or using key 7 of the NUM block.**

*This function is available in all views of Renfert and Expert mode and is not explained again there.*

en



## 4.2.2 File functions

Icon	Name	Description
	New Proj	Creates a new project.
	Open	Loads an existing SIMPLEX sliceware project (*.cwprj file). Alternatively, this can also be used to add an additional print object.
	Save	Saves the project.
	SaveAs	Saves the project under a new name.

*When saving, do not use special characters or umlaut characters in the file name!*

## 4.2.3 Undoing / redoing actions

Icon	Name	Description
	Undo	The last step is undone.
	Redo	The last step is redone.

## 4.2.4 Selecting print objects


When you move the cursor over a print object, the cursor changes into a hand icon.

⇒ Left-click to select an object:

- ♦ The object is shown in color.
- ♦ The space that the object takes up is shown.

To select multiple objects, click them with the Ctrl key pressed.

## 4.2.5 Adding / removing objects

Icon	Name	Description
	Add	Adds a print object.
	Remove	Removes a print object.

### Adding an object:

⇒ Click "Add".

- ♦ A dialog opens where you can open a file.

⇒ Select the directory in which the 3D data (\*.stl or \*.obj) is stored.

⇒ Select the required object.

- ♦ The object is placed in the center of the print bed.

Alternatively, you can also use the "Load" file function to place an object on the print bed or do so directly from Explorer using drag-and-drop.

You can place multiple objects on the print bed.

**All selected print objects must be printed using the same filament!**

### Removing an object:

⇒ Select the object to be removed.

⇒ Click "Remove".

Alternatively, you can also delete a selected object using the delete key on your keyboard.

If you have selected multiple objects, they are all removed.

## 4.2.6 Moving objects

Icon	Name	Description
	Move	The object can be moved to any position.

⇒ Click "Move".

- ◆ An information window opens in which the move action is shown for the 3 spatial directions X, Y, and Z.

⇒ Select the object.

- ◆ Three arrows are shown, for the three possible directions of movement.

⇒ Place the mouse over one of the arrows.

- ◆ The arrow is shown in white.

⇒ Left-click the arrow and keep the mouse button pressed.

- ◆ Moving the mouse allows you to move the object in the direction of the arrow.

Alternatively, you can move a selected object freely in the X / Y plane with the left mouse button pressed.

In addition, the object can also be moved by entering a number for one of the directions in the information window.

If you have selected multiple objects, they are all moved together.


Fig. 1


Fig. 2

**Objects are shown in red if you move them outside of the print range.**

## 4.2.7 Rotating objects

Icon	Name	Description
	Rotate	Rotates objects

⇒ Click "Rotate". An information window opens in which the angle of rotation is shown for the 3 spatial directions X, Y, and Z (see Fig. 1).

⇒ Select the object.

- ◆ Three arrows are shown, for the three possible axes of rotation (see Fig. 2).

⇒ Place the mouse over one of the arrows.

- ◆ The arrow is shown in white.

⇒ Left-click the arrow and keep the mouse button pressed.

- ◆ Moving the mouse allows you to rotate the object around this axis.

In addition, the object can also be rotated by entering a number for an angle of rotation in the information window.

If you have selected multiple objects, they are all rotated together around a shared axis.

## 4.2.8 Dropping objects

Icon	Name	Description
	Drop	The object is dropped onto the print bed.

⇒ Click "Drop".

♦ An information window opens with information on dropping an object (see Fig. 1).

⇒ Place the mouse over the object that you would like to drop.

♦ A small blue arrow appears and shows the layer with which the object is dropped onto the print bed.

⇒ Left-click to perform a drop onto the layer that is shown.

**Before dropping an object, it is helpful to rotate it in such a way that the layer onto which the object is to be dropped can be easily clicked, or to select the corresponding view.**

## 4.2.9 Duplicating objects

There are two ways to duplicate an object:

⇒ Select the object.

⇒ Press and hold the ALT key.

⇒ Move the mouse.

♦ A copy of the object is generated that can be moved and placed as required.

ALTERNATIVELY

⇒ Click "Move".


⇒ Select the object.

⇒ In the information window, click "Multiply".

⇒ Specify how many instances of the same object should be available after multiplication has been completed. If you would like one copy, for example, enter "2".

⇒ Click "OK".

♦ A copy of the object is generated and placed on the print bed.

**If you select several objects, all the selected objects are multiplied.**


## 4.3 Slicing / Creating a G-Code / switching to the G-Code view

Icon	Name	Description
	Slice...	Starts the slicing process / switches to the G-Code view

As soon as a prototype of an object is ready, the tool path information (G-Code) must be generated.

This process is known as slicing.

⇒ Click "Slice...".

⇒ Select the profile of the type of model to be printed by clicking on the corresponding image.


♦ The G-Code is automatically generated.

♦ The G-Code view is automatically shown.

The G-Code view shows how the SIMPLEX 3D-filament-printer will print the object.

## 4.4 G-Code view


In the G-Code view, you can use the icons on the right, top, and left of the screen to perform different actions.

Additionally, information is shown at the top right of the screen about creation of the object in the SIMPLEX 3D-filament-printer.

### 4.4.1 Switching to the Objects view

Icon	Name	Description
	Objects view	Switches back to the Objects view

### 4.4.2 Animation

Icon	Name	Description
	Animate	Switches animation of the print process on / off


During animation, the movement of the tool head is shown.

### 4.4.3 Printing information

Icon	Name	Description
	Stats	Shows / hides print information.

Activates / deactivates the information window for printing the object.

### 4.4.4 Clearing the G-Code

Icon	Name	Description
	Clear	Clears the G-Code


The G-Code generated is cleared.

**Clearing the G-Code cannot be undone.**

To create a new G-Code, you must switch back to the Objects view and perform slicing again.

#### 4.4.5 Loading / saving a G-Code

Icon	Name	Description
	Load	Loads previously saved G-Codes. The object appears again on the print bed.
	Save	Saves the current G-Code.

**When saving, do not use special characters or umlaut characters in the file name!**

#### 4.4.6 Printing using a PC

Before the G-Code can be sent to SIMPLEX print and printed, it must be saved.  
When saving, do not use special characters or umlaut characters in the file name!

Icon	Name	Description
	Print	The G-Code is sent directly to a connected SIMPLEX 3D-filament-printer using the SIMPLEX print program, and the object is printed.

When you click "Print", the SIMPLEX print program is automatically started, and the G-Code is sent to SIMPLEX print.

You can then start the print process in the SIMPLEX print program (see Chapter 6).


#### 4.4.7 Printing using a USB flash drive

Icon	Name	Description
	Pendrive	Saves the print data to a USB flash drive. Do not use umlaut characters or special characters in the file name!

If a SIMPLEX 3D-filament-printer is not connected to the PC, you can also save the G-Code to a USB flash drive. The USB flash drive can then be inserted into the SIMPLEX 3D-filament-printer and the object printed.

The "Pendrive" button is only available when a USB flash drive is connected to the PC.

### 4.5 Options

Icon	Name	Description
	Options	Changes settings

Using the Options window, you can configure different settings in various tabs to adapt the program to your way of working.

#### 4.5.1 General

Changes the language, settings to verify whether updates are available, as well as settings for working with files.

**Changing the language becomes effective the next time you start the program.**

#### 4.5.2 Graphics

Display settings in the Objects view and G-Code view modes.

#### 4.5.3 Keyboard

Key and keyboard assignments for different commands.

#### 4.5.4 Colors

Settings for the colors in which different segment types, extrusion and feed rate, and tool heads should be shown. These colors are used in Expert mode.

#### 4.5.5 About

Information about changes in the different program versions.  
This information is only available in English.

## 5 Expert mode

Among other things, Expert mode offers a variety of settings options when creating a G-Code using the "Slice" function (G-Code view). These settings affect print quality.

Only use Expert mode if you are familiar with the settings options.

**We do not accept any liability for faulty print results when using Expert mode.**

Icon	Name	Description
	To Expert	Switches to Expert mode

In Renfert mode, you can switch to Expert mode in the Objects view using the "To Expert" icon.

You can switch back to Renfert mode from Expert mode using the "To Renfert" icon:

Icon	Name	Description
	To Renfert	Switches to Renfert mode

There are three views in Expert mode:

- Objects view
- Support view
- G-Code view

To switch between the views, click the Object, Support and Slice, or G-Code icons, or press numbers 1, 2 or 3 on the keyboard.

In the "Objects" and "Support" views, the "Undo - Redo" option is unlimited. The "Undo - Redo" function is not available in the G-Code view.

### 5.1 Objects view

The Objects view in Expert mode also contains the same functions as the Objects view in Renfert mode. For details, see Chapter 4.2.

Additional commands include:

- Select
- Scale
- Bisect
- Split
- Autodrop

#### 5.1.1 Select

Icon	Name	Description
	Select	Used for giving individual objects a name and assigning a color.

⇒ Click "Select".

- ♦ An information window is shown at the bottom of the screen.

Individual objects can be given a name and assigned a color here.

#### 5.1.2 Bisecting

Icon	Name	Description
	Bisect	Cuts a plane through the object.

**Before bisecting, it can be helpful to zoom out of the view using the mouse wheel to get a better overview of the position of the cutting plane.**

⇒ Click "Bisect".

- ♦ A cutting plane is shown that can be moved and rotated by clicking on the frame or on the arrow.
- ♦ An information window is shown at the bottom of the screen.

⇒ Click on "Halve" in the information window.

- ♦ The object is split into two parts at the cutting plane, each of which can be moved separately.

### 5.1.3 Split

Icon	Name	Description
	Split	Splits models comprised of multiple parts into several individual parts.

### 5.1.4 Autodrop

Icon	Name	Description
	Autodrop	Automatically drops the model onto the print bed before rotation.

Autodrop can be activated and deactivated.


When activated, the model is dropped onto the print bed before rotation and the axes of rotation are selected so that the model is always touching the print bed.

## 5.2 Support view

Models with large overhang parts must be printed with a support structure.

This support structure can be generated, modified, and shown in the Support view.


Icon	Name	Description
	Support	Switches to the Support view.


An information window is shown at the bottom of the screen in the Support view.

Here you can specify the thickness of the individual supports (the default value is 1.8 mm), as well as the minimum slope of the surface for which a support structure is generated.

"Auto-Generate Support" is used to generate and display the support structure.

Individual supports can be removed using CTRL and the left mouse button. You can add your own supports by left-clicking with the mouse. If a support cannot be inserted at the clicked location, this is indicated. The entire support structure can be deleted using "Clear all the supports".


## 5.3 G-Code view

Use of the "Slice" function and switching to the G-Code view is possible both in the Objects view and in the Support view.

Icon	Name	Description
	Slice...	Starts the slicing process / switches to the G-Code view

When using the "Slice" function you can configure settings in a simple or in an Expert dialog.


*Only configure settings if you are familiar with the effects.*

⇒ In the dialog, click the "Slice!" button.

- ♦ The G-Code is automatically generated.
- ♦ The G-Code view is automatically shown.

The G-Code view shows how the SIMPLEX 3D-filament-printer will print the object.




In the G-Code view, you can use the icons on the right, top, and left of the screen to perform different actions.


### 5.3.1 Switching information windows on/off

Icon	Name	Description
	Legend	Shows / hides the "Legend" information window. The legend refers to the different print process views.
	GCode	Shows / hides the "G-Code" information window.
	Stats	Shows / hides print information.

### 5.3.2 Print process views

Icon	Name	Description
	Type	Colored display of the different segment types. The meaning of the color is provided in the "Legend" information window.
	Printhead	Shows the tool head. Only relevant for printers with multiple tool heads.
	Extrusion	Colored display of the extrusion rate. The meaning of the color is provided in the "Legend" information window.
	Feed rate	Colored display of the feed rate. The meaning of the color is provided in the "Legend" information window.
	Travel	Displays the path taken by the tool head during travel.

### 5.3.3 Display quality

Icon	Name	Description
	Quality	Switches the quality of the screen display of the rendering. Depends on the graphic card used in the computer.

### 5.3.4 Partial view / section

Icon	Name	Description
	Section	A section allows the inside of the printed object to be viewed. The cutting planes are handled in the same way as in the Objects view (see Chapter 5.1.2).

### 5.3.5 Generating a G-Code again

Icon	Name	Description
	Slice...	Generates a G-Code again.

### 5.3.6 Clearing the G-Code

Icon	Name	Description
	Clear	Clears the G-Code. The G-Code generated is cleared.

**Clearing the G-Code cannot be undone.**

To create a new G-Code, you can perform slicing again.

### 5.3.7 Loading / saving a G-Code

Icon	Name	Description
	Load	Loads previously saved G-Codes. The object appears again on the print bed.
	Save	Saves the current G-Code.

**When saving, do not use special characters or umlaut characters in the file name!**

### 5.3.8 Printing using a PC

Before the G-Code can be sent to SIMPLEX print and printed, it must be saved. When saving, do not use special characters or umlaut characters in the file name!

Icon	Name	Description
	Print	The G-Code is sent directly to a connected SIMPLEX 3D-filament-printer using the SIMPLEX print program, and the object is printed.

When you click "Print", the SIMPLEX print program is automatically started, and the G-Code is sent to SIMPLEX print.


You can then start the print process in the SIMPLEX print program (see Chapter 6).

### 5.3.9 Printing using a USB flash drive

Icon	Name	Description
	Pendrive	Saves the print data to a USB flash drive. Do not use umlaut characters or special characters in the file name!

If a SIMPLEX 3D-filament-printer is not connected to the PC, you can also save the G-Code to a USB flash drive. The USB flash drive can then be inserted into the SIMPLEX 3D-filament-printer and the object printed. The "Pendrive" button is only available when a USB flash drive is connected to the PC.

## 5.4 Options

Icon	Name	Description
	Options	Changes settings

In addition to the settings that are also available in Renfert mode, Expert mode also includes a tab with printer settings.

### 5.4.1 Printer

In Expert mode, you can select the printer used and modify the printer parameters.

## 6 SIMPLEX print

If the SIMPLEX 3D-filament-printer is connected directly to a computer using the USB cable, objects can be printed with SIMPLEX sliceware using SIMPLEX print.

**The print process using SIMPLEX print cannot be followed on the display of the SIMPLEX 3D-filament-printer.**

**For this reason, the alternative is available of printing using a USB flash drive.**

### 6.1 Connecting the printer via USB cable

- ⇒ Connect the USB cable to the USB outlet of the SIMPLEX (Fig. 3).
- ⇒ Connect the cable to a USB port on your computer.
- ⇒ Switch on the printer.


Fig. 3

### 6.2 Automatically starting SIMPLEX print

When you click the "Print" command in SIMPLEX sliceware, SIMPLEX print is automatically started.


Fig. 4

- ⇒ Close the error message with "OK".
- ⇒ Click "Printers" and "Detect SIMPLEX 3D printers" to activate the connected SIMPLEX 3D-filament-printer (alternatively, press CTRL+R).
  - ♦ The SIMPLEX is shown in the status bar (1, Fig. 4)
- ⇒ Select and load the G-Code by clicking the load button (3, Fig. 4).
  - ♦ The corresponding file name is shown in the status bar (2, Fig. 4)

You can use the control keys at the bottom right of the program window to remotely control the SIMPLEX 3D printer.

The following actions are possible:

- Adjust the print bed to a specified height.
- Adjust the print bed to the home position.
- Adjust the tool head to the home position.
- Switch off the motors.
- Manual extrusion


Fig. 5

## 6.3 Print process using SIMPLEX print

**Before starting the print process, check that the correct filament has been placed in the SIMPLEX!**

The print process is started in SIMPLEX print by clicking the start / pause button (Fig. 6).

The extruder and print bed are heated to the temperatures specified in the G-Code.

Printing begins as soon as the required temperatures have been reached.

You can pause printing using the pause button and cancel it using the stop button.


Fig. 6

You can use the control keys at the bottom right of the program window to control the print process.

The following actions are possible:

- Modify the printing speed
- Adjust the quantity to be extruded
- Switch the cover fan on / off
- Modify the extruder temperature
- Modify the temperature of the print bed


Fig. 7

**For printing, a larger amount of G-Code data is always sent to the SIMPLEX.**

**There is no response to the pause / stop button or control key commands until the G-Code that has already been sent has been fully processed.**

**In some cases, this can take several minutes.**

## 6.4 Starting SIMPLEX print manually

Alternatively, you can also start and use SIMPLEX print manually using the start menu of your computer.


Fig. 8

It is also the case here that you must activate the connected SIMPLEX 3D-filament-printer and load the G-Code (see description in Chapter 6.2).

**Hochaktuell und ausführlich auf ...**  
**Up to date and in detail at ...**  
**Actualisé et détaillé sous ...**  
**Aggiornato e dettagliato su ...**  
**La máxima actualidad y detalle en ...**  
**Актуально и подробно на ...**

**[www.renfert.com](http://www.renfert.com)**

Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany  
Tel.: +49 7731 82 08-0 • Fax: +49 7731 82 08-70  
[www.renfert.com](http://www.renfert.com) • [info@renfert.com](mailto:info@renfert.com)

Renfert USA • 3718 Illinois Avenue • St. Charles IL 60174/USA  
Tel.: +1 6307 62 18 03 • Fax: +1 6307 62 97 87  
[www.renfert.com](http://www.renfert.com) • [info@renfertusa.com](mailto:info@renfertusa.com)  
USA: Free call 800 336 7422