Fig. 2

Sicherheitsbrenner eco

Quick Start Guide

Made in Germany

Einleitung

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EN Introduction

This quick start guide contains guidelines for use of the unit and relevant safety information. For detailed information on the unit, please see the instructions for use which can be downloaded from our Internet page under www.renfert.com.

Symbols used

Danger

Immediate danger. Observe the accompanying documentation!

Warning

By not following the instructions the burner could become damaged.

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Liquid gas: Burner No. 934-0100

Natural gas E: Burner No. 935-0100

The burner has been certified and released in conformance with the DVGW.

Product identification number: DG-2411CO 0243

Hazard Statements and Warnings

If the burner is not used according to the supplied instructions for use, then the intended safety feature can no longer be guaranteed.

- Always pay attention to the accident prevention regulations stipulated by the respective trade organisations! (DVGW regulations).
- The burner should only be operated under constant supervision.
- Even when the gas valve (3) on the burner is closed, the gas flow is not securely shut off.
- After use, the gas supply should be turned off immediately via the gas shut-off valve at the laboratory workbench (gas isolation valve) or by separating the gas connection safety valve.
- In classrooms, laboratory burners may only be used after turning on all the lined up shut-off systems (central room shut-off, group shut-off system, shut-off valve).
- Should a strong smell of gas be detected, close the gas shut-off valve at the laboratory workbench immediately and initiate the corresponding safety measures (see DVGW law and valid national regulations).

Product Description

Component groups and functioning elements

See Fig. 1 (on cover page)

- | | | | |
|---|-------------|----|---------------------------------------|
| 1 | Mixing tube | 6 | Actuation knob for ignition safeguard |
| 2 | Sensor | 7 | Collection pan |
| 3 | Gas valve | 8 | Sensor holder |
| 4 | Knurled nut | 9 | Clamping screw |
| 5 | Hose spout | 10 | Allen key |

EN

Initial Operation

The initial operation of the burner may only be carried out by qualified personnel according to DVGW stipulations or valid national regulations. Please observe the following:

- ⇒ Check if the type of gas and the connection pressure of the burner correspond to the gas supply.
- ⇒ Check the fit of the connection pipe and that is not damaged.
- ⇒ Close the gas shut-off valve at the laboratory workbench.
- ⇒ Slide the certified gas hose (according to DIN 30664 or DVGW G 5501-P) up to the last groove of the hose connection (5) on the burner. If necessary secure the hose with a clamp or coupling to prevent it from sliding off.

Operation

Igniting the flame

Caution / risk of burning!

When igniting the burner, turn the mixing tube away from the user and other objects. Keep your distance! The flame could be big!

- ⇒ Open the gas shut-off valve at the laboratory workbench (gas isolation valve).
- ⇒ Open the gas valve (3) to approx. 1 turn.
- ⇒ Press the actuation knob (6) in, so that you can hear the gas flowing.
- ⇒ Ignite the flame immediately!
- ⇒ After igniting the flame, keep the actuation knob (6) pressed for another 10 seconds until the flame is burning independently. If the flame goes out, repeat the procedure.

Adjusting the flame

- Soft flame: turn the knurled nut (4) in an upward direction.
Strong flame: turn the knurled nut (4) in a downward direction.

Cleaning / Maintenance

Cleaning

Caution:

Never clean the burner by boiling or steam cleaning! Never use acetone when cleaning!

Do not blow compressed air into the hose connection or the mixing tube on the burner! This could damage the safety valve.

Only the following parts may be cleaned with boiling water or steam cleaner once they have been dismantled:

- ▶ Mixing tube (1)
- ▶ Sensor holder (8)
- ▶ Collection pan (7)

Caution when cleaning:

- Only clean the burner once it is completely cool.
- Before cleaning, close the gas shut-off valve at the laboratory work-bench (gas isolation valve).
- Do not use chemicals for cleaning.
- Never drill or unscrew the burner nozzle.

Cleaning the mixing tube

Disassembly:

- ⇒ Using an Allen key (10), loosen the clamping screw (9) and sensor holder (8).
- ⇒ Unscrew the mixing tube (1).

Cleaning:

- ⇒ Clean the mixing tube (1) with boiling water or the steam cleaner.

Assembly:

See Fig. 2 (on cover page)

- ⇒ The re-assembly is carried out in reverse order.
- ⇒ Position the sensor so that the upper edge of the sensor's brass sleeve and the upper edge of the mixing tube are in line (2).
- ⇒ Carry out a function test.

Cleaning the burner nozzle

If the burner nozzle is blocked, it can be dismantled from the mixing tube and carefully cleaned with a soft cloth.

If cleaning in this manner is not possible, then the burner nozzle can only be cleaned and changed by an authorized specialist.

Maintenance

The burner is maintenance-free.

REACH

Information on REACH and SVHC is available on our website
www.renfert.com, in the Support area.

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Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany