making work easy

Vario basic

FAQ - frequently asked questions

Made in Germany

FAQ - Deutsch

FAQ - English

FAQ - Français

FAQ - Italiano

FAQ - Español

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FAQ - Türkçe

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FAQ - English

For information on operating elements, see [Fig. 1](#) on page 2 of this document.

For chapter indications, refer to the complete instruction manual on the internet.

| Problem | Possible cause | Corrective action |
|---|--|---|
| No air and no abrasive | <ul style="list-style-type: none"> Nozzle blocked. Nozzle blocked at grip and/or blasting hose blocked. Kink(s) in the foot switch connector hose. Switch set to intermediate position. Pressure set too low. | <ul style="list-style-type: none"> Open the blasting nozzle unit and clean the nozzle (Fig. 23, 24, 25). Unscrew the nozzle from the handle and blow it out with compressed air. Without reconnecting the nozzle, activate the foot switch to blow compressed air through the hose. Inspect and correct the connecting hoses. Turn the switch until you feel it click into place. Open the pressure regulator further (recommended pressure: 4 - 5 bar). |
| Air only and no or very little abrasive | <ul style="list-style-type: none"> Dirty or moist abrasive. Insufficient abrasive in the unit. Intake hose blocked. Problems with the additional tank. | <ul style="list-style-type: none"> Replace the abrasive. Add abrasive. Check that the intake hose is not kinked / pinched / blocked. The intake hose must be securely attached to the unit base (Fig. 7). Refer to "Additional tank - Mounting instructions", No. 21-9774. |
| Reduced blasting force | <ul style="list-style-type: none"> Function selector switch not properly set. Nozzle worn. Abrasive used up. Abrasive feed hose is blocked / damaged / kinked. Insufficient abrasive in the reservoir. | <ul style="list-style-type: none"> Make sure the selector switch clicks into place. Replace the nozzle. Replace the abrasive. Drain the abrasive out and inspect / replace the feed hose. Add abrasive. |
| Excessive abrasive consumption | <ul style="list-style-type: none"> Extraction power set too high. | <ul style="list-style-type: none"> Reduce the extraction power. |
| Blasting chamber fills with dust. | <ul style="list-style-type: none"> Insufficient extraction power. Extractor dust bag full. | <ul style="list-style-type: none"> The suction performance of a weak extraction unit can be increased by closing the boreholes on the back side of the unit (Fig. 14, Pos. b) with the supplied magnetic foil. Clean the extractor / replace the dust bag. |
| Abrasive spurts in pulses. | <ul style="list-style-type: none"> Nozzle connectors not vertical. Abrasive used up or too dirty. Kink in the suction hose. | <ul style="list-style-type: none"> Make sure the nozzle is vertical. Replace the abrasive. Inspect the feed hose. |
| Unit continuously blasting and stops when foot switch is activated. | <ul style="list-style-type: none"> Foot switch connected incorrectly. | <ul style="list-style-type: none"> Inspect and correct the foot switch connection. |
| Foot switch fails to switch on. | <ul style="list-style-type: none"> No pressure at the compressed air connection. Kink(s) in the pneumatic hoses. | <ul style="list-style-type: none"> Inspect the compressed air lines and compressor. Realign the hoses. |
| Air leaking from the foot switch. | <ul style="list-style-type: none"> Pneumatic hoses connected incorrectly. | <ul style="list-style-type: none"> Replace the connectors. |

| | | |
|-----------------------------|--|---|
| Lamp fails to go on. | <ul style="list-style-type: none"> • Power cord not connected. • Illumination defect. <p>The device with fluorescent tubes</p> <ul style="list-style-type: none"> • The fluorescent tube is not situated correctly in the socket. • The fluorescent tube is defect. | <ul style="list-style-type: none"> • Plug the power cord into the wall outlet. • Change light using the service or see spare part list. <ul style="list-style-type: none"> • Check that the fluorescent tube is fitted correctly (see chapter 7.7.2). • Changing the fluorescent tubes (see chapter 7.7.2). |
|-----------------------------|--|---|

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