making work easy

FAQ - Deutsch

FAQ - English

FAQ - Français

FAQ - Italiano

FAQ - Español

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FAQ - Türkçe

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FAQ - English

For information on operating elements, see [Fig. 1](#) and [Fig. 2](#) on page 2 of this document.

For chapter indications, refer to the complete instruction manual on the internet.

EN

Trouble	Reason	Solution
When the device is switched on, there is a loud vibration noise for approx. 8 seconds.	The suction unit performs an automatic filter cleaning sequence.	This is functional and requires no intervention.
An acoustic signal is emitted, the suction unit is switched off and a filter cleaning is performed.	- The internal flow velocity limit was not reached. - After 8 hours of operation (turbine air time) without switching off, a filter cleaning is carried out.	- Repeat work after cleaning is finished. - Switch off devices at the end of each working day at the on / off switch (3. Fig. 1) ausschalten.
After switching on, the display shows "Dust drawer" (23. Fig. 2) and a 3 x acoustic signal is heard.	- The time interval for emptying the dust drawer has expired and the dust drawer has not yet been emptied. - The emptying of the dust drawer has not yet been confirmed.	- Empty dust drawer and confirm by pressing the Enter key (30. Fig. 2). - Confirm emptying the dust drawer by pressing the Enter key (30. Fig. 2).
"Cal" is shown in the display.	- The automatic switch-on function cannot be set. At the selected rotational speed on the handpiece, the difference between stand-by current and operating current is too low. - The automatic switch-on function cannot be set.	- Adjust the switch-on threshold to a higher speed on the handpiece. - Operate the device manually (see chapter 5.2).
"Air" is shown in the display.	During self-diagnosis a lack of compressed air was detected.	Check the air pressure connection. When compressed air is available again, clear error message pressing the menu key (20).
"HOT" is shown in the display.	The electronics have become overheated.	- Switch off the device and let it cool down. - Provide adequate cooling, e.g. - Chapter 4.2. - Use external exhaust air duct (see chapter 4.7). - Change the fine filter (see accessories as well as chapter 6.3).
"ELEC" is shown in the display.	- Leakage: Dust drawer not airtight. - Leakage: Second channel open although it should be closed. - Turbine does not run even though electronics have been switched on. - Faults in electronics.	- Check the correct fit of the dust drawer (see chapter 6.2). - Check compressed air pinch valves. - Clean the air filter (see chapter 6.5). - Check that the pinch valve in the unused channel is closed. - Check whether the turbine cable is properly plugged in (especially after changing the turbine). - Contact Renfert / Service.
"FILt" is shown in the display.	The fine filter is so heavily clogged that the filter cleaning does not show sufficient effect.	Change the fine filter (see accessories as well as chapter 6.3).
Dust drawer over full.	Time interval to the "FULL" display of the dust drawer too high.	Set a shorter time interval (see chapter 5.4.2).
The signal for emptying the dust drawer appears, even though it is not yet full.	The set time interval to the "FULL" dust drawer is too low.	Set a longer time interval (see chapter 5.4.2).

Trouble	Reason	Solution
The suction performance is insufficient.	• Suction performance is set too low. • Blockage in the Suction tube. • High loss of suction due to unfavorable tube position. • Leakage: Dust drawer not airtight. • Leakage: Second channel open although it should be closed. • Fine filter blocked	• Select a higher suction level. • Check suction tube. • Read instructions on connection to the extraction point (see chapter 4.5). • Check the correct fit of the dust drawer (see chapter 6.2). • Check pinch valves in compressed air. • Check that the pinch valve in the unused channel is closed. • Switch the device off and on again, so that a filter cleaning is carried out. • Change fine filter (see chapter 6.3) (if filter cleaning does not improve the performance).
Electrical equipment connected to the appliance socket cannot be operated.	• Suction not switched on. • The device fuse has triggered because the power of the connected electrical device is too high.	• Switch suction on. • Observe the maximum connection capacity of the appliance socket. • Connect the electrical appliance to a separate socket and operate the extraction unit manually (see chapter 5.2).
Suction starts even though the electrical appliance connected to the socket is not being used.	• The automatic switch-on function has been set too low.	• Setting the automatic switch-on (see chapter 5.4.3).
Suction does not stop when the electrical appliance is switched off.	• The automatic switch-on function has been set too low.	• Setting the automatic switch-on (see chapter 5.4.3).
Suction does not start even though the electrical appliance connected to the socket is in used.	• The automatic switch-on function has been set too high.	• Setting the automatic switch-on (see chapter 5.4.3).
Fuse triggers when a connected electrical device is switched on.	• Power consumption of the connected device too high.	• Observe the max. connection value (see chapter 8).
Suction stops suddenly in continuous or automatic mode.	• Suction turbine is overheated. • Suction turbine is defect.	• Switch off the device and leave to cool for at least 60 min. • Check whether suction tube is blocked, eliminate blockage. • Switch the device off and then on again to perform a filter cleaning. • Change fine filter (see chapter 6.3) if filter cleaning does not improve the performance. • Change Suction turbine.
After the suction motor was changed, the suction unit stopped functioning or stops suddenly.	• The plug from the suction motor is not sufficiently clicked in place.	• Connect the suction motor plug correctly. Ensure that the plug is fitted correctly and clicked into place by pulling it.

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