
making work easy

- FAQ - Deutsch
- FAQ - English
- FAQ - Français
- FAQ - Italiano
- FAQ - Español
- FAQ - Türkçe
- FAQ - Polski
- FAQ - Українська
- FAQ - 中文

# SILENT V4

FAQ - frequently asked questions

Made in Germany

21-6669 26092016

# FAQ - Deutsch

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# FAQ - English

For chapter indications, refer to the complete instruction manual on the internet.

Error	Possible Cause	Corrective Action
<b>Extractor stops suddenly for no apparent reason and a error code is indicated in the display.</b>	<ul style="list-style-type: none"> <li>An error has been detected.</li> </ul>	<ul style="list-style-type: none"> <li>See table "Error Codes"</li> <li>If the instructions are of no help to you, or the error code is not listed:                             <ul style="list-style-type: none"> <li>- Make a note of the error code.</li> <li>- Send the device in for repair.</li> <li>- Inform the repair service which error code was noted.</li> </ul> </li> </ul>
<b>The suction performance is inadequate at one suction point.</b>	<ul style="list-style-type: none"> <li>The set suction level is too low.</li> <li>The suction tube is either blocked or has a leak.</li> </ul>	<ul style="list-style-type: none"> <li>Select a higher suction level.</li> <li>Check the suction tube.</li> <li>Please observe the notes in section 4.5.</li> </ul>
<b>Dust bag over full.</b>	<ul style="list-style-type: none"> <li>The set value for the dust drawer fill level is too high.</li> <li>The dust drawer was removed for longer than 15 seconds without emptying.</li> </ul>	<ul style="list-style-type: none"> <li>Set a smaller value (see section 5.6.3).</li> <li>If you remove the dust drawer, please always empty it.</li> </ul>
<b>The dust drawer full signal lights up, although the drawer is not yet full.</b>	<ul style="list-style-type: none"> <li>The set value for the dust drawer fill level is too low.</li> </ul>	<ul style="list-style-type: none"> <li>Set a higher value for the fill level (section 5.6.3).</li> </ul>
<b>Suction hose cannot be closed by pressing the suction hose button.</b>	<ul style="list-style-type: none"> <li>An electrical appliance which is connected to the suction hose is still in operation.</li> </ul>	<ul style="list-style-type: none"> <li>Switch the connected electrical appliance off.</li> </ul>
<b>The letter from a suction channel is blinking.</b>	<ul style="list-style-type: none"> <li>The suction channel is closed and the switch off threshold has not yet been reached.</li> </ul>	<ul style="list-style-type: none"> <li>No fault.</li> </ul>
<b>The letter from a suction channel is blinking and the suction channel cannot be opened.</b>	<ul style="list-style-type: none"> <li>In the single use mode only one suction channel can be opened at a time. There is currently a suction channel already open.</li> </ul>	<ul style="list-style-type: none"> <li>Close one open suction channel first and then open the other.</li> <li>Select the multi use mode, but to do so all suction channels must first be closed.</li> </ul>
<b>It is not possible to change between single and multi use mode.</b>	<ul style="list-style-type: none"> <li>A suction channel is still open.</li> </ul>	<ul style="list-style-type: none"> <li>Close all suction channels.</li> </ul>
<b>There is no audible signal.</b>	<ul style="list-style-type: none"> <li>The audible signal is deactivated.</li> </ul>	<ul style="list-style-type: none"> <li>Activate audible signal (Section. 5.6.2).</li> </ul>
<b>Suction hose cannot be closed by pressing the suction hose button.</b>	<ul style="list-style-type: none"> <li>An electrical appliance connected to the QCB is still in use.</li> </ul>	<ul style="list-style-type: none"> <li>Switch the appliance connected to the QCB off, see section 5.4.4.</li> </ul>
<b>The extraction unit requires filter cleaning at frequent intervals.</b>	<ul style="list-style-type: none"> <li>The diameter of the suction tube is too narrow, because a small suction port was connected to the the suction point.</li> <li>The fine filter unit is so dirty, that it is no longer possible to clean.</li> <li>The system has a leak.</li> </ul>	<ul style="list-style-type: none"> <li>Use the original suction tubes and install a suction port adapter if necessary (see accessories).</li> <li>Change fine filter adapter (section. 6.2).</li> <li>Check to make sure that the dust drawer and the fine filter unit have been replaced correctly and that the seals are not damaged.</li> <li>Check to make sure that the dust area cover panel has been replaced correctly and is fixed securely in place.</li> <li>Check to make sure the pinch valve seals are secure.</li> </ul>

Error	Possible Cause	Corrective Action
<b>The suction does not start when an electrical appliance connected to the QCB is activated.</b>	<ul style="list-style-type: none"> <li>• The extraction unit is not switched on.</li> <li>• Device fuse has been triggered.</li> <li>• QCB is not connected with the extraction unit.</li> <li>• On / off threshold of the suction channel is too high.</li> <li>• A connected remote control has priority over the QCB for this suction channel.</li> </ul>	<ul style="list-style-type: none"> <li>• Switch the extraction unit on (section 5.1).</li> <li>• Check fuses, if necessary replace fuse (section 6.3.1).</li> <li>• Connect the QCB to the extraction unit (section 4.6.1).</li> <li>• Adjust the on / off threshold (section 5.6.4).</li> <li>• Press the mouse scroll wheel from the relevant suction channel to give the QCB priority again.</li> </ul>
<b>The electrical appliance connected to the QCB cannot be operated.</b>	<ul style="list-style-type: none"> <li>• QCB is not connected to a power supply.</li> </ul>	<ul style="list-style-type: none"> <li>• Connect the QCB to a power supply (section 4.6.1).</li> </ul>
<b>Suction begins, even though no appliance connected to the QCB has been activated.</b>	<ul style="list-style-type: none"> <li>• The on / off threshold from the suction channel is too low.</li> </ul>	<ul style="list-style-type: none"> <li>• Check which suction channel is open, the relevant letter will be shown in the display.</li> <li>• Adjust the on / off threshold from that particular suction channel (Section 5.6.4).</li> </ul>
<b>The suction does not stop when the electrical appliance is switched off.</b>	<ul style="list-style-type: none"> <li>• The electrical appliance was not connected to the QCB.</li> <li>• The on / off threshold is set too low.</li> <li>• A connected remote control has priority over the QCB for this suction channel.</li> </ul>	<ul style="list-style-type: none"> <li>• Connect the electrical appliance to the QCB (section 4.6.2).</li> <li>• Adjust the on / off threshold (section 5.6.4).</li> <li>• Press the mouse scroll wheel from the relevant suction channel to give the QCB priority again.</li> </ul>
<b>The fault “Err14” appears when programming the on / off threshold for a handpiece.</b>	<ul style="list-style-type: none"> <li>• The rpm you have selected on the handpiece means the difference between stand-by power and operating current is too low.</li> </ul>	<ul style="list-style-type: none"> <li>• Program a higher rpm on the handpiece to activate the on / off threshold.</li> </ul>
<b>It is not possible to find a satisfactory setting for the on / off threshold.</b>	<ul style="list-style-type: none"> <li>• Some older handpieces do not provide sufficient signal for the automatic switch on.</li> </ul>	<ul style="list-style-type: none"> <li>• Switch the suction channel on manually (see section 5.4.3).</li> </ul>
<b>The fuse in the QCB is triggered when a connected electrical appliance is switched on.</b>	<ul style="list-style-type: none"> <li>• The total power consumption of the activated appliances is too high.</li> </ul>	<ul style="list-style-type: none"> <li>• Please observe the max. connected load for the QCB (see section 8.2).</li> </ul>
<b>An audible sound of compressed air is escaping from the deactivated device.</b>	<ul style="list-style-type: none"> <li>• The float for the automatic drainage function in the pneumatic unit has slipped out of its correct position, i.e. due to horizontal transportation.</li> </ul>	<ul style="list-style-type: none"> <li>• Remove the device from the compressed air connection and then reconnect it. If necessary, repeat this step 2 - 3 times.</li> </ul>

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