

*Renfert launches high-performance trimmer with clever features*

### **MT premium: A genius model trimmer**

The new wet trimmer MT premium from Renfert demonstrates masterful expertise in reliable professional equipment. On top of its high performance, the latest member of the MT trimmer family impresses with refined details. Once more, Renfert underlines its commitment to high-quality laboratory equipment that makes work convenient and simple – „making work easy“.

Uncompromisingly powerful and cleverly detailed - with the MT premium wet trimmer, Renfert introduces a paradigm shift in model trimming. The new high-performance trimmer is based on the tried-and-tested technical characteristics of the MT 3 trimmer (e.g. high performance, tool-free cleaning) with the addition of clever features, e.g. the aqua stop function and a spotlight.

And with its superb performance, the MT premium is a real powerhouse too. The construction and design leave nothing to be desired. Slotted angle lines on the trimmer table enable easy guidance of the plaster model. Thanks to a 10° inclination of the unit, an optimal view of the model is ensured. The trimmer table angle can be quickly adjusted using the handy tilt mechanism (90°/98°). Optimally positioned, the spray tube prevents the trimmer disc from clogging and the anti-drop function ensures a clean, dry work surface. In addition, the MT premium is equipped with the aqua stop function that helps saving water. The integrated work light, a spotlight that provides optimal lighting of the work surface on the trimmer table, is a real touch of genius.

### **ORTHO guide for orthodontic models**

With the ORTHO guide retrofit set, the MT premium is easily transformed into a special trimmer. A template and a special trimmer table ensure fast and precise trimming of orthodontic models. An angulation support ensures defined angles. Optimal alignment of the median palatine raphe is eased by a laser guide light.

The MT premium truly triumphs either in the laboratory or in the dental office. Therefore, the high-performance trimmer stands for next-level model trimming. Sheer strength and reliable performance meet sophisticated ease of use and solid expertise. And the result is unmistakably Renfert: „making work easy“.

**1,927 characters**

---

### **More information**

Renfert GmbH | Untere Giesswiesen 2 | 78247 Hilzingen | [www.renfert.com](http://www.renfert.com)

### **Press contact**

Alberta Domnick | Tel. +49 7731 8208 898 | E-Mail: [Alberta.Domnick@renfert.com](mailto:Alberta.Domnick@renfert.com)