



making work easy



POWER steamer

NEW Steam cleaners

NEW POWER steamer 2

Steam cleaning unit with automatic filling

The powerful POWER steamer 2 of the technically optimized new generation stands for maximum reliability and durability, as well as outstanding cleaning performance in the dental laboratory. The POWER steamer 2 is consistently designed for continuous, uninterrupted operation. Any kind of typical soiling can be easily and effectively removed. Using the mains water connection and integrated pump, the pressure tank is filled automatically and ensures availability of steam at all times.

Advantages

- Exceptionally durable heating system: the heating element is permanently protected against water, limescale, and cleaning chemicals thanks to being fully embedded in the solid aluminum boiler base.
- Availability of steam at all times – even during heavy-duty continuous use – thanks to automatic filling using the mains water connection and integrated pump.
- High operational reliability thanks to the specially designed solenoid valve – easily replaceable by the user in just a few simple steps if required.

Details

- Continuous and effective cleaning performance thanks up to 2,000 W and 4.5 bar working pressure in combination with innovative steam nozzle technology.
- High functional reliability and dependability thanks to a real-time calcification indicator.
- Optimal control through display of the current working pressure using a manometer.
- Simple cleaning of the pressure tank thanks to an extra-large service opening.



tool included



easy use



compact



Protected heating element

Technical data

Permissible mains voltage	220-240 V 120 V 100 V
Permissible mains frequency	50/60 Hz
Power consumption	2000 W 1280 W 1550 W
Pressure tank volume	4 l // 1.06 gal
Volume (Recommended fill quantity)	2,8 l // 0.74 gal
Volume (Maximum fill quantity)	3 l // 0.79 gal
Working pressure	4,5 bar // 65.3 psi
Connection pressure	1 – 5 bar // 14.5 – 72.5 psi
Dimensions (W x H x D)	370 x 422 x 350 mm // 14.6 x 16.6 x 13.8"
Weight	~10.1 kg // ~22.3 lbs

Ordering information

POWER steamer 2, 230 V	No. 18460000
POWER steamer 2, 120 V	No. 18461000
POWER steamer 2, 100 V	No. 18462000
Wall bracket	No. 18450200



NEW POWER steamer 1

Steam cleaning unit for manual filling

The powerful POWER steamer 1 of the technically optimized generation sets new benchmarks in reliability and durability and impresses with its outstanding cleaning performance in the dental laboratory. Any kind of typical soiling can be easily and effectively removed. Manual filling allows for flexible placement.

Advantages

- Exceptionally durable heating system: the heating element is permanently protected against water, limescale, and cleaning chemicals thanks to being fully embedded in the solid aluminum boiler base.
- High functional reliability and dependability thanks to a real-time calcification indicator.
- High operational reliability thanks to the specially designed solenoid valve – easily replaceable by the user in just a few simple steps if required.

Details

- Continuous and effective cleaning thanks to up to 2,000 W and 4.5 bar working pressure in combination with innovative steam nozzle technology.
- Flexible location selection thanks to manual filling.
- Simple cleaning of the pressure tank thanks to an extra-large service opening.
- Safe and convenient steam cleaning thanks to an ergonomic handpiece (the housing safely dissipates electrostatic charge).
- Easy filling of water and descaling solution thanks to handy, hopper-shaped fillers.
- Wall mounting possible using an optional wall bracket.



tool included



easy use



compact



Extra-large service opening

Technical data

Permissible mains voltage	220-240 V 120 V 100 V
Permissible mains frequency	50/60 Hz
Power consumption	2000 W 1550 W 1280 W
Pressure tank volume	4 l // 1.06 gal
Volume (Recommended fill quantity)	2,8 l // 0.74 gal
Volume (Maximum fill quantity)	3 l // 0.79 gal
Working pressure	4,5 bar // 65.3 psi
Dimensions (W x H x D)	370 x 422 x 350 mm // 14.6 x 16.6 x 13.8"
Weight	~9.5 kg // ~20.9 lbs

Ordering information

POWER steamer 1, 230 V	No. 18450000
POWER steamer 1, 120 V	No. 18451000
POWER steamer 1, 100 V	No. 18452000
Wall bracket	No. 18450200



STEAM tube

Intelligent steam cleaning station

The functional solution that efficiently masters the challenges of steam cleaning. STEAM tube is designed to specifically collect dirt particles such as wax and residual polish and protects against the loss of small, delicate parts. Thanks to controlled steam discharge, it also improves visibility during steam cleaning.

Advantages

- Maximum cleanliness and rapid cleaning thanks to a unique dirt collection concept and effortless removal of the hopper, filter, and internal collection container.
- Protects valuable small components with the safety collection filter and can catch small parts safely at any time.
- Ensures clear visibility of the work object thanks to targeted steam discharge and enables an unobstructed view of the object being cleaned.



easy use



perfect view



compact

Technical data

Weight	3 kg
Dimensions (in the position foot bar 90° to table) (W x H x D)	270 x 300 – 390 x 400 mm
Tube diameter	200 mm
Drain hose length	1 m
Ø Drain hose exterior	16 mm

Ordering information

STEAM tube (already mounted), ergonomic cleaning spatula, drain hose, 2x sealing plugs, Quick Start Guide	No. 18420000
---	--------------

POWER steamer descaler

Concentrated descaling solution
(for POWER steamer 1 and 2)



POWER steamer descaler is a concentrated descaling solution for effective as well as gentle descaling of both POWER steamers.

Ordering information

POWER steamer descaler	No. 18450100
------------------------	--------------

POWER steamer water softener

Water softening system for POWER steamer 2



silent



compact



Reliable water softening and filter system, especially adapted for POWER steamer 2.

The system provides optimal protection of the POWER steamer, increasing its efficiency, reliability, and durability.

Ordering information

POWER steamer water softener	No. 18460100
------------------------------	--------------

