

Wax-up techniques

Analog mastery for digital precision

Accurate waxing-up

Waxlectric II is the professional device in the electric wax knife range which can be operated with two handpieces simultaneously.



Distortion-free wax copings

Reliable wax dipping pot hotty LED with LED display and membrane keypad for the even and accurate heating temperature of dipping waxes.

The precise die preparation with GEO Dip creates the ideal basis for accurately fitting conventional and all-ceramic restorations.



Tension-free modeling

Thanks to a special wax formulation, the GEO wax wires are tension-free and develop only low restoring forces.





Crown and
bridge techniques



Esthetic
wax-up



All-ceramics and
press ceramics



Special wax

Work flexibly

The GEO Classic range offers high-quality modeling waxes in various degrees of hardness – ideal for every dental requirement, from single crowns to complex bridges. The waxes impress with their excellent modeling, carving and milling properties and burn out without leaving any residue – perfect for full and press ceramics. GEO Classic waxes can be processed precisely and in a controlled manner, whether using an open flame or an electric wax knife.



	No.				
GEO Classic Basic Set	493 0004	■	■	■	■
GEO Classic	497 0100	■	–	–	–
GEO Classic	497 0200	■	–	–	–
GEO Classic	497 0300	■	–	–	–
GEO Classic	497 1200	■	■	–	–
GEO Classic	497 1300	■	■	–	–
GEO Classic Avantgarde	495 0100	■	–	–	–
GEO Classic Avantgarde	495 0200	■	–	–	–

	No.				
GEO Classic Avantgarde	495 0300	■	–	–	–
GEO Classic Natural	499 0300	■	–	■	–
GEO Classic Natural	499 0400	■	■	■	–
GEO Classic Snow-white	499 0101	■	–	■	–
GEO Classic Snow-white	499 0201	■	■	■	–
GEO Classic cervical wax	486 1000	■	■	–	■
GEO Classic cervical and undercut wax	489 1000	■	■	–	■
GEO Classic milling wax	485 1000	■	–	–	■

Work efficiently

The GEO Crowax product range stands for cost-effective modeling without compromise. The medium-hard waxes are ideal for crown, bridge and model casting techniques – easy to use, reliable results.



	No.				
GEO Crowax	475 0100	■	■	–	–
GEO Crowax	475 0200	■	■	–	–
GEO Crowax	475 0300	■	■	–	–

	No.				
GEO Crowax	475 0400	■	–	–	–
GEO Crowax	475 0500	■	–	–	–
GEO Crowax cervical and undercut wax	475 0600	■	■	–	■

Opaque

Transparent



Find out more

www.renfert.com/waxup