



Renfert

making work easy

SILENT XS

NEW The mobile suction unit for precise touch-ups

NEW SILENT XS

The mobile suction unit for precise touch-ups

SILENT XS is a fully mobile suction unit, specifically designed for dental offices and in-house laboratories. Its compact design is impressive, and makes it ideal for minor, small-scale finishing work - for example when working on temporary restorations, splints or dentures. With its high suction power, it collects debris right at the source and ensures a clean working environment. Thanks to battery operation, the device can be moved flexibly between different locations. Its low-noise operation makes it the ideal companion for efficient work - even during mobile applications.

Advantages

- Flexible – cordless and battery-powered for maximum freedom of movement – from chairside to in-house laboratory.
- Simple – one-handed operation and LED lighting for an optimum view of the workpiece. Perfect for quick grinding tasks.
- Clean – clean working air and effective dust protection thanks to H13/HEPA filter with 99.99% collector efficiency. Bag-free.

Details

- Low-noise operation for a calm treatment environment.
- Compact design for space-saving integration and easy storage.
- Long operating times thanks to battery and mains operation.
- Easy to clean with smooth, disinfectable surfaces in a validated process.
- Easy to use with ergonomic design and one-hand operation.

Technical data

Permissible mains voltage (power supply unit)	100–240 V
Permissible mains frequency (power supply unit)	50–60 Hz
Volume flow (max.)	510 l/min
Filter technology	Fine filter
Filter quality	HEPA / H13
Sound pressure level (LpA at max. volume flow)	52,7 dB (A)
Weight (without battery)	1,12 kg // 2.47 lbs
Dimensions (W x H x D) (incl. shield)	223 x 248 x 193 mm // 8.8 x 9.8 x 7.6 "
Battery (not included)	8 x AA (NiMH rechargeable) Attention: Do not use Li-Ion!



easy use




silent



perfect view



compact



Effective dust collection with practical protective grille – for clean and safe work

Ordering information

SILENT XS	No. 29220000
SILENT XS, US / JP	No. 29221000
SILENT XS, ROW	No. 29223000
Working block	No. 29220001
Filter set (3 pieces incl. exhaust air filter)	No. 29220002

