

making work easy

≅**\84.31t**800

Rentert

Rever ENUMER



NEW lay:art evo

Bionic hair brush

lay:art evo is a state-of-the-art premium brush with bionic hair. Thanks to the chemical composition of the newly developed fibers and their specific processing, the brush bristles can absorb moisture and replicate the flame shape that is typical of high-end brushes. Together, both of these aspects demonstrate the characteristic moisture management, which also distinguishes natural bristle brushes. As a result, lay:art evo impresses with incredible dental product detail and pragmatic functionality.

Advantages

- Remarkable durability due to innovative bionic hair.
- Optimum moisture management thanks to the special molecular structure of the fibers.
- Ideal handling of ceramic masses during layering due to the fineness, elasticity, and stability of the brush tip.
- Refills in different sizes enable simple, cost-effective replacement of the brush tip. Also compatible with brush systems supplied by other manufacturers.


Functional details

- Delicate tip with durable tension and elasticity impresses during targeted application of ceramic.
- The tension of the fine bristle tip adapts to usage and is maintained throughout layering.
- Longer service-life without separation or "bunching" due to bionic hair.
- Brush handle in different sizes made from top quality anodized aluminum with smart ergonomics.
- State-of-the-art dental lifestyle thanks to brush tip made of bionic hair and functional design.

Ordering information

-		
lay:art evo fine handle	1 piece, without interchangeable tip	No. 17240000
lay:art evo standard handle	1 piece, without interchangeable tip	No. 17240001
[1] lay:art evo size 2	1 piece, with handle fine	No. 17240002
[2] lay:art evo size 4 flame	1 piece, with handle standard	No. 17240004
[3] ay:art evo size 4 cone	1 piece, with handle standard	No. 17240014
[4] lay:art evo size 6 flame	1 piece, with handle standard	No. 17240006
[5] lay:art evo size 6 cone	1 piece, with handle standard	No. 17240016
[6] lay:art evo size 8 flame	1 piece, with handle standard	No. 17240008
[7] lay:art evo size 8 cone	1 piece, with handle standard	No. 17240018
[8] lay:art evo color	1 piece, with handle fine	No. 17240100
[9] lay:art evo opaque	1 piece, with handle fine	No. 17240200
lay:art evo size 2 refill	2 pieces, without handle	No. 17240012
lay:art evo size 4 flame refill	2 pieces, without handle	No. 17240114
lay:art evo size 4 cone refill	2 pieces, without handle	No. 17241114
lay:art evo size 6 flame refill	2 pieces, without handle	No. 17240116





Choice of two types of tips depending on the size: Flame and Cone

lay:art evo size 6 cone refill	2 pieces, without handle	No. 17241116
lay:art evo size 8 flame refill	2 pieces, without handle	No. 17240118
lay:art evo size 8 cone refill	2 pieces, without handle	No. 17241118
lay:art evo color refill	2 pieces, without handle	No. 17240110
lay:art evo opaque refill	2 pieces, without handle	No. 17240210
lay:art evo build up set	Handle standard with 3 brush sizes (4, 6, 8 flame)	No. 17240300
lay:art evo detail up set	Handle fine with 3 brush sizes (color, opaque, 2)	No. 17240400
lay:art evo refresher	1 Brush stand with glass (without brushes)	No. 17240600

2