

making work easy

Easyclean TEC

NEW Ultrasonic cleaning unit



The Easyclean TEC ultrasonic cleaning unit supports powerful yet simple cleaning of dental objects such as impression trays and instruments. Cleaning is enhanced by using appropriate chemical cleaning agents, increasing the bath temperature, and regulating the application time. What's more, the permanent sweep function enables effective distribution of the ultrasonic frequency, further increasing the impact of cleaning.

Advantages

- An ultrasonic frequency of 37 kHz ensures an optimal cleaning result.
- For effective distribution of the ultrasonic frequency, Easyclean TEC uses the permanent sweep function.
- The integrated heating element means optimized cleaning even with the most stubborn contaminations.

Details

- For efficient and energy-conscious cleaning, both the "eco" and the "dynamic" ultrasonic mode allow adaption to the level of contamination.
- The stainless-steel tank is highly cavitation resistant ensuring an extraordinarily long service life.

• For a disruption-free working environment, the acoustic feedback function indicates when cleaning has been completed.

making work easy

Thanks to a modern, streamlined design, this premium quality ultrasonic cleaning unit blends seamlessly into the dental environment. As well as offering low-noise operation, the unit can also be easily maintained as part of the daily routine. The control elements and the stainless-steel surface of the housing are easy to clean.



Application example: Cleaning of instruments

Technical data

Permissible mains voltage	220-240 V 115-120 V 100 V
Permissible mains frequency	50/60 Hz
Power consumption	280 W
Ultrasonic frequency	37 kHz
Ultrasonic effective power	80 W
Heating power	200 W
Dimensions (W x H x D) (unit, exterior dimensions)	305 x 230 x 190 mm // 12 x 9.1 x 7.5"
Weight (incl. lid)	~3.3 kg // ~7.3 lbs
Volume (tank, max. volume)	2.7 I // 0.71 gal
Volume (tank operating volume)	1.7 l // 0.45 gal

Ordering information

Easyclean TEC, 220 – 240 V	No. 18520000
Easyclean TEC, 115 – 120 V	No. 18521000
Easyclean TEC, 100 V	No. 18522000

Accessories

Easyclean TEC

→ ha 4 1	GO 2011 speed	Suitable for cleaning different instruments and aids quickly and thoroughly in only 50% of the time compared to conventional plaster and alginate solvents.	2 I (0.53 gal)	No. 20120000
	Plastic lid, gray	Quicker heating. Protects against evaporation and dust infiltration.	1 piece	No. 18500001
	Stainless steel cover	Cover for use with the cleaning jar or plastic cup with lid.	1 piece	No. 18500002
	Stainless steel basket	For storage of cleaning items. To protect the floor of the ultrasonic tank.	1 piece	No.18500003
	Stainless steel immersion basket	For cleaning small and very delicate items. For use with the stainless steel cover or with the cleaning jar.	1 piece	No. 18500004
	Plastic acid bath insert	For acids and fluids not suitable for use in the stainless steel tank.	1 piece	No. 18500005
	Cleaning jar	For the use of additional cleaning liquids. For use in the stainless steel cover. 600 ml	1 piece with lid and rubber ring	No. 18500006
	Plastic cup with lid	Ideal for cleaning smaller items and working with acids. Can be also used in conjunction with the stainless steel cover.	1 piece	No. 18500007

