

making work easy



Renfert-Scanspray extra fine

NEW Matting spray

NEW Renfert-Scanspray extra fine

Matting spray

Thanks to its micro-fine, homogeneous structure, Renfert-Scanspray extra fine can be applied very thinly. Enhanced adhesion, especially on polished ceramics/metals and reflecting surfaces.

Advantages


- Thinner layer thickness (< 1 µm) and higher opacity compared to conventional scan sprays. This ensures more homogeneous and more detailed scanning results.
- High separation sharpness with detailed edge representation, especially on delicate structures, for precise scanning results.
- No flake and bubble formation (agglomeration), thus reduction of dimensional errors.
- Particularly good output due to very thin and micro-fine layer application.







Micro-fine application for detailed scanning results



Left: uniform application of Renfert-Scanspray extra fine Right: uneven coating with competitive products

Ordering information

Renfert-Scanspray extra fine, 200 ml

No. 17312000

