

**Renfert**

making work easy



# SYMPRO

Cleaning system for removable  
dental restorations

# SYMPRO

## Cleaning system for removable dental restorations

Compact cleaning unit especially suitable for cleaning removable dental restorations, splints and orthodontic appliances.

### Advantages

- Customer retention – a prophylaxis service for denture wearers.
- Patient satisfaction – cleaning of splints and orthodontic appliances.
- Process reliability – validated, hygienic reprocessing of system components.

### Details

- Positive effect on the patient's overall health due to regular denture cleaning.
- High time savings compared to manual methods thanks to effective and automated cleaning process.
- Externally tested material compatibility with standard dental materials.
- Optimized cleaning pin geometry, gentle and effective at the same time.
- Excellent cleaning performance thanks to optimal bowl inclination.
- Easy cleaning of temporary restorations in combination with SYMPRO mini cup and temp:ex temporary cement remover.

### Making work easy

Uncompromising cleaning with SYMPRO system. The effective help:ex cleaning agents are available for cleaning with the SYMPRO denture cleaning unit.



easy use



silent



compact



Before and after cleaning with SYMPRO and help:ex

### Technical data

Permissible mains voltage	100–240 V
Permissible mains frequency	50/60 Hz
Power consumption	90 VA
Cleaning intensity	low / medium / high
Dimensions (W x H x D)	150 x 240 x 280 mm // 5.91 x 9.45 x 11.02"
Weight (without bowl)	~3 kg // ~6.6 lbs

### Ordering information

SYMPRO, 100–240 V	No. 67001000
-------------------	--------------

### Ordering information accessories

help:ex plaque f	4 x 1 l (4 x 0.26 gal)	No. 67000100
help:ex plaque p	20 x 20 g (0.7 oz.) cleaning powder + 20 x 4,5 g (0.16 oz.) neutralizer	No. 67000000
help:ex discolor f	1 l (0.26 gal)	No. 67000200
Further accessories	<a href="http://www.renfert.com/sympro">www.renfert.com/sympro</a>	

