

Press release

Hilzingen, Germany - October 2023

Renfert launches a premium quality ultrasonic cleaning unit Easyclean TEC: Top Cleaning Performance – Top Design

Renfert, based in Hilzingen, Germany, presents Easyclean TEC, an ultrasonic cleaning unit that is designed with the needs of dentistry in mind. With an ultrasonic frequency of 37 kHz, even stubborn stains can be thoroughly removed. It is also incredibly easy to use and keeps noise to a minimum.

The new Easyclean TEC ultrasonic cleaning unit from Renfert provides powerful cleaning in dental laboratories and the dental office, effortlessly cleaning impression trays, articulation plates, and other equipment. Cleaning is based on the physical action of ultrasonic technology and enhanced by using chemical cleaning agents, increasing the bath temperature, and regulating the application time. What's more, the sweep function enables effective distribution of the ultrasonic frequency, further increasing the impact of cleaning.

- Efficient cleaning thanks to an integrated heating element and continuous sweep function
- Frequency of 37 kHz
- Two pre-configured cleaning programs
- Intuitive operation
- Exceptionally long service life

A choice of two ultrasonic modes is available depending on the level of contamination. "eco" stands for energy-efficient and quiet operation, while "dynamic", with its optimized ultrasonic field distribution, enables comprehensive, thorough cleaning of stubborn stains. The user-friendly control panel provides relevant information such as the operating mode, duration, or temperature. Thanks to its modern design, this premium quality ultrasonic cleaning unit blends seamlessly into the surroundings. The control elements and the stainless-steel surface of the housing are easy to clean too.

A premium quality product "Made in Germany", Easyclean TEC offers an exceptionally long service life. Renfert's commitment to durability is demonstrated by a three-year guarantee. Easyclean TEC impresses with a robust design that is typically Renfert, blending product hygiene with quality you can rely on, ease of use, and a timeless design. The result is unmistakably Renfert: making work easy.

2,137 characters

More information