



making work easy

The background of the advertisement is a close-up photograph showing a dental procedure. A hand holds a small, white, irregularly shaped dental restoration. A dental handpiece with a dark, cylindrical polishing wheel is being used to polish the restoration. The wheel is in contact with the restoration, and a bright, circular reflection is visible on its surface. The background is a soft, out-of-focus light blue.

Renfert Polish ZrO_2 starter kit

NEW Polishing kit for zirconium dioxide restorations

NEW Renfert Polish
ZrO₂ starter kit

Polishing kit for zirconium dioxide res-
torations

Renfert Polish ZrO₂ starter kit combines the es-
sential components for processing zirconia resto-
rations into a coordinated set. The polishing paste
and polishers are tailored to the material and typ-
ical indications, helping users carry out the polish-
ing procedure in a structured, controlled manner to
achieve reproducible results.

Advantages

- Efficient workflow: fast and uniform polishing of zirconium dioxide restorations for a reliably smooth surface.
- Material-friendly polishing: the three-step polishing system ensures consistent surface quality and delivers optimal high gloss.
- Long-lasting high gloss: polishers with perfectly matched geometries enable versatile use on monolithic and partially veneered zirconium dioxide restorations.

Details

- First prepolishing: Prepolisher ZrO₂
- Second prepolishing: Prepolisher twist ZrO₂
- High gloss polishing: Bison polishing brush + Renfert Polish ZrO₂ polishing paste



.....

More polishing kits

Renfert Polish 3D-printed ma-
terials starter kit: a specially
coordinated starter kit for the
diamond polishing paste
Renfert 3D-printed materials.



Renfert Polish all-in-one
starter kit: a specially coor-
dinated starter set for the
Renfert polishing paste
all-in-one.



.....

Ordering information

Renfert Polish ZrO ₂ starter kit		No. 5101200
Renfert Polish ZrO ₂	Polishing paste, 18 g	No. 5101100
Prepolisher ZrO ₂ (Lens)	26 mm, 1 piece	No. 95 0000
Prepolisher twist ZrO ₂ (Lamellenrad)	26 mm, 1 piece	No. 95 1000
Bison	14 mm, 18 mm, 25/100 pieces	No. 76xx00

Scope of delivery

Renfert Polish ZrO₂ polishing paste 18 g (No. 5101100), incl. 1x Prepolisher ZrO₂ (Lens) 26 mm (No. 95 0000) for the initial prepolishing of the entire surface, 1x Prepolisher twist ZrO₂ (Lamellenrad) 26 mm (No. 95 1000) for the second prepolishing of the entire surface and 1x Bison polishing brush 18 mm (No. 766 1000) for high-gloss polishing.