

# Products for All-Ceramic Workflows

The Renfert logo consists of the word "Renfert" in a bold, white, sans-serif font, centered within a red rectangular background. The text is flanked by two white, curved lines that resemble a stylized "R" or a pair of wings.

**Renfert**

making work easy

A close-up, low-angle shot of a blue dental mold. The mold is circular and features several circular cavities. Each cavity contains a small, blue, conical object, likely a wax or ceramic mold for dental work. The lighting is dramatic, highlighting the texture and color of the mold.

Please check  
the QR codes for  
our new videos



[www.renfert.com](http://www.renfert.com)

# EASY blank wax

the blank for  
state-of-the-art hybrid  
workflows

Many have the impression, that dental restorations are now all produced using full digital workflows. However, a day in a dental laboratory paints a very different picture. Yes, the CAD/CAM technology is well established, however, conventional casting and pressing techniques are still being used regularly. For this reason, Renfert has extended its range of products for digital manufacturing with the EASY blank wax for the hybrid processing techniques CAD/Cast and CAD/Press.




The new video for  
EASY blank wax



item number	description	unit
7700014	EASY blank wax 14 mm	piece
7700020	EASY blank wax 20 mm	piece

## Dynex Brillant Grinding discs

Diamond grinding discs with glass fiber reinforcement, created for processing and trimming ceramic materials, especially suited for lithium disilicate.



The short video for cutting and grinding with Dynex



## Dynex Brillant Separating discs

Flexible diamond separating discs with glass fiber reinforcement, created for cutting, separating and trimming veneering porcelain, zirconium oxide and lithium disilicate.




item number	description	unit
560220	Dynex Brillant Separating disc 0,2x20mm, (10 pcs.)	package
560320	Dynex Brillant Separating disc 0,3x20mm, (10 pcs.)	package
560820	Dynex Brillant Separating disc 0.8x20mm, (1 pcs.)	piece
562520	Dynex Brillant Separating disc 0,25x20 mm, (10 pcs.)	package

## Renfert Polish all-in-one

Polishing paste for universal use, with a wide spectrum of diamond grain sizes. Quick, high gloss finish on various different materials such as lithium-disilicate, zirconium oxide, veneering and hybrid ceramics such as high performance polymers and veneering composites.

The short video for Renfert Polish



## Renfert Polish ZrO<sub>2</sub>

Polishing paste for zirconium oxide with specially coordinated abrasive properties due to extra coarse diamond grit.

## Renfert Polish LiSi<sub>2</sub>

Polishing paste specially for lithium disilicate with an exactly coordinated amount of coarse diamond grit.



item number	description	unit
5100100	Renfert Polish all-in-one 18 g	piece
5101100	Renfert Polish ZrO <sub>2</sub> 18 g	piece
5102100	Renfert Polish LiSi <sub>2</sub> 18 g	piece



The short video for Renfert Polish



## Bison Polishing brush

This unique mixture of bristles ensures a powerful and quick high-lustre polish on alloys, porcelain and acrylic facings and hybrid ceramics such as high performance polymers and veneering composites.

item number	description	unit
7631000	BISON polishing brushes, 14mm (25 pcs.)	package
7631100	BISON polishing brushes, 14mm (100 pcs.)	package
7661000	BISON polishing brushes, 18mm (25 pcs.)	package
7661100	BISON polishing brushes, 18mm (100 pcs.)	package



## Silicone Polishers

For pre-polishing porcelain and metal

item number	description	unit
860000	Sil. polishers wheel 22mm, light grey (100 pcs.)	package

## Goat hair brush Polishing brush

For polishing precious metal and acrylics surfaces.

item number	description	unit
2040000	Goathair brushes (12 pcs)	package

## Renfert-Scanspray

Maximum precision thanks to a very thin layer application and homogenous grain size.

item number	description	unit
17310000	Renfert-Scanspray 200ml	piece



## Jacketgrip

### Holding instrument with diamond tips

Tweezers for a secure hold on crowns and bridges.



## Keramogrip

### Holding instrument with diamond tips

Holders for crowns and bridges during porcelain or acrylic.



item number	description	unit
11050300	JACKET-GRIP tweezers	piece

item number	description	unit
11090300	KERAMOG RIP tweezers, angled	piece

## Basic quattro IS

### Fine sandblasting unit



Sandblasting Occlusal Fissures  
with Basic Quattro IS

Professional fine sandblasting unit with 2-4 blasting tanks and the "IS" immediate stop technology, for controlled sandblasting without an after-flow of sandblasting material



## Twister evolution

### Vacuum mixing unit

High speed, vacuum mixing unit, diaphragm pump for the vacuum and eight programmable mixing functions.



item number	description	unit
29590000	Basic quattro IS 25-70 $\mu$ / 70-250 $\mu$	piece
18280000	Twister evolution	piece



## SILENT compactCAM

Extraction unit for CAM systems

Compact, bag-free extraction unit for use with compact, dental CAM units.



## SILENT powerCAM EC

Extraction unit for CAM systems

Powerful, bag-free extraction unit with automatic filter. Cleaning and EC motor technology for heavily used CAM units.

item number	description	unit
29342500	SILENT compactCAM	piece
29390000	SILENT powerCAM EC	piece



## Interface cable

item number	description	unit
29340005	Interface cable Type A for vhf, 5-pin	5m
29340006	Interface cable Type B for Roland DG	5m
29340008	Interface cable Type C for imes-icore	5m
29340009	Interface cable Type D for Amann Girrbach	5m
29340010	Interface cable Type E for Yenadent/Origin + Nema adapter	5m
29340011	Interface cable Type F for Zirkonzahn	5m
29340012	Interface cable Type F for Zirkonzahn + C14 adapter	5m
29340013	Interface cable Type G (6-pin)	5m

# AUTO spin



## Pin drilling unit

Using the AUTO spin pin drill, precise working models with plastic base plates can be fabricated. The drilling unit is also compatible with other model-making systems like Giroform® or Zeiser®.

When uniquely combined with perfectly matched materials of the AUTO spin model-making system, even the most sophisticated model fabrication is simple.



item number	description	unit
18600000	AUTO spin	piece



The new AUTO spin video



95  
years

*Renfert*

95

**Renfert**

**WORKFLOW  
GUARANTEE**

3 year guarantee  
10 year spare parts service  
Activity guarantee

**Renfert**

making work easy



[www.renfert.com](http://www.renfert.com)