

**Renfert**

making work easy



EASYview+

Dental Visual Communicator

# EASY view+

## Dental Visual Communicator

Renfert EASY view+, the Dental Visual Communicator for digital object viewing. Digital visual capture allows you to better visualize, edit and share your daily work. In addition to classic live viewing with three zoom levels, you can spontaneously generate images at the touch of a button. Create videos via the direct interface to your laptop and share e.g. patient cases simultaneously with others. Simply use your familiar programs to edit the image files.

### Advantages

- Three magnification levels (4x, 15x and 20x), for comprehensive or detailed/microscopic viewing.
- Connection to your laptop\* via an online meeting tool makes case discussion with the dentist easy and efficient.
- Free choice of monitor\*\* allows for an individual adaptation to your needs.

Optional accessories: Renfert CONNECT \*\*\* (Renfert CONNECT stick and Renfert CONNECT app for integration into the network and wireless data synchronization between EASYview+ data storage and app).

\* System requirements: Microsoft Windows® 10 and higher; Intel® Core™ i3 and higher | Apple® macOS® 10.15 and higher; Intel® Core™ i5 and higher.

\*\* 3D function is only possible via Renfert monitor; 2D monitor with resolution 1920x1080 (Full HD).

\*\*\* Renfert CONNECT stick is not available in all regions. Please, check availability here: [www.renfert.com/connect](http://www.renfert.com/connect)

### Details

- Free movement in 4x mode: significant height differences (approx. 70 mm) do not lead to any loss of sharpness.
- Image storage in the common JPEG image format.



Additional lower magnification for improved comprehensive viewing

### Technical data

Permissible mains voltage	120–230 V   100–240 V
Permissible mains frequency	50/60 Hz
Power consumption	12 W
Power consumption (Monitor)	25 W
Video and screen resolution	Full-HD (1920 x 1080)
Focus distance (from camera lens to object)	430 mm // 17"
Magnification	4-fold   15-fold   20-fold
Intensity of lighting	40000 lx
Memory capacity USB stick	4 GB
Number of storable images, approx.	600
Max. video sequence length, approx.	~15 min
Weight (Head)	0,7 kg // 1.54 lbs
Weight (incl. swivel arm)	2,1 kg // 4.63 lbs
Weight (Monitor incl. power supply)	3,65 kg // 8 lbs

### Ordering information

EASY view+ 100–240 V	No. 24004000
EASY view+ 3D 100–240 V	No. 24003000
EASY view 3D monitor	No. 900037808
EASY view+ Composite UV filter	No. 24004020
Accessories	<a href="http://www.renfert.com">www.renfert.com</a>