


Renfert

making work easy

Are you still scraping?

Remove temporary
cement easily



temp:ex

Temporary Cement Remover

temp:ex

Temporary Cement Remover

temp:ex is a solvent for zinc oxide based, dental temporary cements.

Advantages

- Working time savings through a significant reduction in manual working steps.
- Avoids damage to the restoration via instruments during manual working steps.
- Tested compatibility with standard dental materials.

Details

- Up to 50% faster cleaning performance using temp:ex in the SYMPRO denture cleaning device.
- Economical use due to dosage aid on bottle.

MD



Simply let it soak for approx. 10 minutes ...



Dental restoration before and after cleaning

Ordering information

temp:ex

250 ml (8.45 fl.oz.) No. 66000000