



Update Info

SIMPLEX slice studio

Type:	SIMPLEX slice studio
Version:	3.1.25245.1
Hardware Kompatibilität:	Windows 10 & 11

Die Software-Aktualisierung enthält folgende Komponenten:

NEU

- Neues Layout / Menü
- Projektbibliothek
- Indikationserkennung und automatische Positionierung auf der Druckplatte
- Layoutfunktion (Positionierhilfe)
- Anzeige & Reparatur kleinerer Fehler im STL
- Kein automatisches Öffnen des G-Code Viewers mehr
- Hinterlegung neuer Materialien
- Integrierte Software-Update Funktion

Um die bestmögliche Leistung zu erhalten, sollten Sie die SIMPLEX slice studio stets auf dem neuesten Stand halten und Software-Aktualisierungen beachten.



Update Info

SIMPLEX slice studio

Type:	SIMPLEX slice studio
Version:	3.1.25245.1
Hardware compatibility:	Windows 10 & 11

The software update contains the following components:

NEW

- New layout / Menu
- Project library
- Indication recognition and automatic positioning on the printing plate
- Layout function (positioning aid)
- Display & Repair of minor errors in the STL file
- No more automatic opening of the G-Code Viewer
- Addition of new materials
- Integrated software update function

To get the best possible performance, always keep SIMPLEX slice studio up to date and consider software updates.