



Renfert Polishing Pastes

Overview of Renfert Polishing Pastes with recommended Polishing Tools

Item			Recommended Polishing Tool		
Polishing Paste	with Dia- mond Grit	Material/ Range of Application	Processing/ Preliminary Polishing	Polishing	High Lustre Polishing
 Renfert Polish 3D-printed materials No. 5104000	✓	3D printed materials			
 Renfert Polish hybrid materials No. 5103000	✓	Hybrid materials			
 Renfert Polish all-in-one No. 5100100	✓	Lithium disilicate, zirconium dioxide, veneering ceramics, high performance polymers			
 Renfert Polish ZrO ₂ No. 5101100	✓	Zirconium dioxide			
 Renfert Polish LiSi ₂ No. 5102100	✓	Lithium disilicate			
 Brinell L No. 5190001	✓	Precious PFM alloys, ceramic			
 Opal L No. 5200001		Veneering composites			
 Saphir No. 5150000		Precious metal alloys			
 Saphir No. 5150000		Cobalt-chrome alloys			
 Universal polishing paste No. 5131000		PMMA, resins			

*Included in Brinell L No. 5190001 scope of delivery

Polishers for Handpiece



Polishers for Lathe



Renfert Polishing Kits

Renfert Polish ZrO₂ starter kit

The Renfert Polish ZrO₂ starter kit combines the essential components for processing zirconium dioxide restorations into a coordinated set. The polishing paste and polishers are tailored to the material and typical indications, helping users carry out the polishing procedure in a structured, controlled manner to achieve reproducible results.

No. 510 1200



Renfert Polish all-in-one starter kit

The Renfert Polish all-in-one starter kit is a specially formulated polishing set for extraoral high-gloss polishing of dental CAD/CAM restorations. It combines the versatile Renfert Polish all-in-one diamond polishing paste with selected polishing instruments for pre-polishing, polishing, and high-gloss polishing. The kit supports material-specific surface finishing of various dental materials, such as veneering ceramics, zirconium dioxide, lithium disilicate, high-performance resins, and other ceramic materials.

No. 510 0103



Renfert Polish 3D printed materials starter kit

The Renfert Polish 3D-printed materials starter kit is specially developed for polishing 3D-printed materials. It contains a polishing paste developed specifically for this purpose, as well as the appropriate polishers, which are perfectly suited to the respective applications to ensure a flawless finish on filament-printed restorations made of PMMA and PETG, as well as on crowns and bridges made of resin.

No. 510 4100

