


Renfert

making work easy



Kolinsky evo

NEW Bionic hair brush

NEW Kolinsky evo

Bionic hair brush

Just like the popular range of Kolinsky brushes, Kolinsky evo, with its Bionic Hair, combines the benefits of a synthetic bristle brush with those of a natural bristle brush. The special molecular structure of the bristles allows moisture to be absorbed and is distinguished by outstanding longevity. The fineness, elasticity, and resilience of the brush tip allow targeted application of ceramic.

Advantages

- Outstanding longevity thanks to innovative Bionic Hair.
- Optimal moisture management thanks to the special molecular structure of the fibers.
- 6 different sizes for a range of individual applications.
- The streamlined design and low weight ensure convenient handling.

Details

- The fine tip with durable tension and elasticity impresses during targeted application of ceramic.
- Long service life without separation or “bunching” thanks to Bionic Hair.

Tip


If the brush tip bends inadvertently (e.g., due to incorrect storage), it can be held briefly in hot water at approx. 90°C. This automatically causes the bristles of the brush to straighten so that the brush tip returns to its original shape.



High-performance brush with Bionic Hair

Ordering information

Kolinsky evo size 1	2 pieces	No. 17243001
Kolinsky evo size 1/0	2 pieces	No. 17243110
Kolinsky evo size 2	2 pieces	No. 17243002
Kolinsky evo size 4	2 pieces	No. 17243004
Kolinsky evo size 6	2 pieces	No. 17243006
Kolinsky evo size 8	2 pieces	No. 17243008
Kolinsky evo set	with 6 brush sizes (1, 2, 4, 6, 8, 1/0)	No. 17243000



Kolinsky evo – in 6 sizes

