



making work easy

Press release

Hilzingen, Germany – April 2025

Small, yet a SILENT! Equipment specialist Renfert enters the dental office

SILENT XS: The mobile suction unit providing a cleaner working environment

Renfert (Hilzingen, Germany), a dental equipment specialist, presents the SILENT XS, the first fully mobile suction solution for dental offices and in-house laboratories. Launched at IDS 2025, this product innovation combines maximum flexibility with high suction power - for enhanced efficiency, clean air and patient comfort.

The SILENT XS is the response to the demand for flexible dust collectors designed for dental offices and in-house laboratories. Featuring LED lighting, H13/HEPA filtration and particularly quiet operation, it enables clean work - exactly where it is needed. This compact suction unit is ideal for milling temporary restorations, splints, crowns or dentures independently of a fully equipped laboratory workstation.

Maximum flexibility meets proven Renfert technology

SILENT XS shows its true potential in practical use. Cordless and battery-powered, this device can be moved flexibly between different locations and captures milling dust right at the source. SILENT XS replaces any improvised solutions with a professional mobile suction unit. The integrated LED lighting provides a clear view of the work area and enhances precision for intricate work. The H13/HEPA filter material is the SILENT XS's centerpiece. High quality glass fiber fleece offers a large filter surface in compact dimensions - for effective filtration and a long service life. The bagless design and smooth, easy-to-disinfect surface also meet modern hygiene standards.

More than a mobile suction unit, the SILENT XS is an investment in efficient workflows and a healthier working environment. Being a specialist in laboratory suction systems, Renfert transfers decades of experience directly into the dental office. With all its details, SILENT XS impresses anyone who values efficient and clean working. The result is unmistakably Renfert: sophisticated technology making everyday work noticeably easier and thus fulfilling the promise: "making work easy".

2,105 characters

More information

Renfert GmbH | Untere Giesswiesen 2 | 78247 Hilzingen | www.renfert.com

Press contact

Alberta Domnick | Tel. +49 7731 8208 898 | E-Mail: Alberta.Domnick@renfert.com